

Ambu Uni-Suction: Portable

Article no. 1002608



250,00 Euro Plus VAT.

Ambu Uni-Suction: Portable suction device specifically for pharyngeal and tracheal suction. Operation as foot pump or via O2-/Air-connector, second-hand

ALLIED L190-GR portable

Article no. 1002835



130,00 Euro Plus VAT.

ALLIED L190-GR portable suction pump, 2 litre glass, 0.7 bar/ 16 litre without battery, second-hand. This item is intended solely for veterinary medicine and export.

ATMOS RP 25 Absaugpumpe mit

Article no. 1000105



490,00 Euro Plus VAT.

ATMOS RP 25 Absaugpumpe mit 1x 4 Liter Absaugglas, 25 Liter/min., -0,8 bar, gebraucht

Atmos RP45 fahrbare Absaugpumpe

Article no. 1002415



590,00 Euro Plus VAT.

Atmos RP 45 fahrbare Absaugpumpe, 1x 2,5 Liter Absaugglas. Regelbare Saugleistung. 30Liter/min, max. -0,8bar, gebraucht

ATMOS S 031 Thorax: Drainage

Article no. 1001219



1.350,00 Euro Plus VAT.

ATMOS S 031 Thorax: Drainage suction device for mobile chest drainage. The device is portable, battery operated and is fitted with an electronic monitoring device. suction output: 3.0 +/- 0.3 l/min., vacuum max -60 mbar, 1-litre secretion container, demonstration device in mint condition, second-hand

BRAUN Vakufix, surgical

Article no. 1001716



190,00 Euro Plus VAT.

BRAUN Vakufix, surgical suction chamber, second-hand, Braun disposable collection bags, N EW

DRÄGER Vakuum-Absaugpumpe, 2x

Article no. 1000119



290,00 Euro Plus VAT.

DRÄGER Vakuum-Absaugpumpe, 2x 5 Liter Behälter, Überlaufglas, zum Anschluß an der Vakuumsteckdose, gebraucht

DRÄGER Vakuum-Absaugung

Article no. 1002953



180,00 Euro Plus VAT.

DRÄGER Vakuum-Absaugung zum Anschluß an eine zentrale Vakuumanlage. Zur Befestigung an einer Einheitschiene, regelbares Vakuum, neuwertig, ge braucht

DRÄGER bronchial suction pump

Article no. 1001222



90,00 Euro Plus VAT.

DRÄGER bronchial suction pump with injector (operated by O2 or compressed air), second-hand

DRÄGER bronchial suction pump

Article no. 1001223



140,00 Euro Plus VAT.

DRÄGER bronchial suction pump with injector (operated by O2 or compressed air) portable, with 1-litre MEDELA FMS bottle (disposable system) with replacement bag, second-hand

DRÄGER bronchial suction

Article no. 1001224



90,00 Euro Plus VAT.

DRÄGER bronchial suction pump, vacuum-operated, second-hand

DRÄGER portable secretion

Article no. 1001225



50,00 Euro Plus VAT.

DRÄGER portable secretion bottle pack, second-hand

Ejector for bronchial

Article no. 1001233



90,00 Euro Plus VAT.

Ejector for bronchial suction, second-hand

DRÄGER Ejector for Suction

Article no. 1002458



180,00 Euro Plus VAT.

DRÄGER Ejector for Suction incl. O2-pressure reducer and stepless flow-controller for O2, second-hand

DRÄGER Ejektorabsaugung NIST, VAC-betrieben

Article no. 1002995



155,00 Euro Plus VAT.

DRÄGER Ejektorabsaugung NIST, VAC-betrieben, neuwertig, gebraucht

FAHL tracheo, same as ATMOS

Article no. 1000091



350,00 Euro Plus VAT.

FAHL tracheo, same as ATMOS LC16, glass: 1,25 litres, suction power 16 l/min., very good condition, used

LAERDAL Suction Unit:

Article no. 1002701



160,00 Euro Plus VAT.

LAERDAL Suction Unit: Portable suction pump for mains and battery operation, secretion container 1 Litre, mains adapter is not original, therefore only for export and veterinary medicine

LAERDAL Premier Portable,

Article no. 1002482



290,00 Euro Plus VAT.

LAERDAL Premier Portable, battery operated suction, adjustable vacuum (max - 0,68 bar). 500 ml disposable container, incl. battery charger and carrying bag, second-hand

MEDAP P120 fahrbare

Article no. 1000088



770,00 Euro Plus VAT.

MEDAP P120 fahrbare Absaugpumpe, 2 x 4 Liter, gebr.

MEDAP P130, mobile suction

Article no. 1002213



790,00 Euro Plus VAT.

MEDAP P130, mobile suction pump with 4 litre suction container (NEW), incl. foot switch and spillover protection, second-hand, 1 litre container pictured

MEDAP P4020, portable

Article no. 1000081




220,00 Euro Plus VAT.

MEDAP P4020, portable aspirator, max. power output -600 mbar, second-hand

MEDAP Venta Carry 26 tragbare Absaugpumpe

Article no. 1003109



690,00 Euro Plus VAT.

MEDAP Venta Carry 26, tragbare Absaugpumpe Netz- und Akkubetrieb. Saugleistung 26 Liter / min., regelbares Vakuum bis -0,9 bar. Absaugbehälter mit 1 Liter, gebraucht

MEDELA Vario. Portable

Article no. 1002644



590,00 Euro Plus VAT.

MEDELA Vario. Portable suction pump, adjustable vacuum up to -0.8bar, 20l/min, 3 litre suction container (new containers and lid), second-hand

MEDELA Basic, mobile suction

Article no. 1000120



490,00 Euro Plus VAT.

MEDELA Basic, mobile suction pump, 2 litres, -0.8 bar, 1.5-litre container, second-hand.

MEDELA 015: milk pump

Article no. 1001228



280,00 Euro Plus VAT.

MEDELA 015: milk pump with 5-legged chasis, without accessories (can be purchased directly from the manufacturer), second-hand

MGE SAM12 suction pump with

Article no. 1001230



190,00 Euro Plus VAT.

MGE SAM12 suction pump with 1.75 litre container, -0.8 bar, portable, but not battery operated, adjustable vacuum, second-hand

MGE SAM12B suction pump with

Article no. 1001231



290,00 Euro Plus VAT.

MGE SAM12B suction pump with 2.0 Liter container, -0,8 bar, portable, but not battery operated, adjustable vacuum, second-hand

OLYMPUS KV-4: Portable

Article no. 1001229



780,00 Euro Plus VAT.

OLYMPUS KV-4: Portable suction pump H 240 x B 400 x T 410 mm, weight 10.5 kg, with 3-litre secretion glass with safety floater; vacuum -0.85 bar, output 20 L/min, second-hand*

OLYMPUS SSU-2, Keymed,

Article no. 1000336



420,00 Euro Plus VAT.

OLYMPUS SSU-2, Keymed, endoskopische Absaugpumpe, gebr.

VALLEYLAB OptiMumm mit Rapid Vac Sensor

Article no. 1002872



980,00 Euro Plus VAT.

VALLEYLAB OptiMumm: Rauchabsauger mit Rapid Vac Sensor E 3600, gebraucht

WEYER Minivac Acu S, portable

Article no. 1001232



290,00 Euro Plus VAT.

WEYER Minivac Acu S, portable battery operated suction pump (new battery) 500-ml suction glas, good condition, second-hand