

MELAG 23 blau, mit zwei

Article no. 1000140



1.590,00 Euro Plus VAT.

MELAG 23 blau, mit zwei Tablett, Baujahr 1999, 2 Programme, überholt, Innenvolumen 19 Liter, Innendurchmesser 23cm, Kammertiefe: 46 cm, gebr.

Melag Euroklav 23S, year of

Article no. 1001239



3.200,00 Euro Plus VAT.

MELAG Euroklav 23S, year of manufacture 2003, chamber size 25x47 cm, 22 litres. 3 programmes: universal programme, prion programme, express programme and 1 test programme. Connection for batch printer (not included in the price), second-hand. The device has been overhauled and is in a very good condition.

## MELAG 24 table-top autoclave

Article no. 1001238



1.290,00 Euro Plus VAT.

MELAG 24 table-top autoclave chamber diameter 23 cm, length 45 cm, 2 programmes 1 bar (121°C), 2 bar (134°C), the device has been overhauled, second-hand

MELAG 29 steam autoclave,

Article no. 1000148



1.190,00 Euro Plus VAT.

MELAG 29 steam autoclave, boiler unit 14.5 litres, chamber diameter: 22cm, length: 40 cm. including tray guidance frame and one tray, 2 bar/134 °C, 25 min. The device has been overhauled, second-hand

## MELAG 30 Autoklav Druck-Dampf Sterilisator

Article no. 1003114



1.490,00 Euro Plus VAT.

MELAG 30 Autoklav/Druck-Dampf Sterilisator, kammergröße 23 x 25cm , Sterilisation bei 2 bis 2,5 bar (135-138°C), Desinfektion unverpackter Güter 0,2-0,4bar (105-109°C), das Gerät eignet sich vor allem zum einsatz in der Dentalmedizin, gebraucht

## Melag Vacuklav 44B Baujahr 2005

Article no. 1003098



4.900,00 Euro Plus VAT.

Melag Vacuklav 44B, Baujahr 2005 5 verschiedene Programme. inkl. Dokumentationseinheit MelaFlash, Tablettführungsgestell und 2 Tablett B-Klasse Autoklav. Zum Betrieb empfehlen wir eine Wasseraufbereitungsanlage MELADem 40 (neu). Preis EUR 390, gebraucht

Aesculap Transporttrolley für Sterilboxen 70x34cm

Article no. 1003039



90,00 Euro Plus VAT.

Aesculap Transporttrolley für Sterilboxen 70x34cm, Höhe 93cm Edelstahl, 2 Rollen feststellbar, gebraucht

W & H Lisa Mdist 4. Water

Article no. 1002827



120,00 Euro Plus VAT.

W & H Lisa Mdist 4. Water Distillation Device 580 W. produced in 5 Hrs. 4 Litre Aqua dist. Year of Manufacture 2004, second-hand

Cartridge desalinator 20

Article no. 1001241



Dr. Wilfried Müller GmbH  
Gebrauchte medizinische Geräte

998,00 Euro Plus VAT.

Cartridge desalinator 20 with guide value meter: basic capacity of approximately 10° d GSG: 2800 L; H 570 mm, with guide value meter 695 mm, diameter 240 mm; with stainless steel quick closures for easy cartridge replacement, N EW

MARTIN instrument sieve,

Article no. 1001242



20,00 Euro Plus VAT.

MARTIN instrument sieve, stainless steel, dimensions: Height: 5 cm x width: 25cm, x depth: 25 cm, second-hand

instrument sieve, stainless

Article no. 1001243



20,00 Euro Plus VAT.

instrument sieve, stainless steel, width: 48 cm, x depth: 25 cm, x height: 5 cm, good condition, second-hand

MARTIN Sterilisiercontainer 1/2 STE 30 x 30 x 30 cm

Article no. 1002951



55,00 Euro Plus VAT.

MARTIN Sterilisiercontainer 1/2 STE 30 x 30 x 30 cm, Filter im Deckel und Boden Edelstahl, gebraucht

**MEMMERT S15 Heissluftsterilisator 220°C**

Article no. 1003063



**390,00 Euro** Plus VAT.

MEMMERT S15 Heissluftsterilisator 220°C, Zeitschaltuhr 120 min, 2 Alutablets, guter Zustand, gebraucht

Steriliser container 1/4 STE

Article no. 1002660



Dr. Wilfried Müller GmbH  
Gebrauchte medizinische Geräte

45,00 Euro Plus VAT.

Steriliser container 1/4 STE (15x30x30cm), sec.h and

**Steriliser container 1 STE**

Article no. 1000157



**90,00 Euro** Plus VAT.

Steriliser container 1 STE (30x30x60 cm), second-hand

**Steriliser container 1/2 STE**

Article no. 1000158



**45,00 Euro** Plus VAT.

Steriliser container 1/2 STE (30x30x60 cm), second-hand

Aluminium Drum Container with

Article no. 1002569



25,00 Euro Plus VAT.

Aluminium Drum Container with filters in lid and base D=26cm Height 23cm

Aluminium Drum Container with

Article no. 1002570



25,00 Euro Plus VAT.

Aluminium Drum Container with filters in lid and base D=32cm Height 24cm

Aluminium Drum Container with

Article no. 1002571



25,00 Euro Plus VAT.

Aluminium Drum Container with filters in lid and base D=21cm Height 20cm

Aluminium Drum Container with

Article no. 1002572



25,00 Euro Plus VAT.

Aluminium Drum Container with filters in lid and base D=21cm Height 15cm

Aluminium Drum Container with

Article no. 1002573



25,00 Euro Plus VAT.

Aluminium Drum Container with filters in lid and base D=16cm Height 14cm

**BELIMED Spülkorb, gebraucht**

Article no. 1002919



**290,00 Euro** Plus VAT.

BELIMED Spülkorb, gebraucht

**BELIMED Spülkorb, gebraucht**

Article no. 1002922



**290,00 Euro** Plus VAT.

BELIMED Spülkorb, gebraucht

**MARTINIT 18/8 Handbürstenständer**

Article no. 1002926



**158,00 Euro** Plus VAT.

MARTINIT 18/8: Handbürstenständer, autoclavierbar, Edelstahl, gebraucht