

MARTIN Elektrotom 40 monopolares HF-Chirurgiegerät

Article no. 1003004



490,00 Euro Plus VAT.

MARTIN Elektrotom 40, monopolares HF-Chirurgiegerät zur Koagulation und Fulguration, Leistung max.: 40 Watt, gebraucht

MIKROMED MK 1 monopoloares HF-Gerät 50W

Article no. 1002879




980,00 Euro Plus VAT.

MIKROMED MK 1 monopoloares HF-Gerät 50W, sehr guter Zustand, gebraucht

ERBE ICC 200 with EndoCut function, monopolar 200W, bipolar

Article no. 1002125



3.190,00 Euro Plus VAT.

ERBE ICC 200 with EndoCut function, monopolar 200W, bipolar 120W. Very good condition, second-hand

ERBE ICC 300: monopolar cutting: 300W, coagulating: 300W, bi

Article no. 1000496



3.300,00 Euro Plus VAT.

ERBE ICC 300: monopolar cutting: 300W, coagulating: 300W, bipolar: 120W, dimensions: 410x152x368mm, weight: 10kg, demonstration device, very good condition, second-hand ALL HF SURGICAL DEVICES ARE DELIVERED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BI-POLAR ELECTRODE, 3 DIFFERENT INSERTS)

ERBE ICC 350, HF surgical

Article no. 1002177



3.600,00 Euro Plus VAT.

ERBE ICC 350, HF surgical device, monopolar cutting 300 W, bipolar 120 W, dimensions 410 x 152 x 368 mm, good condition, second-hand. ALL HF SURGICAL DEVICES ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS)

ERBE Beamer 2: Gas-operated

Article no. 1001309



1.560,00 Euro Plus VAT.

ERBE Beamer 2: Gas-operated argon gas source for argon plasma coagulation. The equipment with a device that combines the activation of the beamer 2 and the high-frequency surgical device. Can be combined with devices ACC 450, 451, 430 and ICC 200 (only model 10128-023), 300, 350, see second-hand

AESCULAP GK 160

Article no. 1002143



1.290,00 Euro Plus VAT.

AESCULAP GK 160 Arthroresektor: Monopolar HF surgical device, output power 170 Watt, frequency 0.5 MHz, weight 4.6 kg, H 150 x W 305 x D 215 mm, second-hand. ALL HF SURGICAL DEVICES ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS)

OLYMPUS PSD-10 HF surgical

Article no. 1001989



980,00 Euro Plus VAT.

OLYMPUS PSD-10 HF surgical device, for endoscopic work, incl. accessories, second-hand. ALL HF SURGICAL DEVICES ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS)

OLYMPUS PSD-20 HF surgical

Article no. 1001990




1.980,00 Euro Plus VAT.

OLYMPUS PSD-20 HF surgical device, for endoscopic work, incl. accessories, second-hand. ALL HF SURGICAL DEVICES ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS)

OLYMPUS PSD-30: HF-Surgery

Article no. 1002549



2.900,00 Euro Plus VAT.

OLYMPUS PSD-30: HF-Surgery Device for endoscopic use, second-hand ALL HF SURGICAL DEVICES ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS)

OLYMPUS UES 30, HF-Surgical

Article no. 1002355



3.200,00 Euro Plus VAT.

OLYMPUS UES 30, HF-Surgical Device for endoscopy and general surgery. Monopolar 300W, Bipolar 120W, second hand ALL HF-SURGICAL DEVICES ARE DELIVERED WITH ACCESSORIES (NEUTRAL PLATES, MONO AND BIPOLAR ELECTRODES, 3 DIFFERENT INSERTS)

SÖRING MBC 601 HF device

Article no. 1002713



2.100,00 Euro Plus VAT.

SÖRING MBC 601 HF device controlled by processor, year of manufacture 2004, monopolar 350W, bipolar Cut 130 W, bipolar Coag. 100W, second-hand

VALLEYLAB Force 2,

Article no. 1002978



1.280,00 Euro Plus VAT.

VALLEYLAB Force 2 CEM, high-frequency surgical device with floating output in microprocessor technology for cutting and coagulating, incl. neutral electrode monitoring REM with automatic limit value adjustment. HF output in wattas can be adjusted steplessly (digital) up to a max. of 100 W - performance optimised via Closed Loop Circuit construction principle. Bipolar coagulation can be adjusted up to 50 watt, second-hand, incl. foot switch. ALL HF SURGICAL DEVICES ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS)

VALLEYLAB FORCE 30: Digitales HF-Chirurgiegerät

Article no. 1002876



1.150,00 Euro Plus VAT.

VALLEYLAB FORCE 30: Digitales HF-Chirurgiegerät, 300 Watt, monopolar, 50 Watt bipolar, mit folgenden Sonderfunktionen: pure cut, blend, desiccate, bipolar: standard precise, gebraucht

VALLEYLAB FORCE 40: Digital

Article no. 1002742



1.250,00 Euro Plus VAT.

VALLEYLAB FORCE 40: Digital HF-Surgical Device, 300 Watt, mono and bipolar, with the following special functions: cut, blend1, blend 2, deiccate, fulgurate, spray standard bipolar, precise bipolar, second-hand

Valleylab Force 300,

Article no. 1002820



2.290,00 Euro Plus VAT.

Valleylab Force 300, HF-Surgical Device, monopolar 300 Watt, bipolar 70W, very good condition, second-hand. ALL HF SURGICAL DEVICES ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS)

Valleylab Force EZ prozessorgesteuertes HF

Article no. 1003012



3.490,00 Euro Plus VAT.

Valleylab Force EZ, prozessorgesteuertes HF Gerät Monopolar Cut 300W, Monopolar COag. 120W, Bipolar 70W , gebraucht

ERBE IES 300, smoke gas

Article no. 1002252



1.790,00 Euro Plus VAT.

ERBE IES 300, smoke gas suction for HF surgery, including a new main filter, for use with all current HF surgical devices, second-hand

Foot switch (2-pedal) for HF

Article no. 1001316




280,00 Euro Plus VAT.

Foot switch (2-pedal) for HF device (connection for MARTIN or ERBE), explosion-protected, second -hand

BOWA ErgoPen handle, narrow,

Article no. 1002298



98,00 Euro Plus VAT.

BOWA ErgoPen handle, narrow, 2 keys for Valleylab devices, second-hand

BOWA monopolar electrode

Article no. 1001317



125,00 Euro Plus VAT.

BOWA monopolar electrode handle with dual keys suitable for ERBE T series including 4.5 m connection cable, new

BOWA Monopolar Electrode

Article no. 1002397



98,00 Euro Plus VAT.

BOWA Monopolar Electrode Handle with dual keys suitable for ERBE ACC/ICC incl. 4.5 m connection cable, new

BOWA monopolar electrode

Article no. 1001318



98,00 Euro Plus VAT.

BOWA monopolar electrode handle for Martin HF devices, 2 keys, shaft 4 mm, cable 4.5 metres, n ew

Bipolar connection cable for

Article no. 1001319



63,00 Euro Plus VAT.

Bipolar connection cable for Martin HF devices, length: 4.5 m, (287040), new

Connection cable for

Article no. 1001320



63,00 Euro Plus VAT.

Connection cable for Martin/Berchtold (295030) neutral electrode, new

Rubber neutral plate (245x145

Article no. 1001321



190,00 Euro Plus VAT.

Rubber neutral plate (245x145 mm), without cable (242003), for Martin/Berchtold, new

Bipolar connection cable for

Article no. 1001322



63,00 Euro Plus VAT.

Bipolar connection cable for Erbe T series /ACC/ICC, (101-040), new

Connection cable for ERBE

Article no. 1001323



63,00 Euro Plus VAT.

Connection cable for ERBE neutral electrode (101-002), new

Rubber neutral plate for ERBE

Article no. 1001324



190,00 Euro Plus VAT.

Rubber neutral plate for ERBE T series ACC/ICC 245x145 mm, without cable, (193-008), new

Spherical monopolar

Article no. 1001325



8,00 Euro Plus VAT.

Spherical monopolar electrode, 4 mm shaft, 4 mm diameter, (500-021), new

Monopolar needle electrode, 4

Article no. 1001326



8,00 Euro Plus VAT.

Monopolar needle electrode, 4 mm shaft, 4 mm diameter, (500-011), new

Monopolar wire sling

Article no. 1001327



8,00 Euro Plus VAT.

Monopolar wire sling electrode, 10 mm diameter, 4 mm shaft (500-015), new

Monopolar knife electrode, 25

Article no. 1001328




8,00 Euro Plus VAT.

Monopolar knife electrode, 25 mm, 4 mm shaft, new

Bipolar tweezers, 195 mm, tip

Article no. 1001329




119,00 Euro Plus VAT.

Bipolar tweezers, 195 mm, tip 1 mm (605-040), for Erbe and Martin, new

Bipolar tweezers, 110 mm, tip

Article no. 1001330



119,00 Euro Plus VAT.

Bipolar tweezers, 110 mm, tip 0.5 mm (605-0420), new

GRIESHABER bipolar eye

Article no. 1001331



98,00 Euro Plus VAT.

GRIESHABER bipolar eye tweezers (titanium), 110 mm, tip 0.5 mm, new

Additional accessories

Article no. 1001333



Dr. Wilfried Müller GmbH
Gebrauchte medizinische Geräte

0,00 Euro Plus VAT.

Additional accessories available on request