

## DRÄGER Transport - Incubator

Article no. 1002619



2.500,00 Euro Plus VAT.

DRÄGER Transport - Incubator 5400 on guidance frame, with O<sub>2</sub> - connection for central gas supply and small bottles, electric connection for 220 Volt and 12 Volt (auto battery), incl. bronchial suction (ejector), second-hand

## DRÄGER Transport-Incubator

Article no. 1002408



3.500,00 Euro Plus VAT.

DRÄGER Transport-Incubator 5400 on Transport-Lafette, for external and internal O<sub>2</sub>- and air supply, bottle connections for O<sub>2</sub>, power supply via mains connection or car battery, with Babylog 2 ,good condition, second-hand

## DRÄGER Transport-Inkubator auf Transport-Lafette

Article no. 1002973



4.800,00 Euro Plus VAT.

DRÄGER Transport-Inkubator 5400 auf Transport-Lafette, für externe und interne O<sub>2</sub>- und Air-Versorgung, Flaschen-Anschlüsse für O<sub>2</sub> + Air, Stromversorgung über Netzanschluß oder Autobatterie, mit Babylog 2000, Höhe: 129 cm x Breite: 115cm x Tiefe: 51cm, inkl. Bronchial Absaugung (O<sub>2</sub> Injektor), Fahrgestell mit zwei Schubladen und 1 Fach; mit Schwanenhals-Leuchte, sehr guter Zustand, gebraucht

**DRÄGER incubator 8000, with**

Article no. 1001351



**2.980,00 Euro** Plus VAT.

DRÄGER incubator 8000, with integrated cold system and temperature (heating output 400 watt), O<sub>2</sub> + humidity measurement, electrical height adjustment, tilted position of the laying area: right and left continuously up to 10° inclination angle, dimensions: H 1498 x W 82 x D 437 mm, weight 104 kg, second-hand

## Dräger Caleo, mobile baby

Article no. 1002720



5.900,00 Euro Plus VAT.

Dräger Caleo, mobile baby incubator. Heat therapy controlled by air intake temperature control or skin temperature control. Air humidification, O<sub>2</sub> concentration. Swivelling lying surface. Incl. monitoring for air temperature, skin temperature relative humidity, O<sub>2</sub> concentration Dräger Caleo, second-hand

**DRÄGER Babytherm 8000:**

Article no. 1003049



1.900,00 Euro Plus VAT.

**DRÄGER PP 800 resuscitation**

Article no. 1000566



**980,00 Euro** Plus VAT.

DRÄGER PP 800 resuscitation unit for newborn infants with heat radiator, second-hand

**DRÄGER Babytherm 8000, mobile**

Article no. 1002485



**1.250,00 Euro** Plus VAT.

DRÄGER Babytherm 8000, mobile thermal bed without acrylic glass cover, second-hand

## DRÄGER Babytherm 4200 thermal

Article no. 1001356



650,00 Euro Plus VAT.

DRÄGER Babytherm 4200 thermal bed (can be converted to a multipurpose intensive care bed by attaching various devices), with continuously controllable mattress heating (180 W), swivels on cross axis, H 1170 x W 1025 x D 555 mm, with adjustable radiator and integrated drawer unit, good condition, second hand

## DRÄGER Babytherm 4200 thermal

Article no. 1002133



450,00 Euro Plus VAT.

DRÄGER Babytherm 4200 thermal bed (can be converted to a multipurpose intensive care bed by attaching various devices), with continuously controllable mattress heating (180 W), swivels on cross axis, H 1170 x W 1025 x D 555 mm, with integrated drawer unit, good condition, second hand

**AIR SHIELDS Resuscitaire,**

Article no. 1002322



**1.900,00 Euro** Plus VAT.

AIR SHIELDS Resuscitaire, Radiant Warmer, RW81-1E: Mobile thermal bed with integrated resuscitation devices, heating module (3 operating modes: Preheating, skin temperature control, manual mode), microprocessor-control heating lamp (can be rotated 180°), integrated Apgar timer, removable side walls, chassis with integrated drawer and bronchial suction, second-hand\*

**DRÄGER Photo therapy 800, on**

Article no. 1001357



**590,00 Euro** Plus VAT.

DRÄGER Photo therapy 800, on chassis, for the light irradiation of premature or newly born infants in the incubator or in the infrared-heated thermal bed to reduce the serum-bilirubin concentration (phototherapy), frequency: 50 Hz, output: 230W, current: 2.2 A, lamp housing 54 x 25.5 x 20.5cm, can be inclined on tripod up to 30°, tripod height adjustable from 107 - 162 cm, tripod base 66x55x10 cm, seco nd-hand

**DRÄGER Photo Therapie 8000, on chassis, for light radiation**

Article no. 1002339




**790,00 Euro** Plus VAT.

DRÄGER Photo Therapie 8000, on chassis, for light radiation of premature infants or newly born children within the incubator or within the infrared heated thermal bed to reduce the serum bilirubin concentration (photo therapy), frequency: 50Hz, output: 230W, current: 2.2A, lamp housing 54 x 25.5 x 20.5cm, inclination on tripod possible up to 30°, tripod height adjustable between 107 and 162cm, tripod foot 66x55x10cm, second-hand

**DRÄGER Photo Therapy 4000,**

Article no. 1002237



**690,00 Euro** Plus VAT.

DRÄGER Photo Therapy 4000, 120 W, 4 tubes and neon lighting that can be switched separately for attachment to an incubator, second-hand

## MEDELA phototherapy

Article no. 1002238



490,00 Euro Plus VAT.

MEDELA phototherapy spotlight, mobile, 4 x 18 Watt, tubes can be switched separately from each other, second-hand

**BEGA Boom 818: Infrared beam**

Article no. 1002421



**160,00 Euro** Plus VAT.

BEGA Boom 818: Infrared beam 150 Watt, high-quality completely closed alu-casing, 360° rotatable 10° swivelling, casing stays cool in operation, second-hand

### Hanulux Hanautherm, red

Article no. 1002175




690,00 Euro Plus VAT.

Hanulux Hanautherm, red light radiant heater 1100 W for ceiling installation. For rooms up to 300 cm in height, second-hand

**DRÄGER baby scales (sliding**

Article no. 1001358



**180,00 Euro** Plus VAT.

DRÄGER baby scales (sliding weight scales) for DRÄGER incubator, second-hand