

**Dr. Wilfried Müller GmbH**

Lechstr. 13-15

D-86931 Prittriching

Phone 0049-(0)8206-9624-0

Fax 0049-(0)8206-9624-29

E-Mail [info@med-geraete.de](mailto:info@med-geraete.de)

Internet [www.med-geraete.de](http://www.med-geraete.de)



## Second-hand equipment for human and veterinary medicine



For all second-hand equipment: 12-month guarantee

For all new equipment: 24-month guarantee

8-year supply guarantee for spare parts and accessories

Visit us in our showrooms (appointment required) in  
Prittriching (between Augsburg and Munich)












For the latest version of our catalogue (mostly with color photos), please see in the Internet at:

**[www.med-geraete.de](http://www.med-geraete.de)**












- All equipment is second-hand, fully-functional and in good condition, unless stated otherwise.
- For all second-hand equipment: 12-month guarantee
- For all new equipment: 24-month guarantee
- Special conditions apply for traders.
- All prices are net prices, subject tot 19% VAT
- If exporting items within the EU, the German VAT is waived when you provide your VAT ID number.
- If exporting items to non-EU countries, the German VAT is waived if we dispatch the products.
- If exporting items to non-EU countries we shall refund the German VAT upon presentation of a valid export certificate that is issued on the name of Dr Wilfried Müller GmbH.
- All offers subject to change and subject to prior sale.
- Transport costs to be shared equally between buyer and seller.
- The acceptance of our terms and conditions must be the basis of any sales agreement.












<b>Suction pumps</b> .....	4	Oxygen bubblers.....	93
<b>Autoclaves and sterilizers</b> .....	6	Vapors (anaesthetic vaporisers).....	94
<b>Blood pressure measurement</b> .....	8	Pressure reducers und central gas supply.....	99
<b>Defibrillators</b> .....	8	Flowmeters.....	102
<b>EEG</b> .....	12	Compressors, concentrators.....	103
<b>ECG</b> .....	12	<b>OP and investigation lights</b> .....	103
<b>Flexible endoscopes</b> .....	13	Bulbs.....	104
Bronchoscopes.....	13	OP lights.....	104
Bronchoscopes for horses.....	13	Transformers and UPS.....	107
Choledochoscopes.....	13	Investigation lights.....	108
Cystoscopes.....	14	<b>Operating tables and patient trolleys</b> .....	108
Duodenoscopes.....	14	Operating tables.....	109
Flexible micro endoscopes.....	14	Accessories for operating tables.....	109
Gastrosopes.....	14	Patient trolleys.....	114
Hysteroscopes.....	14	<b>Osteosynthesis</b> .....	114
Jejunoscope.....	14	<b>X-ray fluoroscopy</b> .....	114
Coloscopes.....	15	x-ray machines.....	114
Laryngoscopes.....	15	Accessories for x-rax machines.....	114
Uretroscopes.....	15	Mobile c-arms.....	115
<b>Endoscopy-cold light sources</b> .....	15	Accessories for mobile c-arms.....	115
Olympus.....	16	<b>Other</b> .....	115
Other.....	17	<b>Ultra sound devices</b> .....	116
Storz.....	18	Ultrasound systems.....	116
Wolf.....	19	Ultrasonic cleaners.....	116
<b>Endoscopy-rigid endoscopes</b> .....	19	Ultrasonic nebulizers.....	116
Rigid endoscopes.....	19	<b>Investigation</b> .....	116
Rigid micro endoscopes.....	19	<b>Video printer and recorder</b> .....	117
<b>Endoscopy-fittings</b> .....	19	<b>Dental treatment</b> .....	117
Endoscopy towers.....	19	<b>General Terms and Conditions of Business</b> .....	119
Video.....	19		
more fittings.....	23		
<b>Chassis and monitor carrier</b> .....	28		
<b>Repository</b> .....	31		
<b>Gynaecology</b> .....	33		
<b>HF surgery</b> .....	36		
<b>ENT</b> .....	41		
<b>Infusors perfusors</b> .....	41		
<b>Incubators</b> .....	45		
<b>Laboratory equipment</b> .....	46		
<b>Laser surgery equipment</b> .....	46		
<b>Medical antiques</b> .....	46		
Lights.....	46		
Furniture.....	46		
Anaesthesia equipment.....	46		
Other.....	47		
<b>Microscope</b> .....	48		
<b>Furniture</b> .....	53		
<b>Monitoring</b> .....	53		
Datex.....	53		
Dräger.....	57		
Hellige / Marquette.....	61		
Nellcor.....	62		
Siemens.....	62		
Other.....	63		
Accessories for monitors.....	66		
<b>Anaesthesia and ventilation</b> .....	70		
Ventilation equipment.....	71		
Gas provision and electricity supply.....	76		
Equipment fittings.....	77		
Breathing circuit and fittings.....	78		
Anaesthesia equipment.....	86		

Category	Article number	Product	Picture	Price
<b>Suction pumps</b>	1003597	OLYMPUS KV-5 suction pump with new suction container (1 liter), used		1.100,00 EUR
	1003235	MGE SAM 12, portable suction pump with 1.75 litre plastic container, max .-0.8 bar, used		350,00 EUR
	1002192	ATMOS LC16 portable suction pump 16 litres/min., 1.0 litre suction glass, adjustable vacuum, second-hand		230,00 EUR
	1001223	DRÄGER bronchial suction pump with injector (operated by O2 or compressed air) portable, second-hand, NIST for O2		230,00 EUR
	1001224	DRÄGER bronchial suction pump, vacuum-operated, second-hand		250,00 EUR
	1004886	DRÄGER holder for secretion bottle package or bronchus suction. The holder can be clamped to a standard rail, second-hand		58,00 EUR
	1002953	DRÄGER ejector exhaust system, air-operated, mount for standard rails, adjustable vacuum, used		155,00 EUR
	1002995	DRÄGER ejector exhaust system, VAC-operated, mount for standard rails, adjustable vacuum, used		155,00 EUR
	1004483	DRÄGER bronchial suction bottle made of polysulfone, NEW		49,00 EUR
	1004484	DRÄGER lid for suction bottle, used		25,00 EUR
	1004425	DRÄGER holder for bronchial suction, second-hand		28,00 EUR

1002458	DRÄGER Ejector for Suction incl. O2-pressure reducer and stepless flow-controller for O2, second-hand		90,00 EUR
1003315	DRÄGER ejector for bronchial suction, used		90,00 EUR
1004817	DRÄGER Sjector 2000: Bronchial aspirator		68,00 EUR
1004924	MEDAP AIR with 030623 safety glass incl. manometer, clamping device for standard rail and pressure hose and plug for wall tapping point Air		85,00 EUR
1003769	HERSILL Eurovac 50 surgical aspirator: mobile, powerful, 70 l/min, 2 x 4 liter secretion container, size 96 x 43 x 34 cm, 25 kg incl, bacteria filter and suction hose, with foot switch, NEW		1.350,00 EUR
1004882	LID for suction container 4 litres, suitable for Eurovac 50 suction container, NEW		39,00 EUR
1003921	Replacement bacteria filter for EUROVAC 50, NEW		9,80 EUR
1003922	Replacement 2.5m suction hose for EUROVAC 50, NEW		34,90 EUR
1004448	Suction container 4 litres for Eurovac, NEW		85,00 EUR
1004629	OLYMPUS MH-481 Filter for KV-5 (and other suction pumps)		19,50 EUR
1004917	MEDELA 077.0521 Fluid Management System rail clamp plastic, for secretion canister, used		22,00 EUR











**Autoclaves and sterilizers**

1005200	Suction cannula 4mm with illumination. Cold light cables from Storz can be attached. Length 220mm, used		125,00 EUR
1003751	MELAG Euroklav 23VS+, chamber size 45 cm x diameter 25 cm, 22 liter, load capacity 4 kg of instruments, 1 kg textiles, dimensions (WxHxD): 42.5 x 48.5 x 69 cm, the unit fits on a 60 cm tabletop, weight 45 kg, power consumption 2600 Watt, used		3.650,00 EUR
1002624	MELAG Vacuclav 40-B, class B autoclave kettle D = 25 cm, length 35 cm, 18 liters 3400W, incl. tray frame and 1 tray, overhauled, including MELADEM 40, used		3.490,00 EUR
1003551	MELAG Vacuklav 24B: Table-top autoclave, B class (corresponds to European standard EN 13060), with 1 tray. The unit requires a fixed water connection and a waste water connection for operation. Dimensions: diam. 25 x D 45 cm, capacity: 6 trays Tray dimensions: 42 x 19 x 2 cm Outer dimensions W 42.5 x H 48.5 x L 66 cm , second-hand		3.250,00 EUR
1003578	MELAG Melaflash CF-Card-Printer, for batch documentation at your Melag autoclave, including CF card, second-hand		199,00 EUR
1004549	MELAG Melanet Box, for connecting Melag autoclaves (Profi class, S class and other series) to the practice network, second-hand		480,00 EUR
1004154	MELAG MELAprint, batch printer for MELAG autoclaves, including power supply and serial interface cable, used		139,00 EUR
1005106	MELAG tray 29x19cm aluminium, suitable for Melag units such as Vacuclav 30B, Euroklav 29VS+, etc., second-hand		29,00 EUR
1005104	MELAG tray 42x19cm aluminium, suitable for Melag units such as Euroklav 23, MelaTronic 23 etc, second-hand		35,00 EUR
1004304	MELAG holding rack for foil-wrapped sterile goods, suitable for autoclaves, e.g. MELAG Vacuklav 23, Vacuklav 44B+ etc., dimensions 36 x 18cm, second-hand		65,00 EUR
1000158	AESCLAP Steriliser container 1/2 STE (30x30x60 cm), second-hand		60,00 EUR


1005153	MELAG ME01081 Boiler cleaning set for powerful cleaning of Autoclaves (kettle, door round, holder) removes greasy and oily residues. Set consists of spray bottle with 500 ml cleaner, microfibre cloth and paper towels, new.		39,00 EUR
1001241	Cartridge desalinator 20 with guide value meter: basic capacity of approximately 10° d GSG: 2800 L, H 570 mm, with guide value meter 695 mm, diameter 240 mm, with stainless steel quick closures for easy cartridge replacement, N EW		998,00 EUR
1004097	MELAG MelaDem 40 with MelaJet spray gun, used		380,00 EUR
1004012	MELAG MelaDem 40, demineralization system for Melag B and S class autoclaves, as new, unused		490,00 EUR
1003828	MELAG MELAdem 47: reverse osmosis system, incl. new filter and cartridges, supplies with expansion tank, used		590,00 EUR
1004829	ERMIS ER200.010E: Sterilisation container with sieve, W 45 x D 28 x H 10 cm, NEW		120,00 EUR
1005127	Instrument sieve stainless steel, width: 45cm x depth: 17cm x height: 6cm, for use in autoclave or WD, good condition, used		30,00 EUR
1005159	HALL E306 autoclave container, W38cm, D22cm, H12cm for compressed air-driven drill, oscillating saw, etc., second-hand		49,00 EUR
1000157	Steriliser container 1 STE (30x30x60 cm), second-hand		90,00 EUR
1003651	TOLZIN Steri-Test- package for Bowie-Dick test, used		58,00 EUR
1005250	MELAG insert rack for steam autoclaves. This rack is used for the safe sterilisation of foil-wrapped sterile goods. Dimensions: 37cmx18cm, height 10cm, used		65,00 EUR













## Blood pressure measurement

1003779	MELAG Steriliser 75, hot air steriliser for 2 trays/cassette 29x19cm, incl. timer and thermostat, second-hand		295,00 EUR
1002387	WELCH ALLYN Spot Vital Signs, measurement of SPO2 and blood pressure, on mobile rolling stand, second-hand		390,00 EUR
1004705	BOSO Roid II: Aneroid sphygmomanometer, NEW		25,00 EUR
1000204	CRITIKON OxyShuttle: Vital Data Monitor, measurement and digital display of blood pressure (NIBP) and pulse, incl. Oxytrak pulse oxymeter, second-hand		390,00 EUR
1003246	CRITIKON Dinamap XL blood pressure measurement device incl. pulse display, supplied with blood pressure cuff and pressure tube, used		490,00 EUR
1001264	CRITIKON Dinamap 8100 Vital Data Monitor, measurement and digital display of blood pressure (NIBP) and pulse, second-hand		390,00 EUR
1004367	DINAMAP Procure 300: Blood pressure monitor, display of systolic and diastolic blood pressure, pulse rate. Adjustable alarm limits. Battery or mains operation possible. Including mains charger, used		320,00 EUR
1001266	GE Dinamap PRO 300V2. Monitor for measurement of NIBP, pulse, SpO2 with battery function and printer, good condition, incl. chassis, used.		590,00 EUR
1001265	CRITIKON Dinamap Compact TS. Measurement of NIBP and SPO2, with battery function, good condition, second-hand		490,00 EUR
1004690	PW-Tourniquet, NEW		5,50 EUR












## Defibrillators





1002401	PHYSIO CONTROL Lifepak 20: Light and compact defibrillator with biphasic ADAPTIV technology semi-automatic/manual; AED mode, colour monitor, ECG, English keyboard, SW German (can be switched to English), delivery incl. Quick Combo electrode cable for disposable adhesive electrodes, delivery without hardpaddles, second-hand		2.890,00 EUR
---------	--	---	--------------



1004633	PHYSIO CONTROL Lifepak 20: Light and compact defibrillator with biphasic ADAPTIV technology semi-automatic / manual, colour monitor, ECG, with Quick Combo electrode cable for use of disposable adhesive electrodes, Spanish keyboard, language of menu guidance can be freely selected, delivery without hardpaddles, second-hand		2.650,00 EUR
1004941	PHILIPS Heartstart MRx, portable, biphasic defibrillator incl. battery. Power up to max. 200J, incl. AED, ECG. . Large intuitive color display, printer. The device is supplied with a connection cable for operation with disposable adhesive electrodes, software and language in English, delivery without hardpaddles, used		2.890,00 EUR
1005204	PHILIPS Heartstart XL+, biphasic defibrillator up to a maximum of 200 joules, delivery incl ECG cable, CO2 module, lead for disposable adhesive electrodes, lithium battery, Software language and keys in English, delivery without hardpaddles, second-hand		3.980,00 EUR
1000224	Philips Heartstart XL, bi-phasic defibrillator with AED function, incl. printer and battery power, English software, used		1.980,00 EUR
1004808	PHILIPS Heartstart XL, biphasic defibrillator with AED function, incl. printer and battery operation, software German, with pacemaker, second-hand		2.290,00 EUR
1004944	PHILIPS intracardiac ladle electrodes for Philips defibrillators (e.g. Heartstart XL). Diameter 45mm, the electrodes are autoclavable, second-hand.		490,00 EUR
1003567	PHYSIO CONTROL Lifepak 12, monophasic defibrillator, printer, table charger,, Including disposable adhesive paddles, Hardpaddles surcharge 450,00 Euro, Keyboard language English, Delivery without batteries, second-hand		1.300,00 EUR
1001282	PHYSIO CONTROL Lifepak 12, monophasic defibrillator with printer and monitor, good condition, incl. defi electrodes (disposable adhesive paddles), 1 pair of hard paddles extra charge 450.00 Euro, delivery without batteries, second-hand		980,00 EUR
1001283	PHYSIO CONTROL Lifepak 12, monophasic defibrillator, printer with writing width 6cm, table charger, SPO2. Tatatur in German, extra charge Hardpaddles EUR 450,00, delivery without battery used		1.400,00 EUR
1002551	PhysioControl Lifepak 12 biphasic, measurement of ECG, SPO2, NIBP, CO2. With printer 10cm width, carrying strap, including disposable adhesive paddles, hardpaddles surcharge 450,00 EUR, delivery without batteries, second-hand		3.200,00 EUR







1004196	PHYSIO CONTROL DC Power Adaptor 12V, under-cabinet charger for Lifepak 12, (for operation with NICD batteries), second-hand		395,00 EUR
1003579	PHYSIO CONTROLL Hard paddles for LP 12, second-hand		450,00 EUR
1004126	AD electronic charger for Lifepak 12 with Li-ion accumulator, including charger cable, NEW		1.200,00 EUR
1004659	PHYSIOCONTROL vehicle mount for defibrillator Lifepak 15, used		490,00 EUR
1002708	SCHILLER Defigard 2002 monophasic semi-automatic & manual defibrillator with ECG module, printer incl. charger kit for 12V or stationary supply , second-hand		980,00 EUR
1005039	SCHILLER BTCD, Schnellladegerät für Schiller Akkus z.B. für Schiller Defiguard, inkl. Netzteil, gebraucht		190,00 EUR
1003150	GE Healthcare Responder 3000 monophasic emergency defibrillator with vehicle charging console 12 V, operation with hard paddles and disposable electrodes, measurement of ECG, SPO2 and CO2 (the CO2 sensor is not included in the price, new price *2600), menu and caption in English, BJ 2002 with CE, used		1.490,00 EUR
1003865	PHYSIOCONTROL Lifepak CRplus, fully automatic AED, including new battery and disposable paddles, used		980,00 EUR
1002626	PHYSIO CONTROL Lifepak 500 defibrillator, 3D biphasic AED, German software and menu navigation, year of manufacture 2001, incl. carrying case, new battery		790,00 EUR
1001286	PHYSIO CONTROL Lifepak 500 defibrillator, biphasic AED, incl. carrying bag, with display for battery charge information, English menu, second-hand		590,00 EUR
1001285	PHYSIO CONTROL Lifepak 500 defibrillator, monophasic AED, English Software and menu, incl. case, used		390,00 EUR








1001275	Philips Heartstream, AED, with English menu, used		590,00 EUR
1001278	PHYSIO CONTROL Lifepak 9, AED-modul, Shock advisory adapter, Defi-adapter for adhering electrodes, good condition, second-hand		90,00 EUR
1003011	SCHILLER FRED, AED, ECG measurement, English software, English menu, with rechargeable battery and charger, used		790,00 EUR
1004711	Weinmann Meducore easy, fully automatic AED, delivery including new primary cell and a set of disposable defibrillation electrodes, used		790,00 EUR
1005248	ZOLL M-Series CCT, biphasic defibrillator incl. pacemaker, measurement of ECG, SPO2, NIBP, printer. The battery is charged in the unit by means of a mains cable, used		2.290,00 EUR
1004153	ZOLL M-Series XL: biphasic defibrillator, incl. pacemaker, ECG, SPO2, NIBP, printer, analysis, battery new, only for 12-24 Volt operation, delivery without hardpaddles, second-hand		2.290,00 EUR
1003907	ZOLL M-Series biphasic defibrillator, incl. pacemaker, ECG, printer, analysis, battery new, incl. SPO2 measurement, capnography (in main current mode), protective bag, delivery without hardpaddles, second-hand		3.450,00 EUR
1004818	ZOLL M-Series CCT, biphasic defibrillator, incl. pacemaker, ECG, printer, analysis, battery new, incl. SPO2 measurement, capnography (in main current mode the sensor is included in the price), NIBP, delivery without hardpaddles, used		3.490,00 EUR
1005003	ZOLL M-Series biphasic defibrillator, including pacemaker, ECG, printer, analysis, battery new, with bag. For use with disposable adhesive electrodes, delivery without hardpaddles, used		1.790,00 EUR
1003955	ZOLL protective bag for ZOLL M series defibrillators, good condition, used		120,00 EUR
1004002	ZOLL Base Power Charger, Quick Charge charger for Zoll M series batteries, used		295,00 EUR

1005249	Zoll CO2 main current sensor for use with Zoll M-series defibrillators, second-hand		780,00 EUR
1004820	ZOLL charging holder , for Zoll M-Series, rotatable, for operation on mains 220V, second-hand		450,00 EUR
1004943	ZOLL Intracardiac ladle electrodes for Zoll defibrillators. Diameter 75mm, The electrodes are autoclavable, used		590,00 EUR
1004277	ZOLL 12-channel ceiling transition for E- and M-series, second-hand		85,00 EUR


## EEG

## ECG



1000244	GE Marquette Mac 1200ST, 6-/12 channel ECG unit, recording width 210 mm, monitor screen, keyboard for patient data, used incl. evaluation option, however measurement without ergometry option, used		1.390,00 EUR
1005225	GE MAC 1600, 6-/12 channel ECG, writing width 210mm. Folded layer paper (thermo), colour display, delivery incl. new battery, ECG cable, second-hand.		1.690,00 EUR
1004727	GE MAC 5500, 12 channel ECG unit, with evaluation function, writing width 210mm, incl. 10-wire patient cable with banana plug. For battery and mains operation. Delivery incl. trolley, monochrome display, used		1.890,00 EUR
1004728	GE MAC 5500HD, 12 channel ECG device with evaluation function, writing width 210mm incl. 10 core patient cable with banana plug. Mains and battery operation possible. Delivery incl. device trolley color display, used		1.990,00 EUR
1001292	HP Pagewriter XLE, 6 channel ECG unit with keyboard for input of patient data, battery powered, used		790,00 EUR
1004306	NIHON KOHDEN Cardiofax 6851G portable, battery-operated 1-channel ECG device in carrying case, incl. ECG cable and accessories, second-hand		490,00 EUR











1004327	PHILIPS Pagewriter Trim II, 6-/12 channel ECG device with monochrome display, battery and mains operation, second-hand		850,00 EUR
1003444	Schiller AT-102 6-12 channel ECG with keyboard and screen, battery operated, writing width 20cm, delivery incl. ECG cable used		1.390,00 EUR
1002360	SECA CT 3000, 3/6-channel ECG device, writing width 90mm, battery and mains operation, incl. protective bag and ECG cable, second-hand		690,00 EUR
1005252	SECA CT6i, 6-/12-channel ECG recorder, writing width 210mm, battery and mains operation possible, delivery incl. ECG cable, second-hand		790,00 EUR
1004301	Oesophageal probe small (for animals up to 5kg), NEW		275,00 EUR
1004302	Oesophageal probe medium (for animals up to 15kg), NEW		275,00 EUR
1004303	Oesophageal probe large (for animals from 15kg), NEW		275,00 EUR

## Flexible endoscopes

<b>Bronchoscopes</b>	1004726	OLYMPUS BF-P160: video bronchoscope L 60 cm, diam. 5,5 mm, diam. of working channel 2.0 mm very good condition, second-hand		4.900,00 EUR
----------------------	---------	---	---	--------------

## Bronchoscopes for horses










<b>Choledochoscopes</b>	1002930	OLYMPUS CHF-P10: Paediatric fibre choledochoscope, 0 fractures in the fibres, L: 30 cm, diameter: 5.6 mm, used		1.500,00 EUR
	1004174	OLYMPUS CHF-P20: paediatric fibre choledochoscope, 2 fibre breaks, L: 30 cm, diameter: 5.6 mm, very good condition, second-hand		2.400,00 EUR

<b>Cystoscopes</b>	1003925	OLYMPUS CYF-240: Video cystoscope, diameter 5.3 mm, L 38 cm, diameter of working channel 2.0 mm, used		2.890,00 EUR
<b>Duodenoscopes</b>	1003200	OLYMPUS JF-V2: Video duodenoscope, diameter 11 mm, length: 123.5 cm, very good condition, used Please note: All duodenoscopes have a 90° viewangle		1.300,00 EUR
	1002697	OLYMPUS TJF-M20: Fibre Duodenoscope, Diameter 13mm, Length 124 cm, Side view 90°, very good condition, 0 fibre breakages, second-hand, therapeutic mother scope! Please note: All duodenoscopes have a 90° viewangle		980,00 EUR
<b>Flexible micro endoscopes</b>				
<b>Gastrosopes</b>	1003504	OLYMPUS GIF-Q165 video gastroscope, L 103 cm, diameter 9.5 mm, diameter of working channel 2.8mm, used		5.600,00 EUR
	1004546	OLYMPUS GIF-XP260N: video gastroscope length: 110 cm, diameter 5.0 mm, diameter working channel 2.0 mm, second-hand		5.800,00 EUR
	1002468	OLYMPUS GIF-V2: video gastroscope, diameter 9.8mm, length 105cm, good condition, second-hand		2.900,00 EUR
	1002746	OLYMPUS Gastro-Vet 240K (Video-Gastroscope for small pets), good condition, L: 155 cm, Ø 4,8 mm, - This item is intended solely for veterinary medicine and export - NEW		4.980,00 EUR
	1003937	OLYMPUS Gastro-Vet 240 P: video gastroscope for horses, field of view 140°, depth of field 5 to 100mm, angling 4x diameter 13.2 mm, length 320 cm, working channel diameter: 3.7 mm, suitable for Olympus CV-240 and CV-260 processors, NEW		6.800,00 EUR
	1000407	OLYMPUS Gastro-Vet 240 H (Video-Gastroskop für Kleintiere), L: 135 cm, Ø 9,5 mm, NEUAUFBAU - Dieser Artikel ist ausschliesslich für die Veterinärmedizin und den Export bestimmt - NEU		4.980,00 EUR
<b>Hysteroscopes</b>	1004063	WOLF 7325.122: fibre hysteroscope 2.5mm/7.5 Charr., working length 20 cm, diameter of working channel 1.5 mm, no fibre breaks, second-hand		1.300,00 EUR

## Jejunoscope



## Coloscopes

1002240	OLYMPUS CF-1T 20L, good condition, fully insertable, 8 fibre breaks fibre breaks, L:168 cm, diameter 13 mm, working channel 3,2mm, second-hand		2.800,00 EUR
1003503	OLYMPUS CF-Q165I video colonoscope, L 133, diameter 12.9 mm, diameter of working channel 3.7 mm, used		6.800,00 EUR
1004717	OLYMPUS CF-Q180AI: video colonoscope L 133 cm diam. 12.8 mm, diam. of working channel 3.7 mm, very good condition, second-hand		7.900,00 EUR
1004982	OLYMPUS CF-HQ 190i, flexible videocoloscope, length 133cm, diameter 13,2mm, diameter of working channel 3,7mm, very good condition, used		9.800,00 EUR
1004861	OLYMPUS CF-Q260Di: video colonoscope, L 138 cm, diameter 12.4 mm, diameter working channel 3.2 mm, used		6.900,00 EUR
1005130	OLYMPUS CF-H260DL, diameter 13.2mm, biopsy channel diameter: 3.7mm, working length 168cm, angulation up/down 180°, right/left 160°, second -hand.		6.900,00 EUR
1003508	OLYMPUS PCF-240S video sigmoidoscope, L 100 cm, diameter 13.2 mm, diameter of working channel 3.2 mm, used		1.500,00 EUR
1004141	OLYMPUS PCF-Q260AL: pediatrice video coloscope, diameter 11,3 mm, length 168 cm, working channel diameter 3.2mm, used		6.900,00 EUR
1004862	OLYMPUS OSF-4: Fibre sigmoidoscope L 60 cm, diam. 12.6 mm, good condition, 0 fibre breaks, used		1.450,00 EUR







## Laryngoscopes











## Uretroscopes











## Endoscopy-cold light sources



**Olympus**

1004275	Olympus CLE-E, halogen light source 150 Watt with reserve lamp, suitable for flexible endoscopes of the E-series, second-hand		980,00 EUR
1001968	OLYMPUS CLE-10 cold-light source, 2 x 150 watt halogen, used		780,00 EUR
1004258	OLYMPUS CLE-165: Halogen cold light source 2x150Watt, second-hand		2.200,00 EUR
1000305	OLYMPUS CLK-4 150 Watt halogen, with light insufflation, portable unit, used.		580,00 EUR
1001961	OLYMPUS CLV-10, OES Xenon cold light source 300 W, second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		780,00 EUR
1001963	OLYMPUS CLV-U20, Xenon cold light source 300 W with insufflation, only (OES mode, no OES 100 or OES 200), second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		980,00 EUR
1003083	OLYMPUS CLE 145 EXERA 2 x 150 watt cold-light source, halogen, for endoscopy, good condition, used		1.100,00 EUR
1001962	OLYMPUS CLV 160, xenon cold light source 300W, good condition, new xenon lamp, used The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		2.900,00 EUR
1004163	OLYMPUS CLV-180: Cold light source, xenon 300 Watt, including new bulb, second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		3.900,00 EUR
1003541	OLYMPUS CLV-260: Cold light source 300 Watt Xenon, used The cold light source is equipped with a new xenon bulb. The average lifetime is about 500 hours. The lifetime indicator of the bulb is currently not available. This is only used for approximate determination of the service life.		2.800,00 EUR

1002332	OLYMPUS CLV-S cold light source for rigid endoscopes, 220 V, 4 A, second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		580,00 EUR
1002333	OLYMPUS CLV-S 20 cold light source , 220 V/300 W,xenon, second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		680,00 EUR
1002646	OLYMPUS CLV-S30, xenon cold light source for rigid endoscopy, new xenon lamp, new lamp life meter, used The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		780,00 EUR
1002474	OLYMPUS CLV-S40 Visera: Xenon cold-light source for rigid endoscope, 300 watt, used The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.		980,00 EUR
1004889	OLYMPUS AC-10S: Light source adapter		68,00 EUR
1004890	OLYMPUS AR-SX12E: Camera adapter for rigid endoscopes		120,00 EUR
<b>Other</b>	1002740 ERBE/WOLF: Cold Light Source 150 Watt Halogen (with integrated replacement lamp), the device is able to be combined with ERBE-devices in the T-Series, second-hand		290,00 EUR
	1000295 LEMKE HL752 Endospot, 2 x 75 Watt halogen, second-hand		290,00 EUR
	1001958 LEMKE XL 754-2 Endobeam: Xenon cold light source 35 watt, with 1 replacement lamp, second-hand		390,00 EUR
	1001955 CIRCON MV 9083 automatic video light source Xenon 300W, with Storz connection, second-hand		480,00 EUR

1001954	DYONICS Auto Brite Illuminator II 250 W metal vapour discharge lamp, second-hand		380,00 EUR
1003221	DYONICS Xenon light source 200 watt with change adapter for ACMI/Wolf/Olympus/Storz, used		680,00 EUR
1000299	FUJINON Light Source FIL-150, with facility for insufflation and rinsing, used		320,00 EUR
1003960	GIMMI Alpha Xenon Light 180 W: STORZ connection, lamp life 500 hours, used		980,00 EUR
1002766	LEMKE L100 Arthrolight 75 Watt Halogen Cold Light Source for arthroscopy, second-hand		390,00 EUR
1001976	SCHÖLLY FX.650.2 halogen cold-light source with spare bulb, 150 watt, used		390,00 EUR
<b>Storz</b>	1004594 STORZ 20132620 Xenon cold light source 100 Watt, adjustable brightness, including insufflation pump, for use with flexible Storz endoscopes, second-hand		1.550,00 EUR
	1002145 STORZ 20131520 XENON NOVA 175: Xenon cold-light source 175 watt, colour temperature 5600K, choice of STORZ or OLYMPUS connection, management of the brightness of the light is carried out using a mechanical dimmer, instabilities of the light arc are therefore avoided and a maximum service life of the bulb is ensured, used		1.550,00 EUR
	1004088 STORZ 487-OES: Adapter for connection of flexible endoscopes to STORZ light sources, used		180,00 EUR
	1004887 STORZ 487-U: Adapter for Storz cold light sources. Rotatable adapter for use with light guide cables ACMI, WOLF, STORZ, PILLING, second-hand.		150,00 EUR
	1005094 STORZ 487-UO: Adapter for Storz cold light sources. Rotatable adapter for use with light guide cables ACMI, WOLF, STORZ, Olympus, second-hand.		150,00 EUR








## Wolf

### Endoscopy-rigid endoscopes

### Rigid endoscopes












### Rigid micro endoscopes

### Endoscopy-fittings



<b>Endoscopy towers</b>	1004387	OLYMPUS Emergency Trolley, mobile endoscopy trolley, incl. 3 removable disinfection trays, 2 wide shelves 95cmx50cm, used		380,00 EUR
	1004925	Olympus endoscopy trolley with 3 shelves, bottle holder, holding bar for forceps, brushes etc., second-hand		350,00 EUR
	1005049	OLYMPUS endoscopy trolley with 4 floors, 2 braked castors, holding bar for brushes etc. Dimensions WxHxD 72cmx106cmx70cm, used		350,00 EUR
	1002888	OLYMPUS WM-30 endoscopy trolley with monitor support, 2 adjustable bases 60 x 60 x 150 cm, used		780,00 EUR
	1002338	OLYMPUS WM-60 mobile workstation, second-hand		680,00 EUR
<b>Video</b>	1003823	Olympus OVC-140, attachable camera for flexible endoscopes of Olympus 140 series, used		480,00 EUR
	1002034	OLYMPUS OVC-200 camera adapter to connect Olympus fibre-optic endoscopes to CV-200 processors, good condition, the CV-200 and OVC-200 unit gives a video camera for fibre-optic endoscopes, second-hand		380,00 EUR












1001191	OLYMPUS OTV-S5 video camera with camera head for rigid endoscopy, second-hand		1.180,00 EUR
1005122	OLYMPUS VISERA Pro OTV-S7Pro, video camera for endoscopy, delivery includes camera head Olympus OTV-S7P HD, second-hand		3.900,00 EUR
1004011	OLYMPUS OTV-S6 adapter for rigid endoscopes NEW		190,00 EUR
1004420	OLYMPUS OTV-S6 PAL system, camera head with cable and connection, second-hand		680,00 EUR
1005015	OLYMPUS EMC camera adapter, used		490,00 EUR
1003939	OLYMPUS OTV-SI: video processor with integrated cold light source 150 Watt halogen (MD-151), for connection of video endoscopes: ENF-V2, ENF-VT2, ENF-VQ, or of camera heads: OTV-S7H-N, OTV-S7H-1N. OTV-S7H-1D, used		2.900,00 EUR
1004274	OLYMPUS CV-E, video processor for Olympus e-series flexible endoscopes, second-hand		980,00 EUR
1000365	OLYMPUS CV-140 EVIS video processor, good condition, second-hand		1.650,00 EUR
1003084	OLYMPUS CV-145 video processor for flexible endoscopy, good condition, used		1.980,00 EUR
1004846	OLYMPUS MAJ-845 keyboard for Olympus CV-160 series processors, second-hand.		270,00 EUR
1004257	OLYMPUS CV-165: video processor, second-hand		2.500,00 EUR












1002657	OLYMPUS EVIS EXERA CV-160 Video-Processor, good condition, second-hand		2.550,00 EUR
1004172	OLYMPUS CV-180 Video-Prozessor, gebraucht		3.500,00 EUR
1003542	OLYMPUS CV-260SL: video processor, used		2.500,00 EUR
1002431	OLYMPUS CV-240: Video Processor, good quality, second-hand		1.650,00 EUR
1004843	Olympus MAJ-658 Keyboard for Olympus processors of the CV-240 series, second-hand		220,00 EUR
1004844	OLYMPUS MAJ-844, keyboard for Olympus processors of the CV-145 series, second-hand		220,00 EUR
1004847	OLYMPUS MAJ-1428 keyboard for Olympus processors of the CV-180 series, second-hand		750,00 EUR
1004848	OLYMPUS MAJ-1463, keyboard for Olympus processors of the CV-165 series, second-hand		350,00 EUR
1004849	OLYMPUS MAJ-1536, keyboard for Olympus processors of the CV-260SL series, second-hand		350,00 EUR
1002028	OLYMPUS Keyboard, keyboard for Olympus CV-200 processor, second-hand		120,00 EUR
1002029	OLYMPUS PigTail (for the systems 100, 130, 200, 230) connecting cable from the processor to the video endoscope, second-hand		580,00 EUR












1002030	OLYMPUS PigTail (for the systems 140, 240) connecting cable from the processor to the video endoscope, second-hand		580,00 EUR
1002031	OLYMPUS PigTail (for the system 145, 160) connecting cable from the processor to the video endoscope, second-hand		880,00 EUR
1003807	OLYMPUS MD-680: Pigtail, used		580,00 EUR
1004851	OLYMPUS MH-237, pigtail cable for connecting an endoscope to the processor, second-hand		580,00 EUR
1003809	OLYMPUS MH-976: Pigtail, used		580,00 EUR
1003810	OLYMPUS MAJ-843: Pigtail, used		580,00 EUR
1003814	OLYMPUS MAJ-1154: Pigtail, used		580,00 EUR
1003811	OLYMPUS MAJ-1430: Pigtail, used		2.250,00 EUR
1002663	OLYMPUS MH 199 Character Generator for connection to an Olympus OTV S5 NTSC Version, second-hand		290,00 EUR
1002032	OLYMPUS adapter for universal connection, with this adapter a camera head for a rigid endoscope can be adapted to a flexible endoscope (OES series), new		198,00 EUR
1004165	STORZ Telecam SL II 202 130 20, for rigid endoscopy, very good condition, used		1.790,00 EUR

























1001192	STORZ Telecam SL Pal 202120 20, endoscopy camera for endoscopy, incl. camera head, good condition, used		1.450,00 EUR
1000372	WOLF 5501 CCD ENDOCAM with camera 5501.901, second-hand		980,00 EUR
1004481	STORZ 20042020 Medi-Pack PAL. Camera processor with cold light source and screen. Including camera head 20212030 and input keyboard, second-hand		3.900,00 EUR
1004837	OLYMPUS OEV 191, flat screen for endoscopy, 19 inch LCD display, used		790,00 EUR
1005093	SONY PVM-14N5E, 14 inch video monitor for endoscopy, second-hand		290,00 EUR
1002023	SONY PVM 14L1 MDE video monitor, 500 lines, 14 inches, second-hand		280,00 EUR
1002600	Sony PVM 14N5MDE endoscopy monitor 14 inches, good condition, second-hand		280,00 EUR
1002181	SONY Triniton PVM 1442 QM video colour monitor 14 inches, second-hand		280,00 EUR
1003158	Sony PVM 1443 MD medical monitor with BNC, RGB, Y/C connections, used		280,00 EUR
<b>more fittings</b>			
1001187	WOLF 2232 LAPARO CO2-Pneu: CO2 insufflator with a max. flow rate of 30 litres/min., supplied without high-pressure hose for the CO2 cylinder, second-hand		980,00 EUR
1003704	Olympus UPD: device for positioning and visualization of the colonoscope in the intestine, second-hand This device is compatible with the following endoscopes: CF-240DL/I, CF-Q260DL/I and CF-H260D L/I, second-hand		2.980,00 EUR

1003800	OLYMPUS MH-869: Remote control, used		68,00 EUR
1002080	OLYMPUS LS-10 Teaching, this device is fitted to the ocular of a fibrescope and enables a second person to observe the endoscopic image. Good condition, no fibre ruptures, second-hand		298,00 EUR
1002007	OLYMPUS MB-155 leakage tester for all flexible Olympus endoscopes, to connect to MU-1 or cold light source, second-hand		98,00 EUR
1000077	OLYMPUS MU-1: Maintenance pump for flexible Olympus endoscope, second-hand		330,00 EUR
1004657	Leak tester with adapter for Olympus endoscopes, NEW		169,00 EUR
1005029	STORZ leak tester with Storz adapter for flexible endoscopes, second-hand		88,00 EUR
1004292	OLYMPUS ECR-VET: CO2 insufflator for gastroscopy and colonoscopy, including trolley and CO2 pressure reducer. Flow between 0,2 and 6l/ min, max. pressure 280 mmHG or 0,4 bar, used This item is exclusively designed for being used in the field of veterinary medicine or for export purposes.		980,00 EUR
1004495	OLYMPUS MH-317 foot switch insufflator, used		280,00 EUR
1004412	OLYMPUS MAJ-210: valve for working channel (bronchoscopes etc.), second-hand		9,80 EUR
1001951	OLYMPUS MD-431: rinsing bottle for endoscope series 10, 20, 30, E-Line, 100, 130, 200, 230, V-class, new		259,00 EUR
1004132	OLYMPUS MAJ-901: Rinsing bottle for endoscope from the 140 range, NEW		230,00 EUR

1004693	OLYMPUS MB-399: Footswitch, used		89,00 EUR
1002927	STORZ 263100 20 CO2 Hydromat: price applies for the basic device with no accessories, used		520,00 EUR
1002928	STORZ 263100 20: Complete accessories, NEW, without basic device		520,00 EUR
1004256	AESCLAP Fluid 3: powerful irrigation pump for laparoscopy, pressure-controlled, for reusable and for disposable tubing sets, second-hand		1.480,00 EUR
1002131	Wisap CO2 Aqua-Purator 1611. The device provides continuous flushing of the abdomen and aqua dissection during a diagnostic laparoscopy. The accessory (bottle pin for NACL) is not included in the delivery and must be ordered directly from the manufacturer. Technical data: continuously adjustable pressure up to 750 mmHg. NACL instillation quantity: selectable 1 litre or Hi-Flo 3 l/min, second-hand		480,00 EUR
1001985	WOLF 2221 Hystero-CO2-Pneu: Max. intrauterine pressure 200 mm HG, max. gas flow 99 ml/min, intrauterine working pressure 10-150 mm HG can be continuously selected in advance, second-hand		680,00 EUR
1004074	WOLF 2211 Endo-Irrigator, used		680,00 EUR
1001986	OLYMPUS MD-777 Pulse Unit, second-hand		250,00 EUR
1001999	Biopsy excision forceps, oval spoon without spur, used		98,00 EUR
1002000	Biopsy excision forceps, oval spoon with spur, used		98,00 EUR
1002001	Biopsy forceps, alligator, new		235,00 EUR












1002693	Probe Excision Forceps, tripple gripper, second- hand		194,00 EUR
1004187	PAULDRACH 155-929-01: Disposable biopsy forceps with mandrel (coated), L 230cm, diameter 2.3 mm, / piece, sterile packed		45,00 EUR
1002005	Suction/irrigation probe with syringe connection length 2300mm, diameter 2.3mm, new		48,00 EUR
1003756	OLYMPUS MB-156 ETO Cap, used		15,00 EUR
1003758	MB-107: AW Channel Cleaning Adapter, NEW		25,00 EUR
1005026	OLYMPUS MAJ-420, silicone connection tube, used		24,50 EUR
1005119	OLYMPUS Flushing Adapter Set VET. With this adapter set, the air/water channel of our equine gastroscope Gastrovet 240P can be cleaned. Consisting of a cleaning valve and an adapter hose. This article may only be used in veterinary medicine, second-hand.		69,00 EUR
1003760	Olympus BW-15B Channel Cleaning Brush, used		55,00 EUR
1004994	Olympus BW-7B, Channel Cleaning Brush, reusable, new		69,00 EUR
1005028	OLYMPUS MH-856, cleaning adapter for the suction channel, good condition, used		23,00 EUR
1003572	OLYMPUS MB-196 (also MB-627). Air / Water valve, NEW		59,00 EUR









1003573	OLYMPUS MH-438: Air / water valve for endoscopes OES 40 and higher and EVIS 140 (240), NEW		169,00 EUR
1004411	OLYMPUS MH-081: suction valve for bronchoscopes, second-hand		9,80 EUR
1003574	OLYMPUS MB-197: suction valve, NEW		59,00 EUR
1003575	OLYMPUS MH-443: Suction valve for endoscopes OES-40 and higher and EVIS 140 (240), NEW		139,00 EUR
1004281	OLYMPUS MAJ-207: Suction valve, NEW		12,90 EUR
1003854	OLYMPUS MH-946, NEW		70,00 EUR
1003855	OLYMPUS MH-944, NEW		98,00 EUR
1004143	OLYMPUS MH-948: AW Channel Cleaning Adapter, second-hand		35,00 EUR
1004262	Olympus MD-252: protective cap waterproof, second-hand		98,00 EUR
1002683	Olympus MH-553: Watertight protection cap, second-hand		125,00 EUR
1002786	OLYMPUS Rubber Cap for Work Channel, NEW		15,00 EUR









	1005030	STORZ 13391 SAV Suction valve reusable, good condition used		69,00 EUR
	1005031	STORZ 13991 SSV, air-water valve Silverscope reusable, good condition, used		89,00 EUR
	1005032	STORZ 13991 SS, irrigation tube reusable for endoscopy, new		25,00 EUR
	1005033	STORZ Biopsy valve, new		9,80 EUR
	1005034	KEYSURGICAL Stopcon, disposable valve set sterile for use with Olympus endoscopes: OES 40/Evis 140/145/160/165/180/185/190/240/260/290, new		11,50 EUR
<b>Chassis and monitor carrier</b>	1004609	GCX WMM-0001-01B: monitor holder for anaesthesia equipment, max. load 27 kg, with adapter plate for Philips monitors, second-hand		125,00 EUR
	1004608	GCX WMM-0001-018: monitor holder for anaesthesia equipment, max. load 27 kg (please specify required adapter plate for monitor), second-hand		180,00 EUR
	1004610	BYTEC Monitorhalter für Narkosegeräte, max. Tragkraft 5-9 Kg, für TFT-Monitor und Tastatur, gebraucht		320,00 EUR
	1001252	Instrument supply table (silent nurse) made entirely from stainless steel, mobile, height-adjustable approximately 80-130 cm, stainless steel plate 42.5 x 62.5 cm with margin, 2 rollers adjustable, new		350,00 EUR
	1004542	Mobile equipment trolley with 4 drawers and sloping storage surface (60 x 46 x 34cm), second-hand		240,00 EUR
	1000353	OLYMPUS TC-C1 Clinical Trolley, mobile endoscopy trolley with 4 surfaces, H: 111 x W: 48 x D: 56cm, good condition, second-hand		580,00 EUR

1004599	Sideboard trolley with holders for cups etc., used		245,00 EUR
1001799	DATEX wall support for DATEX Light Monitor with wall rail for continuous height adjustment, second-hand		135,00 EUR
1001934	DRÄGER chassis for Babylog, second-hand		380,00 EUR
1004715	DRÄGER chassis for Babylog 8000, second-hand		580,00 EUR
1001929	DRÄGER chassis for EVITA 2 / EVITA 2 / EVITA XL, very good condition, second-hand		440,00 EUR
1004225	Chassis for Trajan (please specify which model when ordering: Trajan, Trajan 800 or Trajan 808), with two bottle holders, NEW This item is exclusively designed for veterinary medicine and export-.		740,00 EUR
1001933	DRÄGER chassis (original) for Titus, with monitor support and drawer unit with 2 drawers, second-hand		780,00 EUR
1005280	Waste bin with trolley for use in operating theatres and treatment rooms, capacity approx. 15 litres, second-hand		75,00 EUR
1003775	DRÄGER small cabinet for DRÄGER modular system, H28 x W 21 x D 30 cm, second-hand		80,00 EUR
1004860	DRÄGER small cabinet for DRÄGER modular system, H24 x W 21 x D 30 cm, second-hand		80,00 EUR
1004643	DRAEGER carrier / wall mount) for DRAEGER Sulla 808 V and others, used		320,00 EUR













1001936	DRÄGER monitor support for standard rail B43 (44) x D 30 (36)cm, ideal for supporting DRÄGER monitor PM series and/or Capnodig, Oxysat, Iris, etc. second-hand		120,00 EUR
1000510	ERBE chassis for ERBE HF devices, with bottle holder for argon bottle, used.		280,00 EUR
1003316	ERBE device trolley for HF devices e.g. ACC 450, ICC 350, used		180,00 EUR
1004583	ERBE Trolley for ERBE ICC series and for ERBE APC 300, used		280,00 EUR
1003872	GE equipment trolley, suitable for GE MAC 5000 and 1200; 2 braked castors, used		125,00 EUR
1005007	Valleylab chassis for HF surgical device Valleylab Force FX, the middle floor offers space to accommodate a smoke evacuation unit, used		175,00 EUR
1005111	CONMED chassis, for mounting an HF surgical unit from Conmed. Dimensions of platform: 31x31cm, height 84cm, incl. 2 storage baskets, second-hand.		290,00 EUR
1004085	KLS MARTIN maXium-Cart II: equipment trolley, used		490,00 EUR
1003970	Mobile infusion stand on 5-legged chassis, height-adjustable, 2 lockable castors used		30,00 EUR
1001345	Infusion support, stainless steel, mobile, adjustable rollers, stepless height adjustment including dropper, new		69,00 EUR
1004920	DRÄGER 2M86464 Quickstop Hined arm, AS NEW		98,00 EUR











1003902	Hose holder, 3-piece, used		68,00 EUR
1003901	Hose holder, 2-piece, used		68,00 EUR
1005168	DRÄGER holding arm, 3-part, for patient tubes and ECG cables. Suitable for anaesthesia devices such as Julian, Cicero EM, used.		79,00 EUR
1003773	TRUMPF chassis with 2 shelves, used		180,00 EUR
1003774	TRUMPF trolley without drawer, used		180,00 EUR
1000734	DRÄGER Julian mobile anaesthesia workstation, with PM8060 Vitarä (ECG monitor with module box, recording of ECG, IBP, NIBP, temperature, SPO2, good condition, second-hand) comprising bellows-in-bottle ventilator with electrical control, respiratory volume of 50 ml-1400 ml and for pressure-controlled ventilation (PCV). Electronic gas dosing for O2, N2O and air. Safety device: failure warning device for O2,N2O,Air,nitrous oxide block, oxygen-ratio controller, monitoring of CO2 and the volatile anaesthetics Halothane, Enflurane, Isoflurane, Sevoflurane and Desflurane with automated anaesthetic gas identification, second-hand The buyer can return the equipment offered in the lost and found within 14 days after receipt. Any further warranty is excluded. For more information, please see our terms and conditions.		4.900,00 EUR
1000733	DRÄGER Julian mobile anaesthesia workstation comprising bellows-in-bottle ventilator with electrical control, respiratory volume of 50 ml-1400 ml and for pressure-controlled ventilation (PCV). Electronic gas dosing for O2, N2O and air. Safety device: failure warning device for O2,N2O,Air,nitrous oxide block, oxygen-ratio controller, monitoring of CO2 and the anaesthetics Halo.,Enfl.,Isofl.,Sevofl. and Desfl. with automated anaesthetic gas identification, second-hand The buyer can return the equipment offered in the lost and found within 14 days after receipt. Any further warranty is excluded. For more information, please see our terms and conditions.		3.900,00 EUR
1004480	OLYMPUS EU-M30, Endoscopic Ultrasound Center, used The devices offered here in the lost and found section can be returned by the buyer within 14 days after receipt. Any further warranty is excluded. For further information please see our terms and conditions		250,00 EUR

1005131	OLYMPUS BF-260, video bronchoscope, length 60cm, diameter insertion tube 4.9mm, diameter biopsy channel 2.0mm. The device has a defective CCD chip that needs to be replaced. The buyer can return the equipment offered here in the treasure trove within 14 days after receipt. Any further warranty is excluded. You will find further information in our terms and conditions		690,00 EUR
1005167	Olympus GIF-E: Fibre gastroscope: diameter 9.8 mm, 103 cm, 22 fibre breaks, the insertion tube has 2 kinks, the supply tube has 1 kink at the light connector, second-hand. The buyer can return the equipment offered in the lost and found within 14 days after receipt. Any further warranty is excluded. For further information, please see our general terms and conditions		750,00 EUR
1003390	OLYMPUS (?) mini scope, diameter 2.1 mm, L 99 cm, with no eyepiece, 0 fractures in the fibres, good condition, used. The devices offered here in the repository can be returned by the purchaser within 14 days of receipt. Any further warranty is excluded. Please see our general terms and conditions for further information.		890,00 EUR
1002076	OLYMPUS HYF-P, diameter: 3.6mm, length: 20cm, working channel diameter: 1mm, 2 fiber breaks, can be angled in 2 directions. The buyer can return the equipment offered here in the lost and found within 14 days after receipt. Any further warranty is excluded. For more information, please see our terms and conditions.		680,00 EUR
1003383	FRIATEC LQ 8250 IR microprocessor-controlled cold-light source, 250 watt halogen with low level light, used. The devices offered here in the repository can be returned by the purchaser within 14 days of receipt. Any further warranty is excluded. Please see our general terms and conditions for further information.		390,00 EUR
1004989	DRÄGER C700 patient monitor. For integration into an existing monitoring system, second-hand. The devices offered here in the repository can be returned by the purchaser within 14 days of receipt. Any further warranty is excluded. Please see our general terms and conditions for further information		490,00 EUR
1004477	Pentax EPM 3500, processor with Xenon cold light source, second-hand The devices offered here in the lost and found may be returned by the buyer within 14 days of receipt. Any further warranty is excluded. You will find further information in our general terms and conditions		980,00 EUR
1004856	GEUDER G-28252ECV Megatron Cool" III-PV" (version ext. sterilizable pressure transducer) ophthalmosurgical system, microprocessor-controlled unit with I/A system, ultrasound incl. Cool Flash" mode for MICS" (bimanualPako), vitrectomy (magnetic and pneumatic) and diathermy; incl. pressure transducer G-29850 With G-28540V combination foot switch FS4 Both units are tested by the manufacturer. The devices offered here in the repository can be returned by the purchaser within 14 days of receipt. Any further warranty is excluded. Please see our general terms and conditions for further information.		2.900,00 EUR











**Gynaecology**

1004857	GEUDER G-29850 Pressure transducer 3 pcs. GEUDER G-28430 Piezoceramic ultrasonic handle 4 pcs. The devices must be serviced by the manufacturer!		1.980,00 EUR
1004858	GEUDER Box with accessories for MEGATRON G-28252 4 pcs. The devices must be serviced by the manufacturer!		1.980,00 EUR
1004964	STEPHAN AG-Monitor, Multigasmonitor zur Messung von Halothan, Isofluran, Sevofluran und Desfluran. Einstellbare Alarmgrenzen, gebraucht. Die hier in der Fundgrube angebotenen Geräte kann der Käufer innerhalb von 14 Tagen nach Erhalt zurückgeben. Jede weitere Gewährleistung ist ausgeschlossen. Weitere Informationen erhalten Sie in unseren AGB		490,00 EUR
1005045	SMITH & NEPHEW 660HD Image Management Center. Documentation unit for endoscopy. The device is functional, but the front USB input is defective. The devices offered in the treasure trove can be returned by the buyer within 14 days of receipt. Any further guarantee is excluded. Further information can be found in our terms and conditions		680,00 EUR
1005084	ZEISS camera Medilive for mounting on the beam splitter, delivery without processor. The devices offered in the lost and found may be returned by the purchaser within 14 days of receipt. Any further warranty is excluded. For further information, please refer to our general terms and conditions. used		390,00 EUR
1005008	SCHMITZ Medi-Matik series 118, gynaecological examination chair with height-adjustable leg holders and collecting tray, electric height adjustment controllable with foot switch, back plate manually adjustable, the back and seat upholstery are newly upholstered. second-hand.		1.200,00 EUR
1004037	COROMETRICS 120 Cardiotocograph with 1x ultrasound and 1x Toco transducers (new), additionally measurement of SPO2 and blood pressure, second-hand		1.890,00 EUR
1005091	KARL KAPS ViCo S, mobile video colposcope, zoom with autofocus, working distance 250-350mm, freeze function, one-hand operation. Illumination by LED's, video output with BNC connection, electronic green filter, second-hand.		3.590,00 EUR
1005285	LEISEGANG 3DLF, mobile colposcope on mobile stand. Straight binocular view, working distance 400mm. Magnification from approx. 5-20x. Halogen illumination 50Watt. Zoom 7.5 minimum, 30 maximum, second-hand.		1.100,00 EUR
1005217	LEISEGANG 103, mobile colposcope on rolling stand. Straight binocular view, working distance. Magnification of approx. 5-20x. Halogen illumination 20Watt. Zoom 7.5 minimum, 30 maximum, used.		1.100,00 EUR

1002442	LEISEGANG 8 BW: colposcope, binocular, magnification 20 x, green filter, light source 6 V, 15 Watt, for attachment to a gynaecological chair, second-hand		890,00 EUR
1005206	LEISEGANG colposcope for mounting on a gynaecological chair, incl. light source 6V, colour filter, second-hand		490,00 EUR
1003172	KAPS SOM 52, mobile colposcope, working distance 300mm, magnification from 3x to 21x, zoom 0.4minimal, 2.5maximal eyepieces WF 10x, halogen lighting 150Watt, color filter, used		2.300,00 EUR
1003522	GE Corometrics 170 series, cardiotocograph, US and toco probe, used		1.690,00 EUR
1003470	PHILIPS Series 50 A: Cardiotocograph with 1 US and 1 Toco probe, second-hand		1.590,00 EUR
1004581	SONICAID Team duo, cardiotocograph incl. thermal printer and 2x ultrasound probe and 1x Toco probe, second-hand		1.450,00 EUR
1005025	Huntleigh Sonicaid FM800 Encore, Cardiotocograph, with the possibility to operate a 2nd ultrasound transducer for twin imaging. Delivery incl. 1x Toco probe, 1x ultrasound probe, device chassis with 2 shelves, 2 castors lockable, second-hand.		1.875,00 EUR
1004706	RIESTER 5201 Graeve vaginal speculum 75 x 20 mm, NEW		25,00 EUR
1004707	RIESTER 5232 Cusco vaginal speculum 2,85 x 35 cm, NEW		25,00 EUR
1005178	MEDIELLE gynaecological chair, new (3 motors), Medielle is equipped with 3 electric motors to adjust the height, the seat angle as well as the angle of the backrest. The minimum height at 0° is 590 mm, the maximum height is 830 mm. The minimum height at maximum seat angle is 810 mm, the maximum height is 1050 mm. The gynaecological chair is operated with a hand control. The colour of the upholstery is at your choice. Without anatomical leg supports (must be ordered separately), new		4.725,00 EUR

1005179	Medielle gynaecological chair, new, (2 motors) Medielle is equipped with 2 electric motors with which the height, and the angle of the backrest can be adjusted. The minimum height at 0° is 590 mm, the maximum height is 830 mm. The minimum height at maximum seat angle is 810 mm, the maximum height is 1050 mm. The gynaecological chair is operated with a hand control. The colour of the upholstery is at your choice. Without anatomical leg supports (must be ordered separately), new		4.250,00 EUR
1005180	Medielle leg plate B101_40-00-00, used for horizontal examinations of the patient. With this accessory the chair is transformed in to an examination table, new		385,00 EUR
1005181	Anatomical foot rest (Goepel) - pair of 2, new		1.100,00 EUR
1005182	Memory functions (3 memory slots) and 3 predefined position (home, horizontal and trendelenburg) included in hand control, gebraucht		275,00 EUR
1005183	Paper roll holder at the back rest for 600 mm paper roll, new		50,00 EUR
1005185	Medielle foot switch, adjustment of the height of the seat angle and the backrest angle is done via the foot switch, new		250,00 EUR
1005186	Medielle comfort footrest (1 pair), stainless steel, new		770,00 EUR
1005187	Colposcope holder for systems such as Olympus, Kaps, Zeiss or Leisegang, can be adapted to the left or right, new		275,00 EUR
1005188	Medielle infusion stand for attachment to the multifunctional rails, new		220,00 EUR
1005189	Medielle Urological drainage, new		275,00 EUR























	1005190	Medielle multifunctional support rail (1 pair) for attachment to the seat, new		85,00 EUR
	1005191	Medielle multifunctional support rail (1 pair) for mounting on the backrest, new		95,00 EUR
	1005192	Medielle retractable castors, operated via the hand control. By pressing the two buttons, the castors are extended and retracted, new		1.250,00 EUR
	1005193	Castors, the chair is on the castors - not retractable. The height increases by 25 mm. Four swivel castors with brake, new		650,00 EUR
<b>HF surgery</b>	1005223	ERBE Vio 300 S HF surgical unit, maximum power in CUT mode: 300 watts at 500 ohms, maximum power in COAG mode: 200 watts safety system: NESSY, frequency: 350 kHz, delivery includes new double pedal footswitch, second-hand .		5.500,00 EUR
	1005086	ERBE APC 2, argon plasma coagulator for use on ERBE Vio units. The control takes place on the ERBE Vio module, used		2.790,00 EUR
	1004905	ERBE ERBEJET 2, basic module for waterjet surgery (hydrosurgery). Particularly gentle on nerves and blood vessels. The unit is supplied without accessories (these are only available from ERBE), secondhand.		8.900,00 EUR
	1005220	ERBE IES 2, smoke evacuator for HF surgery, including new main filter, for use with HF units of the VIO series, second-hand		1.580,00 EUR
	1002177	ERBE ICC 350, HF surgical unit, monopolar cutting 300W, bipolar 120W, dimensions 410 x 152 x 368 mm, good condition, used ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS).		2.900,00 EUR
	1001214	ERBE ICC 200 with EndoCut function, monopolar 200 Watt, bipolar 120 Watt. Very good condition, used ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS).		2.980,00 EUR








1002125	ERBE ICC 200 VET with EndoCut function, monopolar 200 watts, bipolar 120 watts, good condition, used ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS). This item is exclusively intended for veterinary medicine and export-.		2.750,00 EUR
1005171	ERBE ICC 200 with EndoCut function, German version, monopolar 200 watts, bipolar 120 watts. Very good condition, used ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS).		2.980,00 EUR
1000505	ERBE ICC 200 VET without EndoCut function (button present on operating panel, but function not connected). Monopolar 200 watt, bipolar 120 watt. Very good condition, used. ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS). - This item is exclusively designed for being used in the field of veterinary medicine or for export purposes -		1.980,00 EUR
1000496	ERBE ICC 300: monopolar cutting: 300 W, coagulation: 300 W, bipolar: 120 W, dimensions: 410x152x368 mm, weight: 10 kg, very good condition, used. ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS).		2.700,00 EUR
1002426	STORZ 26021 Bipolar Coagulator 50 Watt with foot switch, second-hand. HF-SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (, BIPOLAR FORCEPS, CABLE). This article is exclusively intended for veterinary medicine and export-.		890,00 EUR
1005175	ERBE Erbotom T400 VET, mono- and bipolar, 400 Watt, second-hand. ATTENTION: This unit is intended exclusively for use in veterinary medicine. The neutral plate socket has been replaced by a non-original plug socket. Therefore, the unit must not be operated in human medicine. The function and safety are fully guaranteed. ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS).		1.290,00 EUR
1000499	ERBE Erbotom T175VET, monopolar and bipolar, HF power: monopolar cutting 175 watts, monopolar coagulation 100 watts, bipolar 50 watts, both monopolar and bipolar functions can be activated without foot switch, branded ATTENTION: This device is intended exclusively for use in veterinary medicine. ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODES, 3 DIFFERENT INSERTS)		1.100,00 EUR






1002689	ERBE APC 300, Argon Beamer for argon plasma coagulation. Only for use in combination with ERBOTOM ACC 450, 451, ICC 200 Argon Coag, 300, 350. For argon plasma coagulation (coagulation method in which argon gas which has been ionized by HF current is passed onto the tissue to be coagulated), without pressure reducer, used.		1.590,00 EUR
1003780	ERBE argon pressure reducer, used		395,00 EUR
1002252	ERBE IES 300, smoke gas suction for HF surgery, including a new main filter, for use with all current HF surgical devices, second-hand		980,00 EUR
1004045	KLS MARTIN Maximum beamer, Argon beamer for MARTIN MAXIMUM HF unit (without pressure reducer and accessories), used		1.450,00 EUR
1003148	Olympus PSD-10 VET. Monopolar HF device, the handle can be activated using the foot switch, CUT 80 W, blend 55 W, coag 40 W, used. THIS HF SURGICAL DEVICE IS SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR ELECTRODE, 3 DIFFERENT INSERTS and FOOTSWITCH) This item is exclusively designed for being used in the field of veterinary medicine or for export purposes		780,00 EUR
1001990	OLYMPUS PSD-20 HF surgical device, for endoscopic work, second-hand. delivery without accessories		980,00 EUR
1003012	Valleylab Force EZ, processor-controlled HF device, monopolar cut 300 W, monopolar coag. 120 W, bipolar 70 W, used ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS).		2.590,00 EUR
1002705	Valleylab Force FX, processor-controlled HF unit monopolar Cut 300W, monopolar COag. 120W, bipolar 70W, used, ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS).		3.350,00 EUR











1000516	BERCHTOLD Elektrotom 621: HF surgical device with the following functions: Monopolar type of cutting current with low sloughing (max 200 W) and moderate sloughing (max 150 W), monopolar haemostatic type of cutting current (polypectomy, papillotomy) (max 100 W), monopolar contact-coagulation (max 120 W), forced coagulation (max 150 W) and spray coagulation (max 90 W). Bipolar coagulation (max 95 W). Digital performance display, service display, memory (last setting saved at shut-down), used ALL HF SURGICAL DEVICES ARE DELIVERED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BI-POLAR ELECTRODE, 3 DIFFERENT INSERTS)		2.390,00 EUR
1004494	OLYMPUS MH-551: Fußschalter, gebraucht		320,00 EUR
1002298	Monopolar electrode handle for VALLEYLAB / ERBE, 2 keys, shaft 4 mm, 89214040, new		108,00 EUR
1004968	OLYMPUS ESG-100, endoscopic HF surgical unit, particularly easy and intuitive operation. Cutting power max. 120Watt, coagulation min. 50Watt, max. 120Watt, delivery incl. double pedal foot switch, second-hand.		2.450,00 EUR
1001317	Monopolar electrode handle with dual keys suitable for ERBE T series including 4.5 m connection cable, new 89213041		108,00 EUR
1002397	Monopolar Electrode Handle with dual keys suitable for ERBE ACC/ICC incl. 4.5 m connection cable, new 89213040		108,00 EUR
1001318	Monopolar electrode handle for Martin / Berchtold devices, shaft 4 mm, 2 keys, 89216040, new		108,00 EUR
1001319	Bipolar connection cable for MARTIN / BERCHTOLD, length: 300cm, (80100019), new		74,00 EUR
1001322	Bipolar connection cable for Erbe T series /ACC/ICC, 80100017, new		74,00 EUR
1003668	Bipolar connection cable, 5.00m for VALLEYLAB / ERBE int 80100013, new		81,00 EUR

1005213	Neutral plate VET made of rubber for connection to HF units from Valleylab (Force2, Force EZ, Force FX). The neutral plate has a 2-piece connection cable. Dimensions: 26x16cm, new. As this article does not comply with the regulations for human medicine, we would like to point out that it may only be used in veterinary medicine.		215,00 EUR
1005214	Neutral plate VET made of rubber for connection to Erbe and Martin HF units (Erbotom T series, Erbe ACC and ICC series, Martin Elektrotom series). The neutral plate has a one-piece connection cable. Dimensions: 8x15cm, new. As this item does not comply with the regulations of human medicine, we would like to point out that this item may only be used in veterinary medicine.		175,00 EUR
1005215	Neutral plate VET made of rubber for connection to Erbe and Martin HF units (Erbotom T series, Erbe ACC and ICC series, Martin Elektrotom series). The neutral plate has a one-piece connection cable. Dimensions: 26x16cm, new. As this item does not comply with the regulations of human medicine, we would like to point out that this item may only be used in veterinary medicine.		189,00 EUR
1004522	Neutral electrode cable for ERBE and MARTIN for use with disposable adhesive electrodes, NEW 89101 050		74,00 EUR
1001325	Spherical monopolar electrode, 4 mm shaft, 4 mm diameter, new		9,00 EUR
1001326	Monopolar needle electrode, 4 mm shaft, 4 mm diameter, new		9,00 EUR
1001327	Monopolar wire sling electrode, 10 mm diameter, 4 mm shaft, new		9,00 EUR
1001328	Monopolar knife electrode, 25 mm, 4 mm shaft, new		9,00 EUR
1001329	Bipolar tweezers, 180 mm, tip 1,0 mm for Erbe and Martin, 81518201 new		172,00 EUR
1004280	Bipolar forceps, 150 mm, tip 0.6 mm, 81515201, new		172,00 EUR












1003525	HF loop, oval, length 3600mm, diameter sling 30mm, tube 2,3mm, without handle, new		129,00 EUR
1004324	HF sling with removable handle. When ordering, please indicate the diameter of your working channel and the desired length, NEW		225,00 EUR
1003526	HF loop handle, NEW		135,00 EUR
1003527	HF cable, NEW (connection cable for loop handle)		65,00 EUR
1005216	HF-Test Mate, autoclavable test unit for HF accessories. This unit can be used intraoperatively for testing monopolar and bipolar HF accessories. new		135,00 EUR

## ENT











<b>Infusors perfusors</b>	1002518	BRAUN Infusomat fms infusion pump, for mains and battery operation. 0.1 ml - 999.9 ml in 0.1ml-levels, very good condition, with English software, new battery, second-hand The device already has the current BRAUN Spaceline infusion set.		690,00 EUR
	1002848	BRAUN Infusomat FMS infusion pump, mains and battery operated. 0.1 ml-999.9 ml in 0.1 ml increments, very good condition, comes with German software, new battery, used. The device already has the current BRAUN Spaceline infusion set.		790,00 EUR
	1002190	BRAUN Infusomat P infusion pump, delivery volume 0.1 ml - 999.99 ml/h, new battery, good condition, second-hand.		650,00 EUR
	1001215	BRAUN Infusomat fm infusion pump, mains and battery operated, 0.1-999.9 ml/h in 0.1 ml increments, new battery, used. The unit already has the current BRAUN infusion set.		490,00 EUR
	1003307	B. BRAUN FM controller with FM computer and option to include 3 Perfusors FM or Compact, used		590,00 EUR

1002986	BRAUN fm workplace, mobile, with 3 intakes Infusomat fm, fmS and P, and 6 intakes for Perfusor fm and compact, with 3 infusion holders and data box with connecting cable, RS 232 interface, used		980,00 EUR
1001335	BRAUN Perfusor Compact, very good condition, delivery rate 0.1-99.9 ml/h, bolus rate 800 ml/h, new battery, second-hand This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please! For retroactive changement (after sale) we have to charge EUR 55,-		690,00 EUR
1004688	BRAUN Perfusor Compact Plus, syringe pump with intuitive menu navigation, internal lithium-ion battery, continuous delivery rate 0.01-999.9 ml / hour, second hand		1.190,00 EUR
1005041	B.Braun 871167 KFZ-Ladekabel für Braun Perfusor compact, gebraucht		70,00 EUR
1001337	BRAUN Perfusor fm, very good condition, delivery rate 0.1 to 200 ml/hour, with new battery and integrated mains component, second-hand This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please! For retroactive changement (after sale) we have to charge EUR 55,-		490,00 EUR
1001338	BRAUN Perfusor segura ft, very good condition, delivery rate 0.1 to 99.9 ml/hour, new battery, second-hand This perfusor is ready for use with B. BRAUN perfusor syringes.		390,00 EUR
1003137	BRAUN Perfusor Space, mains/battery operated perfusor (new battery), incl. charging cable, delivery rate for 50 ml syringes: 0.1 - 200 ml/h, battery life approx. 8 hours. Also suitable for use in space stations, second-hand B. Braun perfusor syringes are approved as standard for this perfusor. If you would like to have other manufacturers syringes approved, please let us know when ordering. For a subsequent release of other syringe types we charge EUR 55,00 flat rate		950,00 EUR
1002838	BRAUN Infusomat Space vol. infusion pump 0.1 - 1200 ml/h, NIMH battery new (operation 4 H at 100 ml) incl. power supply unit, second-hand		790,00 EUR
1005226	Braun 8713142 Space Station with integrated SpaceCom, holder for up to 4 space pumps, second -hand		680,00 EUR
1005235	BRAUN 8713140 Space Station, extension module without cover and handle, second-hand		450,00 EUR





















1005238	BRAUN 8713140 Space Station with comfort lid, holder for up to 4 space pumps, second-hand		650,00 EUR
1005172	B.Braun 38245302 Stand clamp for Braun space station for attachment to an infusion pole or a standard rail, second-hand		69,00 EUR
1004927	BRAUN drop sensor space. This drop sensor is exclusively for use with Braun Infusomat space devices. The drop sensor is not absolutely necessary for operation with this model series.		98,00 EUR
1005237	BRAUN 8711747 Universal adapter for drop sensor Infusomat FMS and Infusomat P, as good as new		25,00 EUR
1005158	BRAUN 8713120D, Plug-in power supply for infusomats and perfusors of the Braun space series, new and original packaging.		139,00 EUR
1005253	B.Braun drop sensor 3450578A, suitable for infusion pumps Braun Infusomat FMS and FM, new		154,00 EUR
1004972	B.Braun plug-in power supply UK version for Perfusor / Infusomat space, used		65,00 EUR
1001346	BRAUN holding clamp for BRAUN Infusomat or perfusor (for attachment to the infusion supports, etc.) for segura, segura 7, f and ft models, second-hand		25,00 EUR
1003209	BRAUN holding clamps rotatable for BRAUN Infusomat FM /FMS or Perfusor FM and compact, used		55,00 EUR
1003841	BRAUN holding clamp for BRAUN Infusomat FM, FMS or Perfusor FM and compact, used		40,00 EUR
1002906	BRAUN holding clamps to secure Perfusor or Infusomat Space to a normal splint or an infusion stand, used		98,00 EUR



1004783	BRAUN 8701644 Short infusion pole for Infusomat fm and fms, NEW		45,00 EUR
1004462	BRAUN Spaceline Infusomat line, NEW		2,40 EUR
1004986	FRESENIUS Pilote Anesthesia 2, syringe pump for anesthesia. ATTENTION: the device has a French menu guide and a French keyboard, used.		490,00 EUR
1004987	FRESENIUS Pilote A2, syringe pump for anesthesia. ATTENTION: the device has a French menu guide and a French keyboard, used.		490,00 EUR
1003805	FRESENIUS Volumat MC Agilia, volumetric infusion pump with battery and mains operation, 1-999ml / h, battery new, used		690,00 EUR
1003631	FRESENIUS Injectomat MC Agilia. injection pump with battery and mains operation, including mounting clamp for infusion stand. Approx. 50 possible syringe types. Flow rate 0.1 to 99.9ml/h in 0.1ml steps. Bolus rate 50-1200ml, used		690,00 EUR
1005251	FRESENIUS Injectomat Agilia. Syringe pump for battery and mains operation, incl. retaining claw for infusion stand. Approx. 50 usable syringe types. Delivery rate 0.1 to 99.9ml/h in 0.1ml steps. Bolus rate 50-1200ml, second hand		640,00 EUR
1005117	CAREFUSION Alaris plus GH, syringe pump, compatible with many standard disposable syringes 50ml (Braun Omnifix, BD, Codan, Fresenius, Terumo). Max. Infusion rate with 50ml. Syringes 0.1ml/h-1200ml/h. Max. Bolus rate 10ml/h-1200ml/h. System accuracy +/-2%. Rechargeable NIMH battery to operate the device for approx. 4 hours, used		690,00 EUR
1005170	BAXTER EVO IQ LVP VET, volumetric infusion pump for mains or battery operation. Rate 0.1-2000ml/h, bolus rate 0.1-2000ml/h. We recommend the use of B.Braun Intrafix Safeset infusion lines, used. This article is exclusively intended for use in veterinary medicine or export.		590,00 EUR
1004442	STIHLER Astotherm plus, infusion warmer, for attachment to the infusion stand, second-hand		380,00 EUR












## Incubators









1002322	AIR SHIELDS RW81-1: Resuscitaire, Radiant Warmer, mobile thermal bed with integrated reanimation units, heater module (3 operating modes: pre-heating, skin temperature control, manual mode), microprocessor-controlled heating lamp (180° rotation), integrated Apgar timer, hinged side walls, trolley with integrated drawer and bronchial aspiration, used		1.450,00 EUR
1004833	DRÄGERMU10841 Breathing hose (for Autobreath) disposable, NEW Accessory for DRÄGER/AirShields Isolette 2000		85,00 EUR
1004834	DRÄGER MU11225 Skin temperature Sensor Premi-Probe 3, reusable		95,00 EUR
1004835	DRÄGER MU11228 Skin temp.-sensor3, Dual, Premi-Probe 2, disposable		4,50 EUR
1004000	DRÄGER Babytherm 8000, mobile reanimation unit for neonates, heat radiator 8010 with illumination, electric height adjustment, used		1.980,00 EUR
1003411	DRÄGER Caleo, mobile baby incubator, warming treatment regulated by air temperature control or skin temperature control, air humidification, O2 enrichment, pivotable lying surface, incl. monitoring of air temperature, skin temperature, relative air humidity, O2 concentration, used		4.890,00 EUR
1004113	DRÄGER/AirShields Isolette C2000, mobile incubator for newborns, temperature controlled, monitoring of humidity and O2 content, second-hand		2.980,00 EUR
1004996	ARDO Amelie star, mobile incubator for newborns, monitoring of temperature and humidity in the incubator. Electric height adjustment, inclinable lying surface. Including 2 drawers, braked castors, used.		2.750,00 EUR
1002974	DRÄGER 5400 transport incubator on transport cradle, for external and internal O2 and air supply, bottle connections for O2 + air, power supply via mains connection or car battery, with Babylog 2000, dimensions without cradle height: 86 cm x width: 150cm x depth: 56cm, incl. bronchial suction (O2 injector), dimensions of cradle H 86 x W 189 x D 61cm, very good condition, used		3.200,00 EUR
1002485	DRÄGER Babytherm 8000, mobile thermal bed without acrylic glass cover, second-hand		1.250,00 EUR

	1001357	DRÄGER Photo therapy 800, with chassis, for the light irradiation of premature or newly born infants in the incubator or in the infrared-heated thermal bed to reduce the serum-bilirubin concentration (phototherapy), frequency: 50 Hz, output: 230W, current: 2.2 A, lamp housing 54 x 25.5 x 20.5cm, second-hand		590,00 EUR
	1002237	DRÄGER Photo Therapy 4000, 120W, 4 tubes and separately switchable neon lighting, without chassis, used		790,00 EUR
	1004150	DRÄGER PTS 68-1C, phototherapy lamp on chassis, second-hand		480,00 EUR
	1003588	VICKERS Air-Shields Fluoro Light phototherapy unit with new tubes, used		480,00 EUR
	1005267	EICKEMEYER laboratory centrifuge, rotor for 8 sample tubes 15ml, infinitely variable up to 4500rpm, timer 0-60min. used. This article may only be used in veterinary medicine !!!		220,00 EUR
<b>Laboratory equipment</b>				
<b>Laser surgery equipment</b>				
<b>Medical antiques</b>				
<b>Lights</b>	1001368	HERAEUS mobile examination light with battery function, 4 x 25 watt, battery and mains operated, rollers with brakes, height adjustable side arm, body of lamp can be swivelled, second-hand - This item is intended solely for veterinary medicine or export -		590,00 EUR
<b>Furniture</b>				
<b>Anaesthesia equipment</b>	1004909	DRÄGER Vapor: antiquarian halothane evaporator, on a carrying frame, good condition, device must not be used, only for decorative purposes, used		125,00 EUR
	1004601	DRÄGER Romulus: fahrbares Narkosegerät mit DRÄGER Halothan Verdunster mit integrierten Schubladen und Flaschenhalterung für zwei Flaschen, mit integrierten Druckminderern und Druckmessern für O2 und N2O, mit O2 Flush und Bronchial-Absaugung, gebraucht Medizinische Antiquität! Darf nicht zu medizinischen Zwecken eingesetzt werden.		1.250,00 EUR










**Other**

1004547	DRÄGER circular part for anaesthesia equipment, second-hand Medical antique! Must not be used for medical purposes.		149,00 EUR
1004548	DRÄGER Volumeter for circular part, used Medical antique! Must not be used for medical purposes.		98,00 EUR
1004528	Tube adapter set in wooden stand, used Medical antique! Must not be used for medical purposes.		23,00 EUR
1004529	Test tube stand 2-tier, metal, used Medical antique! May not be used for medical purposes.		15,00 EUR
1004531	MGF spreading speculum size 1, used Medical antique! Must not be used for medical purposes.		18,00 EUR
1004532	GSZ speculum, used Medical antique! May not be used for medical purposes.		18,00 EUR
1004533	Orthopedic Hammer Stainless Germany, used Medical antique! May not be used for medical purposes.		25,00 EUR
1004534	Reflex hammer, used Medical antique! May not be used for medical purposes.		9,00 EUR
1004536	Retractor, used Medical antique! Must not be used for medical purposes.		9,00 EUR
1004535	REFILA glass syringe 5 ml, used Medical antique! Must not be used for medical purposes.		12,00 EUR
1004537	REFILA glass syringe 10 ml, used Medical antique! Must not be used for medical purposes.		12,00 EUR

1005140	Castration forceps antiquarian. This article may not be used in human medicine, second-hand.		75,00 EUR
1005141	Enema syringe antiquarian, approx. 50ml. This syringe must not be used in human medicine, second-hand.		35,00 EUR
1005142	HAUPTNER tuberculin syringes antiquarian, 3 pieces. These articles are antiquarian and must not be used in human medicine, second-hand.		25,00 EUR
1004541	Autoclave container made of aluminium, diameter 16 cm, height 15cm, used Medical antique! May not be used for medical purposes		18,00 EUR
1004550	Urine sampler according to Dr. Vogel, used Medical antique! May not be used for medical purposes.		25,00 EUR
1005143	HAUPTNER instrument or ampoule containers, 3 pieces (15x8cm, 13x8cm, 18x8cm), chrome-plated metal. These containers must not be used in human medicine, used		25,00 EUR
1005144	HAUPTNER Shearing machine hand-operated, antiquarian. This item is exclusively intended for veterinary medicine, second-hand		29,00 EUR
1004606	BROWN glass bottle with cap Acid.acetic. 99-100%, used Medical antique! May not be used for medical purposes.		12,00 EUR
1004658	Pharmacist's scale max. 1kg, dismountable with wooden container, second-hand		168,00 EUR
1005145	HAUPTNER antiquarian wooden instrument case. The instrument set is used for castrating chicken birds, incl. original instructions, second-hand. These instruments may not be used in human medicine.		79,00 EUR
1004252	ZEISS foot switch for OPMI, second-hand		780,00 EUR












1004492	Surcharge for conversion of our surgical microscope from halogen to LED illumination: High-contrast imaging by LED, stepless brightness adjustment, light color 5000 to 5600 Kelvin. As the conversion is a modification of the design, this article is only approved for technical and veterinary purposes		980,00 EUR
1004729	ZEISS OPMI, mobile operating microscope, eyepieces 10x, binocular view 45°, zoom min. 0,5- 2,5 max. electrically adjustable by hand switch. Working distance 150mm magnification 6,7-33x, the device is especially suitable for technical applications because it has no own light source, used		1.500,00 EUR
1004616	ZEISS OPMI 1: Binocular (10X / 22B) surgical microscope, f = 170, straight ahead optic, zoom 0.4 to 2.5 times (tap changer), working distance 200 mm, field of view in the eyepiece: 60 to 10 mm, halogen lighting 100 watts, with reserve lamp and color filter, brightness infinitely controllable, used		4.500,00 EUR
1002410	ZEISS OPMI 99: Mobile surgical microscope, binocular (19 x) tube f = 80, 19 x, objective f = 250, magnification factors: 0.4 - 0.6 - 1.0 - 1.6 - 2.5, halogen lighting 100 W, used		3.600,00 EUR
1003045	Zeiss OPMI 1 mobile operating microscope Eyepieces 12.5x, magnification 0.4/0.6/1/1.6/2.5, illumination 30 W, working distance 200mm, magnification 6.4x-42x. second-hand This item is exclusively designed for veterinary medicine and export-		1.980,00 EUR
1004630	ZEISS OPMI 6-CFC LED: Binocular (10x22) Surgical Microscope Straight view, zoom and focus electronically controlled via footswitch, zoom 0.4 to 2.4, focal length tube f = 170 working distance 200 mm, field of view in eyepiece 10 to 55 mm, Magnification 3.4 to 21.3-fold; Lighting: LED; ATTENTION: It is NOT an ORIGINAL ZEISS light source. Therefore this device is not approved for human medicine.; used		7.800,00 EUR
1005098	ZEISS OPMI, mobile operating microscope on universal stand S2, eyepieces 12.5x, binocular view straight, zoom minimum 6, maximum 40, focal length tube 160mm, working distance 200mm, field of view in eyepiece minimum 50mm, maximum 8mm, magnification 6-40x. Halogen light source 100W infinitely variable, incl. colour filter. second-hand. ATTENTION: not all parts on this device are original ZEISS products - therefore the device may only be used in veterinary medicine or in the technical field.		3.980,00 EUR
1005205	ZEISS OPMI Tec, mobile microscope for technical applications such as fine deburring. This unit must not be used in the medical field, as it has retrofitted LED illumination. Eyepieces 12.5x, straight binocular tube, zoom min. 6, max. 40. Focal length tube 160mm, working distance 200mm, field of view in eyepiece 40mm minimum zoom, 6mm maximum zoom. Magnification 6.4 to 42.5. The unit is supplied with an external power supply to power the LED light source, second-hand.		1.980,00 EUR





1005218	ZEISS OPMI 1F TEC, mobile operating microscope for veterinary medicine or technical use. Eyepieces 12.5x, binocular view, straight, zoom 0.4minimum, 2.5maximum adjustable with step switch, working distance 200mm, focal length tube 125mm, illumination 100W halogen, second-hand.		1.980,00 EUR
1001416	Zeiss OPMI 6-CFC: Binocular (12.5 x) surgical microscope, magnification factor 0.4-2.4 steplessly adjustable, ZOOM and focus controllable by hand switch (can also be operated with nose or mouth), double beam splitter for connecting two observation tubes, Objective f = 200, working distance 400mm, magnification 5-31 times, field of view in the eyepiece 80mm to 15mm maximum. Lighting 100 watt halogen, continuously adjustable, used		5.900,00 EUR
1003692	SEILER SSL-202 ENT microscope for wall mounting, ocular 12.5x, zoom with step switch, magnification approx. 20-40x, 150W, halogen illumination used		1.800,00 EUR
1004487	WILD M650: mobile operating microscope binocular (10 x 21) straight view, zoom 6 to 40 (step switch), magnification 3 to 20x, field of view 130 to 20 mm, focal length 200, working distance 400 mm, illumination 2 x 50 Watt halogen, incl. monocular co-observer tube, second-hand		4.900,00 EUR
1004615	ZEISS 1270-341: Touchscreen 15 Zoll auf ZEISS Fahrgestell CZ-MNR308202-920, gebraucht		980,00 EUR
1004932	ZEISS 304995 Handset for ZEISS OPMI, second-hand		480,00 EUR
1004933	ZEISS foot switch for OPMI, second-hand		600,00 EUR
1005082	Swivel castor for ZEISS S8 column, new		115,00 EUR
1005095	Heavy duty castor for Zeiss tripod S8. Diameter 200mm with thermoplastic polyurethane tread, axle 20mm, new. Please note: this product is NOT an original spare part from Zeiss. The User is responsible for observing the legally valid standards. We recommend the use in the veterinary or technical field.		329,00 EUR
1002681	ZEISS Illuminant 6V/30W with adjustment ring, NEW		17,90 EUR





1001423	Zeiss Photo-/video adapter, second-hand		980,00 EUR
1001425	ZEISS beam separator , double for OPMI, second -hand		480,00 EUR
1005079	ZEISS bridge, for use Wirth a second binocular tube, second-hand		390,00 EUR
1005083	ZEISS XY coupling electric, for use on ZEISS tripod columns S2, S3, S4, second-hand		780,00 EUR
1003353	ZEISS double beam splitter with video adapter and video camera, used		1.980,00 EUR
1003608	ZEISS beam splitter, double, ZEISS photo adapter with mechanical aperture setting, USB HD web cam C310, modified, software with facility for video or still image recording, video up to HD quality, photos up to 5 megapixels, automatic or manual lighting control via software, used		1.980,00 EUR
1003361	ZEISS video adapter (f60/340) for connection to a beam splitter, with Sony BSC 20Q video CCD camera, used		1.980,00 EUR
1001428	ZEISS monocular (12.5 x) observation tube, second-hand		195,00 EUR
1002864	ZEISS binocular (12.5 x) observation tube f = 125/16, second-hand		490,00 EUR
1004950	ZEISS Adapter für Binokulartubus. Mit diesem Adapter kann ein weiterer Tubus z.B für den Assistenten an einen Strahlteiler angeschlossen werden, gebraucht		395,00 EUR
1003351	ZEISS pivotable (0 to 50°) binocular tube, f=170, movable eyepieces (12.5 x ), used		1.870,00 EUR


1001429	ZEISS binocular tube f=125, ocular 12.5 x, second-hand		680,00 EUR
1002278	ZEISS binocular tube (f 170 T*) for threaded ocular, second-hand		650,00 EUR
1004953	ZEISS binocular (12.5x) observation tube f=125/16, second-hand		390,00 EUR
1002266	ZEISS threaded ocular 12.5x/ 18B T*, second-hand		560,00 EUR
1003301	LEICA video equipment suitable for LEICA/WILD surgical microscope, consisting of double beam splitter, photo adapter, digital video camera with mains adapter, used		1.690,00 EUR
1002379	WILD Heerbrugg 384000 Monocular Observer Tube with Ocular 10/22, second hand		490,00 EUR
1002770	WILD HEERBRUGG Binocular Tube with ocular, second-hand		890,00 EUR
1005155	LEICA/WILD holding arm for surgical microscopes, second-hand		180,00 EUR
1005161	LEICA/WILD 182526 Optics suspension for mobile operating microscopes, second-hand		290,00 EUR
1002366	LEICA WILD 319449 Beam Separator with adjustable diaphragm, second hand		850,00 EUR
1005100	LEICA binocular tube for Leica laboratory microscope. Delivery without eyepieces, used		320,00 EUR


1003822 LED illumination for all of the surgical microscopes which we offer NB: Installation of this illumination is only intended for veterinary and technical purposes. No type approval is available for use in human medicine.  450,00 EUR


1004898 Sterilisable handle for surgical microscope WILD M 690 and others, 1 pair, new  78,00 EUR


1005118 HALOGEN LIGHT FITTLE for Leica and Wild surgical microscopes, 12V/50Watt, new  45,00 EUR

1005080 ZEISS OPMI 9 optics, for use on an existing Zeiss tripod or for technical use. Illumination 6V/30W, power consumption via 2 4mm banana plugs, eyepieces 12.5x, working distance 200mm, magnification approx. 10x, incl. colour filter. (Delivery without power supply unit, which can be purchased optionally), second hand.  690,00 EUR


1005147 ZEISS OPMI 6-DF optics for use on an existing Zeiss column or for technical use. Zoom electrically adjustable 0.5-2.5, focus electrically, working distance 200mm, field of view in eyepiece 8mm-40mm. A control unit is required to operate the optics, this is not included in the scope of delivery, second-hand.  850,00 EUR

1005146 ZEISS Colposcope Plus optics, for use on an existing Zeiss column or for technical use. Zoom 0.6/1.0/1.6, eyepieces 12.5x, field of view 80mm-30mm, working distance 400mm, connection for cold light cable, beam splitter for e.g. camera adapter, second-hand.  890,00 EUR

1005148 ZEISS OPMI 9-F optics, for mounting on an existing Zeiss column or for technical use. Eyepieces 10x, working distance 225mm, field of view in eyepiece 18mm, coaxial illumination 6V/30Watt, second-hand.  690,00 EUR


1005157 ZEISS OPMI optics unit, for use on an existing Zeiss tripod or for use in technical applications. 3 articulated axes, cold light cable connection with mirroring (the cold light cable is not included). Front lens 200mm, eyepieces 20x, tube 160mm, zoom electric 0.5-2.5. field of view 20mm-5mm, magnification 8-40x, second-hand.  780,00 EUR


#### **Furniture**












1004125 BAISCH desk with door and 3 drawers, H83cm x W 120cm x D 60cm, moss green, very good condition, second-hand  780,00 EUR












#### **Monitoring**












#### **Datex**

1005272 GE Carescape B650, network-capable multi-parameter monitor with touchscreen, delivery incl. E-sCaio module for measuring anaesthetic gases, oxygen and CO2, E-PSMP-01 module for measuring ECG, SPO2, temperature, pressure, NIBP, year of manufacture 2013, second-hand  6.400,00 EUR












1003386	DATEX Cardiocap 5: Anaesthesia monitor with the functions: ECG, NIBP, SpO2, spirometry, O2, CO2, and N2O, including new battery, ECG cable, SPO2 sensor, NIBP cuff with tube, spirometry sensor, second-hand		3.300,00 EUR
1003389	DATEX Cardiocap 5: ECG monitor with SpO2 and NIBP, used		1.600,00 EUR
1004672	DATEX Ohmeda F-FM 00 ECG monitor with color LCD, incl. module for ECG, SPO2, NIBP, TEMP measurement, battery operated, incl. DATEX N-FC-00 (Co2 Modul), used		3.900,00 EUR
1003632	DATEX Ohmeda F-FM 00 ECG monitor with color LCD, incl. module for ECG, SPO2, NIBP, TEMP measurement, battery operated, used		3.300,00 EUR
1002781	DATEX S/5 Light, Portable ECG-Monitor with battery operation, ECG, Resp, SpO2, NIBP, temp, second-hand		1.050,00 EUR
1004400	DATEX S/5 Light, Portable ECG-Monitor with battery operation, ECG, Resp, SpO2, NIBP, PRESS, temp, second-hand		1.200,00 EUR
1004461	DATEX S/5 Light, Portable ECG-Monitor with battery operation, ECG, Resp, SpO2, NIBP, temp and CO2 (without sensor) second-hand		1.450,00 EUR
1004587	DATEX zusätzlicher Unterbau-Akku (für DATEX S/5 light) incl. 2 Stk. Ni-Cd Akkus (NEU), Betriebsdauer wenn vollständig geladen: 2 Std., gebraucht		250,00 EUR
1004854	DATEX-OHMEDA Light Monitor Power Adapter Type N-LPOWT. Car charger for Datex light series monitors, used.		165,00 EUR
1004965	DATEX OHMEDA N-LCHGR external battery charger for 2 batteries. Suitable for batteries of the Datex S5 light series, second-hand.		255,00 EUR
1003374	Datex S/5 monitor with flat screen incl. function module M-NESTR, for ECG (3-channel) SPO2, NIBP, temp., Press, used		1.590,00 EUR










1004838	DATEX OHMEDA S5 Compact, portable monitor (mains and battery operation) with M-NESTPR, ECG, SpO2, PRESS, NIBP, temperature and M-C module for capnography, second-hand.		1.980,00 EUR
1001131	DATEX OHMEDA S5 Compact, portable compact monitor (mains and battery operated) with M-ESTPR, ECG, SpO2, PRESS, Temp., second-hand anesthesia monitor		2.400,00 EUR
1003954	DATEX OHMEDA S/5, portable compact monitor (mains and battery operated) with function modules for NIBP, ECG, SpO2, PRESS, Temp., second-hand intensive care monitor		1.600,00 EUR
1004041	DATEX OHMEDA S/5, monitor with function modules ESTPR: NIBP, ECG, Resp., SpO2, Temp., M-CO: Co2, M-Rec: recorder and flat screen monitor, used		2.400,00 EUR
1000850	DATEX SATLITE SpO2 monitor with programmable alarm limits, second-hand		320,00 EUR
1001808	DATEX Satlite Plus, ECG, SPO2, second-hand		420,00 EUR
1003853	DATEX M-ESTPR Module (suitable for AS/3 and S/5) with the following functions: ECG + RESP, SpO2, P1, P2, T1, T2, second-hand		450,00 EUR
1003375	Datex M-NE12STPR, module for measurement of ECG (12-channel), SPO2, NIBP, temp., 2x Press, used		690,00 EUR
1003652	DATEX M-MEM for memory card, used		130,00 EUR
1004364	DATEX M-PP: print module, second-hand		235,00 EUR
1003093	DATEX M-ESTR module: ECG + resp., Spo2, T1, T2, used		350,00 EUR











1005050	DATEX module M-NESTR (suitable for AS / 3 and S / 5) with the functions EKG, NIBP, SpO2 T1, T2, used		580,00 EUR
1003819	DATEX E-PSMP-00: module for ECG, SpO2, Temp. (T1, T2), PRESS (D1, D2) and NIBP, used		1.200,00 EUR
1004918	GE E-INTPSM, with this plug-in module for the Datex S5 series, a GE E-PSM module can be operated on the S5 monitors by means of a mountable adapter, used		490,00 EUR
1004949	GE N-FREC-00 printer module for use gei GE F-FM monitors, second-hand		350,00 EUR
1004769	DATEX Ohmeda E-P, plug-in module for measuring invasive blood pressure, second-hand		390,00 EUR
1005236	GE E-PP, pressure module with 2 inputs for invasive blood pressure measurement, second-hand		420,00 EUR
1003095	DATEX module E-PRESTN: SpO2, T1, T2, NIBP, ECG, PRESS (D1, D2) used		980,00 EUR
1003936	DATEX M-MRI: module for ECG, SPO2 and NIBP, MRI-compatible, used		650,00 EUR
1001802	DATEX M-NIBP Module with the NIBP function, second-hand		350,00 EUR
1004991	DATEX E-CAIO plug-in module for DATEX S5 monitors for measuring CO2, halothane, isoflurane, enflurane, sevoflurane, desflurane, O2, used		2.300,00 EUR
1003719	DATEX E-CAIOV plug-in module for DATEX S5 monitors for measuring CO2, halothane, isoflurane, enflurane, sevoflurane, desflurane, spirometry O2, second-hand		2.600,00 EUR

1002905	DATEX M-CAIO plug-in module for DATEX AS3 and S5 compact for measurement of 5 anaesthesia gases, CO2 and O2, used		1.900,00 EUR
1003092	DATEX M-P module: Invasive blood pressure, P3, used		180,00 EUR
1002800	DATEX M-COP MODULE: Cardiac output, second-hand		380,00 EUR
1004671	GE N-FC, CO2 plug-in module suitable for the GE monitor series GE-F FM measurement of the CO2 concentration in the bypass mode, second hand		980,00 EUR
1004670	GE E-miniC, CO2 plug-in module suitable for the monitors of the S5 series. Measurement of the CO2 concentration in the bypass mode, used		980,00 EUR
1001803	DATEX M-REC modul (Printer), second-hand		350,00 EUR
1003091	DATEX hollow module, used		25,00 EUR
1003720	DATEX blank module, white, used		30,00 EUR
<b>Dräger</b>	1005136 DRÄGER MS13836, Bezel Lens Delta XL 12 Pin, bezel for monitor Dräger Infinity Delta, as new		395,00 EUR
	1002724 DRÄGER Infinity Delta, Multi-Parameter Monitor for ECG, SPO2, NIPB, IPB, Temp., charging station, inclusiv rechargeable battery (NEW), second-hand		2.980,00 EUR
	1005128 DRÄGER Infinity Delta, multi-parameter monitor for ECG, SPO2, NIPB, IPB, temperature, charging station, incl. battery (new) and chassis with power supply, second-hand		3.100,00 EUR

























1002348	DRÄGER Infinity Gamma XL, network-compatible multi-parameter monitor (ECG, respiration, SPO2, temp., NIPB) incl. 1 battery (NEW), good condition, keyboard in English, second-hand		2.890,00 EUR
1004983	DRÄGER IDS docking station, for use with Dräger monitors e.g. Infinity Gamma and Infinity Delta, second-hand		380,00 EUR
1005121	DRÄGER 7256931 MIBII Proto CONV-LS, with this converter a Dräger patient monitor of the INFINITY series can be connected to an anaesthesia or ventilation device. The data are then displayed on the patient monitor, used		190,00 EUR
1004785	DRÄGER MS16356 SPO2 POD Masimo X8. With this pod a Dräger Infinity monitor with X8 interface can be converted to the Masimo SPO2 system, as new.		175,00 EUR
1005009	DRÄGER 5589663 MultiMed 12 Pod, 12-channel ECG module for the Infinity Gamma and Delta monitor series, incl. patient cable and connection cable to the monitor, second-hand		265,00 EUR
1003540	DRÄGER 7862696E EtCO2 Microstream module for use with DRÄGER Infinity Delta monitors, used		680,00 EUR
1004883	DRÄGER Infinity HemoMed, module for measuring invasive blood pressure. Suitable for Dräger monitors of the Infinity series, second-hand.		490,00 EUR
1004757	DRÄGER Vamos plus: anaesthetic gas monitor for monitoring the CO2 concentration, the concentration of N2O and the volatile anaesthetics halothane, enflurane, isoflurane and sevoflurane, with automatic gas detection, second-hand		3.800,00 EUR
1001812	DRÄGER PM 8010 patient monitor (ECG, SpO2, NIBP, temp), second-hand		980,00 EUR
1001813	DRÄGER PM 8014 patient monitor, measurement of: ECG, SpO2, NIBP, temperature, used		850,00 EUR
1003251	DRÄGER PM 8050 respiratory gas monitor, airway pressure, expir. breathing volume, insp. O2 conc., anaesthesia conc. in the bypass procedure for halothane, enflurane, isoflurane, desflurane, sevoflurane, CO2, with SPO2, very good condition, second-hand		3.400,00 EUR

1001827	DRÄGER PM 8050 respiratory gas monitor, airway pressure, expir. breathing volume, insp. O2 conc., anaesthesia conc. in the bypass procedure for halothane, enflurane, isoflurane, desflurane, sevoflurane, CO2, very good condition, second-hand		2.900,00 EUR
1001814	DRÄGER PM 8060 Vitara ECG monitor with module box, recording of ECG, IBP, NIBP, temperature, SPO2, good condition, second-hand		1.890,00 EUR
1005222	DRÄGER PM 8060 respiratory gas monitor with module box, measurement of ECG, IBP, NIBP, temperature, SPO2, insp. O2, exp. airway pressure, CO2, anaesthetic gases (halothane, enflurane, isoflurane, desflurane, sevoflurane). respiratory volume, airway pressure, CO2, anaesthetic gases (halothane, enflurane, isoflurane, desflurane, sevoflurane), used		2.980,00 EUR
1001144	DRÄGER PM 8060 respiratory gas monitor with module box, recording of ECG, IBP, NIBP, temperature, SPO2, insp. O2, exp. respiratory volume, positive airway pressure, CO2, anaesthesia gases (Halothane, Enflurane, Isoflurane, Desflurane, Sevoflurane), very good condition, second-hand		1.580,00 EUR
1004959	DRÄGER PM 8040, monitor for installation in a Dräger anaesthesia device, e.g. Dräger Cato. Delivery includes module box. Measurement of ECG, SPO2, NIBP, pressure, temperature, used.		890,00 EUR
1001832	DRÄGER Oxysat 2, pulse oxymeter, second-hand		510,00 EUR
1001817	DRÄGER IRIS: anaesthetic monitor to measure the concentration continuously and monitor the inhalation anaesthetics halothane, enflurane and isoflurane in the fresh gas supply from inhalation anaesthesia devices. Measurement range 0 to 9.9 Vol. %, measurement limits can be adjusted individually, second-hand		580,00 EUR
1001818	DRÄGER Barolog A, electronic device firstly to test the anaesthesia device incl. the circuit component for leakages (strongly recommended), secondly to measure and monitor the ventilation pressure continuously and thirdly to measure the ventilation frequency, second-hand		290,00 EUR
1001819	DRÄGER Spirolog, electronic device to determine and display the ventilation parameters: respiratory minute volume, breath volume, ventilation frequency, technical data: respiratory minute volume display: 0-25l/min, frequency display: 0 - 60 min <sup>-1</sup> , breath volume display 0-2l, minimum required breath volume: 0.15 l, warning display: visible and acoustic - The device is designed for use with ventilators and during spontaneous respiration, second-hand		290,00 EUR












1001840	DRÄGER Empty module with the size of Iris, Spirolog, Capnolog ecc, used		120,00 EUR
1005169	DRÄGER Leer-Modul VET, size like Iris, Spirolog, Capnodig etc., used. This article is exclusively intended for veterinary medicine and export !		89,00 EUR
1001822	DRÄGER Anemone: (Ventilation Monitoring) Measuring instrument and monitor for ventilatory parameters: Airway pressure, expiratory flow, inspiratory O2 concentration, inspiratory breathing gas temperature, used		550,00 EUR
1001826	DRÄGER PM 8030 respiratory gas monitor, airway pressure, expir. breathing volume, insp. O2 conc., anaesthesia conc. in fresh gas for halothane, enflurane, isoflurane, second-hand		1.250,00 EUR
1001823	DRÄGER PM 8030 SD respiratory gas monitor, airway pressure, expir. breathing volume, insp. O2 conc., second-hand		750,00 EUR
1001821	DRÄGER Oxydig: O2 measurement and warning device used in oxygen tents, incubators, ventilators and anaesthetic devices, etc. Measuring range 0-100 Vol.% O2, spiral sensor cable can be extended to a maximum length of 1.5 m, dimensions of the display instrument 170 x 93 x 65 mm, weight 0.5 kg, sensor dimensions: diameter 40 mm, H 45 mm, alarm is activated visually and acoustically when the individual adjustable warning limits are exceeded, second-hand		620,00 EUR
1004923	DRÄGER holder for DRÄGER Oxydig with clamping device for standard rail		45,00 EUR
1001833	DRÄGER Volumeter 2000, to determine the respiratory minute volume and respiratory stroke volume during assisted and controlled ventilation, during anaesthesia and to determine the vital capacity 3.5 to 120 l/min., second-hand		198,00 EUR
1001834	DRÄGER Volumeter 3000, to determine the respiratory minute volume and respiratory stroke volume during assisted and controlled ventilation, during anaesthesia and to determine the vital capacity 3.5 to 120 l/min., the location of the bodies of rotation in the Volumeter 3000 can be lubricated with the regreasing set 2M 18180, second-hand		265,00 EUR
1001835	DRÄGER regreasing set for Volumeter 3000. Comprises oil pipette and refill bottle, new		9,80 EUR

**Hellige /  
Marquette**












1001839	DRÄGER volumeter heating system, second-hand		90,00 EUR
1001837	DRÄGER pressure gauge, a pressure gauge is recommended as a monitoring instrument for anaesthesia (indispensable during artificial respiration) and to check your anaesthesia device incl. circuit component for leakages, se cond-hand		350,00 EUR
1003942	Pressure gauge für DRÄGER, a pressure gauge is recommended as a monitoring instrument for anaesthesia (indispensable during artificial respiration) and to check your anaesthesia device incl. circuit component for leakages, se cond-hand		198,00 EUR
1001838	DRÄGER PRECOM pressure gauge with acoustic alarm, second-hand		380,00 EUR
1001169	GE Dash 4000 multi-parameter monitor, battery operated, ECG, SPO2, temp. NIBP, 2x invasive blood pressure, very good condition, second-hand		2.850,00 EUR
1002465	MARQUETTE Dash 3000 Multi-Parameter Monitor with battery operation, ECG, SPO2, Temp. NIBP, second-hand		1.900,00 EUR
1004934	GE Dash 2500, portable battery/power supply patient monitor for the parameters ECG, SPO2, NIBP, used		3.450,00 EUR
1004593	GE Dash 5000, portable mains/battery powered ECG monitor with colour display, measurement of ECG, SPO2, Masimo, Nibp, temperature, PRESS, second-h and		3.200,00 EUR
1005002	GE Capnoflex LF CO2 module 2013427-001,Module for CO2 measurement in the bypass process for the GE Dash 3000, 4000, 5000 monitors, second-hand		680,00 EUR
1005246	GE Dash Port 2 Docking station. This quick docking station is used to mount a Dash series monitor (3000/4000/500) to the anaesthesia machine or wall mount. Mounting by means of GCX base plate. Connections for mains power, Ethernet, used		195,00 EUR
1004911	GE Marquette SAM, smart anaesthesia module for GE TRAM Rac. Anesthesia gas module for use in Marquette Solar 8000 monitors, used		490,00 EUR

<b>Nellcor</b>	1001874	NELLCOR N-550 pulse oximeter, adjustable alarm limits, battery operated, good condition, second-hand		590,00 EUR
	1002748	NELLCOR N-560 Portable Pulse Oxymeter with battery function, year of manufacture 2009, adjustable alarm limits, second-hand		850,00 EUR
	1004034	NELLCOR NPB-295, pulse oximeter with graphic display, battery powered, adjustable alarm limits, incl. finger clip sensor for adults, used		490,00 EUR
	1001872	NELLCOR N-100 anaesthesia monitor, for pulse oxymetry, adjustable alarm limits, second-hand		590,00 EUR
	1003115	NELLCOR N-3000 pulse oximeter with PSS1 charging device, battery/mains operated, adjustable alarm limits, incl. finger clip sensor, used		390,00 EUR
<b>Siemens</b>	1002992	SIEMENS Hemomed module box for 4 PRESS connections, used		480,00 EUR
	1002310	SIEMENS HEMO2: Transducer for SC9000 with the following connections: 2x PRESS, Temp A, C.O./Temp B, IBP sc 9000, second-hand		480,00 EUR
	1002308	SIEMENS remote control for SC monitors, second-hand		225,00 EUR
	1001158	SIEMENS R50 alarm recorder for SC6000, 7000, 9000 series, writing width: 50mm, second-hand		120,00 EUR
	1005261	SIEMENS 4319567 E530U Docking station. This docking station can be used to connect a single-channel recorder to a fully-fledged IDS, as good as new		320,00 EUR
	1003605	SIEMENS Capnostat Co2 Sensor with adapter for adults and children, used		1.200,00 EUR












**Other**












1003795	CRITICARE Vital Care 506DXN, vital data monitor for blood pressure and SPO2 incl. rolling stand, second-hand		390,00 EUR
1004852	OHMEDA TruSat pulse oximeter. Portable device with battery function. Delivery incl. rubber protection, charging power supply unit, finger clip sensor for adults, second hand		680,00 EUR
1005090	Masimo Radical 7 Rainbow, pulse oximeter incl. SPO2 finger clip sensor for adults, second-hand		790,00 EUR
1005129	HEYER Vizor 10, anaesthesia monitor for measurement of ECG, SPO2, etCO2 in main current mode, NIBP, Press. temperature. Clear colour display, integrated battery, second-hand.		2.980,00 EUR
1003660	OHMEDA Biox 3740 portable pulse oximeter, for battery and mains operation, adjustable alarm limits, used		350,00 EUR
1001880	OHMEDA Biox 3700e pulse oximeter, second-hand		390,00 EUR
1003929	Masimo Radical 7 Rainbow, pulse oximeter incl. SPO2 finger clip sensor for adults, used		980,00 EUR
1004380	NONIN 2500 Series, portable, battery-operated and rechargeable pulse oximeter, delivery includes charging cradle, power supply unit and adult sensor, second-hand		295,00 EUR
1002389	Agilent Viridia V24CM1205A, monitor with module rack, including ECG/SPO2/NIPB module, second-hand		1.450,00 EUR
1003492	Agilent Viridia V24CM1205A, monitor with attached module rack, including ECG/SPO2/NIPB module, used		1.550,00 EUR
1003269	Critikon Dinamap PRO 1000, portable, mains/battery-operated monitor for measurement of ECG, SPO2, NIBP, with printer, new batteries		980,00 EUR

























1005113	WELCH ALLYN Propaq 244, transport monitor with battery operation for measuring ECG, SPO2, invasive blood pressure, temperature, second-hand		980,00 EUR
1005064	PHILIPS IntelliVue MP 20, patient monitor on a mobile stand, incl. Module M3001A for measuring ECG, SPO2, NIBP, press, mains and battery operation, keyboard and language in English, used		1.650,00 EUR
1003225	HP Omnicare CMS 24, M1204A monitor for measurement of ECG, SPO2, NIBP incl. rolling stand, used		980,00 EUR
1002414	HP Viridia CMS M 1095A year of manufacture 1999, flat screen monitor with fastening clamp and module rack, modules for measuring ECG, SPO2, NIBP, PRESS, used		1.450,00 EUR
1005065	PHILIPS IntelliVue MP 30, patient monitor, incl. M3001A module for measuring ECG, SPO2, NIBP, press, mains and battery operation, keyboard and language in English, used		1.650,00 EUR
1004469	Philips IntelliVue MP 50 Anesthesia, portable anesthesia monitor with function modules for ECG, SPO2, Temp, Press, language and keyboard in English, used		3.400,00 EUR
1004468	Philips IntelliVue MP 50 Anesthesia, portable anaesthesia monitor with function modules for ECG,, SPO2, Press, Temp, CO2. Delivery includes IntelliVue G5 M1019A gas measuring bench for measuring anaesthetic gases Language and keyboard in English, second-hand		4.900,00 EUR
1004767	PHILIPS IntelliVue MP 70 patient monitor incl. ECG module M3000A for measuring ECG, SPO2, NIBP, temperature. Menu navigation and buttons are in English, second-hand		3.900,00 EUR
1004560	PHILIPS IntelliVue X2 Model M3002A Multi-Parameter Module: ECG monitoring with three to ten electrodes, 12-lead ECG with five or ten electrodes, arrhythmia detection, SpO2 Masimo, NIBP, NPRESS and Temp., integrated accu (max. 3 hrs.), operation with external power supply (not included in price) or operation with IntelliVue MP5 to MP90 possible, with new accu, used		780,00 EUR
1004579	PHILIPS External power supply for IntelliVue X2, as new		880,00 EUR
1004841	BLT M9000A, portable patient monitor with colour display. Measurement of ECG, SPO2, NIBP and temperature, second-hand		750,00 EUR
























1004839	PHILIPS IntelliVue G5 M1019A, gas measuring bench suitable for the monitor series Intellivue MP50, for measuring the anaesthetic gases halothane, enflurane, isoflurane, sevoflurane, as well as CO2 and O2 (paramagnetic), as well as automatic anaesthetic gas detection, second-hand		1.980,00 EUR
1003973	Philips M3046A/M3 transport monitor with battery operation, including the module M3000A: Measurement of ECG, SPO2, NIBP, IBP, temp., used		1.100,00 EUR
1004928	Philips M3046A/M3 transport monitor with battery operation. Incl. module M3000A, measurement of ECG, SPO2, NIBP, IBP, temp., used		950,00 EUR
1003524	PHILLIPS Sure Signs VM 4 monitor for SPO2, NIBP measurement, incl. battery, used		780,00 EUR
1004907	NIHON KOHDEN BSM-3763 LifeScope, multi-parameter monitor for measuring ECG, SPO2, NIBP, temperature, modular smart cable technology (The multi-connector socket automatically recognizes the parameter type: BIS, IBP, HZV, CO2). Extra large color touch display, network compatible, used		2.980,00 EUR
1004926	NIHON KOHDEN flat screen monitor 19 inch, for operation on a Nihon Kohden patient monitoring system, second-hand		590,00 EUR
1003132	HP module M1008B for measurement of non-invasive blood pressure. Used, supplied with tube and cuff		390,00 EUR
1003133	HP module M1002B for ECG. Used, supplied with ECG cable		390,00 EUR
1003134	HP module M1020A for measurement of SPO2. Used		490,00 EUR
1004716	PHILIPS M1029 Temp. module, NEW		250,00 EUR
1003135	HP module M1006B for measurement of invasive blood pressure. Used, supplied with press master cable		340,00 EUR

1004914	HP M1016A CO2 module. Plug-in module for capnography. Delivery without sensor, used		690,00 EUR
1003938	Phillips M3000A: module for ECG/Resp, SpO2, NIBP, Temp and PRESS, used		780,00 EUR
1005107	Philips M3001A: Module for ECG/Resp, SpO2, NIBP, Temp and PRESS, for operation with Philips monitors of the IntelliVue MP series, second-hand		1.100,00 EUR
1003928	NIHON-Kohden BSM-2303K, network-capable patient monitor for measurement of ECG, SPO2, Temp. with color touchscreen, portable with battery, used		1.980,00 EUR
1001882	PULSION Picco: Monitor for low-invasive monitoring of haemodynamic parameters: cardiac output (CO), cardiac function index (CFI), intrathoracic blood volume (ITBV), extravasal water on the lung, cardiac output from pulse correction, stroke volume (SV), second-hand		980,00 EUR
1005120	MINDRAY 115-011185-00, CO2 module for measuring CO2 in the bypass method, as new.		1.450,00 EUR
1003873	DRÄGER / SIEMENS remote control for DRÄGER Infinity monitors for control of NIBP, SPO2, used		180,00 EUR
1001867	Siemens ECG cable for Sirecust series 400, 341, 960, 1281. Complete cable with 3-pin clamp adapter, red, yellow, green, new. The cable is also available with crocodile clamps for use in veterinary medicine, NEW		105,00 EUR
1004865	DRÄGER compatible ECG cable, 3-wire with clip attachment for adhesive electrodes, NEW		89,00 EUR
1003261	Temperature probe for SIEMENS and DATEX monitors, used		85,00 EUR
1003595	DRÄGER Y-piece with mounting for temperature sensor, used		12,00 EUR

1003596	DRÄGER temperature sensor for PM8030 and 8050, used		89,00 EUR
1002934	DATEX ECG cable for Cardiocap etc., complete 3-wire cable. The cable is also available with crocodile clamps for use in veterinary medicine (addition 9.00 euros), NEW		112,00 EUR
1002403	GE Capnostat Co2 Sensor with adapter for adults and children, used		980,00 EUR
1002935	NELLCOR compatible SPO2 finger clip sensor for adults, NEW		125,00 EUR
1002936	DATEX compatible SPO2 finger clip sensor for adults, NEW		164,00 EUR
1004868	DATEX compatible adult SPO2 finger clip sensor, for connection to monitor trunk cable, NEW		108,00 EUR
1002937	OHMEDA compatible SPO2 finger clip sensor for adults, NEW		139,00 EUR
1001875	NELLCOR replacement veterinary sensor for Nellcor NPB-40 pulse oxymeter with 2 clips, new, - This item is intended solely for veterinary medicine-		239,00 EUR
1002938	DATEX Vet. sensor with tongue clip for horses, dogs and cats, NEW suitable for Cardiocap, Capnomac, Satlite etc.		199,00 EUR
1004870	OHMEDA compatible veterinary sensor for use on dogs and cats (for use with OHMEDA pulse oximeters e.g. Biox etc.), NEW		125,00 EUR
1004692	Skin temperature probe, second-hand		65,00 EUR










1004804	Temperature probe reusable for DATEX Cardiocap 5 and Cardiocap 2. Especially rigid probe for nasal application in horses, length 300cm, diameter 4mm, new		165,00 EUR
1002939	Blood pressure cuff for adults or children, single tube, suitable for SIEMENS, DRÄGER, HP, etc. including adapter (please state type), NEW		32,00 EUR
1002940	Blood pressure cuff for adults or children, double tube, suitable for DATEX, CRITIKON, etc. including adapter (please state type), NEW		32,00 EUR
1002941	DATEX NIBP pressure tube for Cardiocap II etc. adults, NEW		75,00 EUR
1002942	SIEMENS PRESS master cable for inclusion of BRAUN Combitrans monitoring set, NEW		90,00 EUR
1004159	DRÄGER PRESS master cable for inclusion of BRAUN Combitrans monitoring set, NEW		90,00 EUR
1002943	DATEX PRESS master cable for inclusion of BRAUN Combitrans monitoring set, NEW		90,00 EUR
1001854	BRAUN Combitrans Monitoring set venous: Monitoring set with disposable transducer for physiological pressure measurement with blue pressure hoses, NEW		32,00 EUR
1002944	CO2 sample tube 2 metres male/male, NEW		3,70 EUR
1004902	LUER-LOCK adapter female/female. For connecting 2 CO2 sample tubes with male/male connectors, New		2,00 EUR
1004794	CO2 sample hose NAFION male/male, This hose serves as a drying stage for the gas mixture to be measured (before it reaches the actual sample hose and the water trap) to prevent possible clogging of the system. -this item is intended exclusively for veterinary medicine-		55,00 EUR

1004816	SURESTREAM Sure VentLine H Set: Adult-pediatric CO2 sampling line and airway adapter for humid environments, NEW		59,00 EUR
1004558	DATEX D-fend reusable moisture separator BLACK for DATEX Capnomac, Cardiocap 5 etc., NEW		31,00 EUR
1004877	DATEX mini D-fend, water trap for CO2 modules (e.g. GE E-miniC), NEW		28,00 EUR
1002966	DRÄGER WaterLock: Moisture separator for PM8050, NEW		22,00 EUR
1002967	DRÄGER WaterLock2: Moisture separator for PM 8050, Vamos NEW		33,00 EUR
1003436	DRÄGER moisture separator for PM 8050 (old model) or Cicero, NEW		38,00 EUR
1005109	DRÄGER water trap for CO2 monitor Capnodig, used		119,00 EUR
1004894	Upgrade kit for DRÄGER PM 8050, PM 8060, Cato (if the fast cell" DRÄGER 8650930 is "no" longer available) consisting of: DRÄGER sensor" cable with housing, DRÄGER valve cap with adapter for above mentioned housing, ENVITEC sensor D		384,00 EUR
1004611	DRÄGER 6871500: CO2 main current sensor for Evita, second-hand		180,00 EUR
1004543	Dräger 5740738E etCO2 Pod incl. main current sensor. Module for the measurement of etCO2 using the main current method on the patient's tube. This pod is suitable for operation with Dräger Infinity Gamma XL, second-hand.		980,00 EUR
1002393	DRÄGER 8403735 Flow-Sensor for Spirolog, Evita, NEW		33,00 EUR










1002946	ENVITEC Sensor D replacement product for DRÄGER 6850645 for Oxydig, Babylog, EVITA etc., NEW		86,00 EUR
1002382	ENVITEC O2 Sensor for SIEMENS Servo ventilators ENVITEC V300, new		95,00 EUR
1002947	DATEX d-lite adult cuvette for spirometry, reusable, NEW		46,00 EUR
1002948	DATEX Pedi-lite cuvette for children for spirometry, reusable, NEW		46,00 EUR
1002950	DATEX dust filter for Cardiocap/Capnomac, NEW		4,90 EUR
1001298	Veterinary adapter for ECG patient cable: For patients, crocodile clamp, with snap-on connection for ECG cable, NEW		9,50 EUR
1001299	Crocodile clamp (with screw) for ECG cable		3,00 EUR
1004872	VET ECG clamps, suitable for most ECG cables: 3mm banana plug, 4mm banana plug, clip. Gentle wide jaws with large contact surface, set of 5, NEW		42,00 EUR
1004647	NUOVA S-D13 / P disposable pulse oximetry sensor suitable for Nihon Kohden appliances pediatric 10kg-50kg, new		5,50 EUR
1004648	NUOVA S-D13 / I disposable pulse oximetry sensor suitable for Nihon Kohden infant 1kg-20kg, new		5,50 EUR









## Anaesthesia and ventilation








**Ventilation equipment**









1005245	DRÄGER Evita Infinity V500, intensive care ventilator for adults, children and neonates. The Evita Infinity V500 offers different mandatory ventilation modes as well as different modes to support spontaneous breathing. Airway monitoring is displayed on the Dräger Medical Cockpit Infinity C500, used		9.500,00 EUR
1002432	DRÄGER Evita XL: Long-term ventilator for intensive care. For adults, children and neonates, basic unit with hoses for ZGV (central gas supply) for O2 and Air. Disposable patient tubing. Ventilation modes: IPPV, SIMV, MMV, SB, CPAP, ASB, BIPAP, BIPAPAssist, APRV, PPS, ILV. used		4.900,00 EUR
1005242	FISHER & PAYKEL MR850AEK, respiratory humidifier for long-term ventilators. Delivery incl. 1 set of accessories, second-hand		450,00 EUR
1005243	FISHER & PAYKEL HC550AGU, airway humidifier for long-term ventilators, delivery includes 1 set of accessories, second-hand		450,00 EUR
1004756	SIEMENS Servo-i long-term ventilator for adults, with holder for standard rail (both for the device and for the monitor) ventilation modes: VK, DK, PRVC, SIMV (VK) + DU used with NEW battery		9.800,00 EUR
1004720	SIEMENS Servo Screen 390: monitor for servo 900 C/E or 300/300A, with monitor holder, second-hand		890,00 EUR
1001554	SIEMENS Servo 900 C ventilator, ventilation pattern: volume control rt, volume control with sighs, pressure-controlled, pressure-supported, SIMV, SIMV + pressure support, CPAP manual, min. vol. 0,5-40 l/min, respiratory rate: 5-120 breaths/min, similar display: ventilatory pressure, exp. minute vol., dig. display: respiratory rate, O2 concentration, insp. tidal volume, exp. tidal volume, exp. minute volume, peak, intermittent, moderate, max. working pressure: 120 mbar, 50x23x32cm, weight : 19kg , from bodyweight of 10 kg , second-hand		1.980,00 EUR
1005211	DRÄGER 5704281 Wall bracket for Oxylog 2000plus, second-hand		245,00 EUR
1001559	DRÄGER Oxylog 2000 transport ventilator with accessories, time-controlled, constant volume for patients with a tidal volume from 1000 ml ventilation mode: IPPV \ SIMV \ CPAP. Variables 1:E, also IRV, with integrated PEEP valve 0-15 mbar. Air-Mix \ No-Air-Mix switch, integrated expiratory volume measurement with display, second-hand		2.450,00 EUR
























1001556	DRÄGER Oxylog, time-controlled, constant volume ventilator for controlled respiration in emergency medicine. Including ventilation accessories comprising: silicone breathing tube, respiratory valve and carrying strap, second-hand		690,00 EUR
1005208	WEINMANN Medumat standard, portable emergency ventilator with oxygen module and carrying frame. A 2-litre O2 bottle can be attached to this carrying frame, including pressure reducer, second-hand.		1.680,00 EUR
1005125	DRÄGER carrying frame for Dräger Oxylog 2000. In addition to the ventilator, a 2-litre oxygen cylinder can also be attached. (the pictured Dräger Oxylog 2000 is not included in the price), second-hand.		240,00 EUR
1004745	DRÄGER carrying bag for Oxylog 3000: 1 large bag with carrying strap and 2 small removable bags, second-hand		320,00 EUR
1004773	DRÄGER carrying bag for Oxylog with accessories, second-hand		280,00 EUR
1005115	BLEASE 2200 MRI, anaesthesia ventilator for use at the MRI scanner, incl. patient valve, second-hand		590,00 EUR
1001542	DRÄGER Evita 2 time-controlled, constant volume ventilation device starting with a respiratory volume of 50 mL. controlled assisted ventilation IPPV, IPPVAssist - synchronised intermittent mandatory ventilation SIMV, - mandatory minute volume ventilation MMV, - spontaneous ventilation under continuous positive airway pressure CPAP, BIBAP, second-hand		1.980,00 EUR
1002396	DRÄGER Expiration Valve, patient system for Evita2, Evita 2 dura and Evita 4, used		395,00 EUR
1004998	Expiration valve for Evita 2 dura, 4 edition, XL, disposable, new		22,00 EUR
1005229	DRÄGER MP01061 Expiration valve one-way for long-term ventilator Evita Infinity V500, new		36,00 EUR

1001527	DRÄGER Ventillog ventilation device for small animals and foals,, incl. accessories, powered with oxygen or compressed air, for controlled ventilation, constant volume, time-controlled, tidal volume 50 - 150 ml (small bellows) and 150 - 1600 ml (large bellows), frequency 6-60 min, used, Inspirat: flow 20-80 l/min, compressed gas supply: 2-6 bar, compressed gas consumption: 30 l/min, used. The small bellows are not included in the price		780,00 EUR
1001528	DRÄGER Ventillog 2 ventilation device, incl. accessories, powered with oxygen or compressed air, for controlled ventilation, constant volume, time-controlled, tidal volume 50 - 150 ml (small bellows) and 150 - 1600 ml (large bellows), frequency 6-60 min, used, Inspirat: flow 20-80 l/min, compressed gas supply: 2-6 bar, compressed gas consumption: 16 l/min, used, selectable I : E ratio = 1:1, 1:2, 1:3, used. The small bellows are not included in the price		980,00 EUR
1003087	DRÄGER Ventillog 3: Anaesthesia ventilation device, pneumatic driver, time-controlled ventilation rate: 6 - 60 / min I:E ratio: 1:1, 1:2, 1:3 Inspiration flow: 20 to 80 L/min at 20 mbar stroke volume 30 to 150 ml (bellows for children) 100 to 1400 ml (bellows for adults) minute volume max 25 L/min, compressed gas supply up to 2 bar (O2 or oil-free compressed air), 16 L/min, used. The small bellows are not included in the price		1.380,00 EUR
1003341	DRÄGER Ventillog C: Anaesthesia ventilation device, pneumatic driver, time-controlled ventilation rate: 6 - 60 / min I:E ratio: 1:1, 1:2, 1:3 Inspiration flow: 20 to 80 L/min at 20 mbar stroke volume 30 to 150 ml (bellows for children) 100 to 1400 ml (bellows for adults) minute volume max 25 L/min, compressed gas supply up to 2 bar (O2 or oil-free compressed air), 16 L/min, used. The small bellows are not included in the price, only for operation with COSY SYSTEM, used		1.480,00 EUR
1001529	DRÄGER Ventillog ML with circuit system ML, improved Ventillog 2 with modified circuit system 9, switching valve incorporated in the Ventillog ML patient system, circuit system with passive unidirectional valve RV and PEEP valve, second- hand		1.900,00 EUR
1003965	DRÄGER patient component for Ventillog (2), used		280,00 EUR
1003089	DRÄGER switching valve for Ventillog and Ventillog 2, used		180,00 EUR
1003088	DRÄGER FGE valve: Switching valve for Ventillog 3, used		580,00 EUR

1004317	DRÄGER Bellows for children (250 ml) for Ventilog or RA2, second hand		350,00 EUR
1004181	Scale for child bellows in Dräger Ventilog, NEW		29,00 EUR
1004551	Bellows for children (250 ml) for Ventilog or RA2, NEW CAUTION nor original DRÄGER product, This item is exclusively designed for being used in the field of veterinary medicine		298,00 EUR
1001531	DRÄGER Bellows for adults(1.600 ml) for Ventilog or Pulmomat, second-hand		45,00 EUR
1001525	DR. MÜLLER horse ventilator, large animal ventilator, ventilator for horses: pressure-controlled ventilator for horses based on the bag-in-bottle system, with BIRD Mark 7 (sec. hand) ventilation unit that generates pressure and vacuum in an acrylic glass cylinder. Manual, assisted or controlled ventilation. The device is mobile, extremely robust (4-wheel trolley with acrylic glass bin), can be easily dismantled, stand-alone device (for every large animal - anaesthesia device), new, Suction capacity: 370l/min - intended solely for veterinary medicine - information: MyClean® FD: Disinfectant for Dr. Müller horse ventilator (100g contain: 2,5g Polyhexamethylenbiguanidhydrochlorid, 0,8g Didecylimethylammoniumchlorid) We recommend an annual service of the complete machine. This concerns in particular the BIRD Mark 7		4.900,00 EUR
1001532	DRÄGER Babylog N anaesthesia ventilation device for baby animals weighing 50 g - 15 kg, Continuous-Flow-Respirator with time control and/or manual release, max. ventilatory pressure: 10-60 mbar, expiration pressure - PEEP/CPAP 0-10 mbar, insp. time: 0.3-2 sec., exp. time: 0.5-6 sec, driving gas O2 or oil-free compressed air 2-6 bar, max. gas consumption 2l/min., H: 176 mm, x W: 212 mm x D: 300 mm, weight: 5 kg, second-hand		680,00 EUR
1001535	DRÄGER Babylog 2, short-term ventilation device for premature babies and newborn infants to carry out controlled ventilation (IPPV) with PEEP or spontaneous respiration with CPAP during transportation or for resuscitation, inspiration time-controlled or triggered manually, pressure-limited expiration is time-controlled, inspiration pressure 10-60 mbar plateau pressure, expiration pressure 0-10 mbar, ventilatory rate 10-60/min., driving gas O2 or dried, oil-free compressed air 2-6 mbar, H:195 mm x W:120 mm x D:200 mm, mit DRÄGER Babymix 2 O2 mixer, second-hand		780,00 EUR

1000805	DRÄGER Babylog 2000: Ventilator for premature and newborn infants with a bodyweight of up to 6 kg, with integrated gas supply mixer function - oxygen and compressed air. For the following ventilation modes: controlled ventilation IPPV, intermittent mandatory ventilation IMV, spontaneous ventilation under continuous positive airway pressure CPAP With monitoring of positive airway pressure PAW, electricity supply to the monitoring module and the gas supply, second-hand		1.980,00 EUR
1000806	DRAEGER Babylog 8000: respirator for premature, newborn and Children up to 20kg body weight, controlled ventilation (IPPV), synchronized ventilation (SIPPV), intermittent mandatory ventilation (IMV), synchronized intermittent mandatory ventilation (SIMV), spontaneous breathing under positive airway pressure (CPAP), monitoring for insp. O2 concentration FiO2, airway pressure Paw, Flow V, minute volume MV, respiratory volume VT, dimensions basic device: (WxHxD) 212x280x390mm, 14.5kg, used		4.900,00 EUR
1001536	DRÄGER Babylog 8000 plus: Ventilator for premature babies, newborn infants and children up to 20 kg in bodyweight, controlled ventilation (IPPV), synchronised ventilation (SIPPV), intermittent mandatory ventilation (IMV), synchronised intermittent mandatory ventilation (SIMV), spontaneous respiration with continuous positive airway pressure (CPAP), monitoring for insp. O2 concentration FiO2, positive airway pressure PAW , Flow V, minute volume MV, respiratory volume TV, dimensions of basic device: (WxHxD) 212x280x390 mm,14.5kg,second-hand		5.200,00 EUR
1004160	DRÄGER monitor for operation on Babylog 8000, second-hand		580,00 EUR
1004992	DRÄGER 8408950 Patient part for Dräger Babylog 8000, incl. silencer and hose connection, second-hand		690,00 EUR
1001537	DRÄGER accessories for Babylog 8000: Flow sensor, Y piece and connection cable, NEW		465,00 EUR
1004892	DRÄGER flow sensor cable to accommodate the flow sensor for Babylog 8000, New		319,00 EUR
1005138	DRÄGER 8412935, medical drug nebuliser white for use with Dräger EVITA 4, XL, and Savina, as new		139,00 EUR












1002798	External battery for DRÄGER Evita 4: 24V, max. running time 2 Hrs., H 13 x W 30 x D 40,second-h and		680,00 EUR
1002244	DRÄGER CF 800 Continuous Flow CPAP System: Flow meter for oxygen and air, incl. bellows for manual ventilation, with bracket for fastening on a standard rail. For CPAP and O2 therapy with adults and children from 3 kg weight, for the support of spontaneous respiration, second-hand		680,00 EUR
1003875	DRÄGER 5730460 Converter DC/DC 12 for DRÄGER Respicare, as-new , used		165,00 EUR
1001548	DRÄGER Aquapor: Respiratory gas humidifier with continuous temperature adjustment (39 °C to 66 °C), respiratory minute volume up to 30 l/min, very good condition, second-hand		450,00 EUR
1003362	SIEMENS Servo Ultra Nebuliser 345 for Servo 300, NEW,		340,00 EUR
1004675	Ambu pouch for toddlers, with peep valve		45,00 EUR
1004788	AMBU Silicone Resuscitator for adults with Ambu demand valve, used.		198,00 EUR
1002307	Ambu bag set incl. PEEP valve, 2 masks, mask band, gas leads and cables etc., NEW		82,00 EUR
1002477	Earth Cable, second-hand		25,00 EUR
1004743	DRÄGER medical socket strip 4-fold for attachment to an anesthesia device Dräger or to a ventilator, used, the strip can also be used elsewhere		45,00 EUR
1002478	Earth Cable with triple socket outlet, second-ha nd		35,00 EUR












<b>Equipment fittings</b>	1000867	DRÄGER secretion bottle package consisting of carrying frame, 2x polysulfone container 700ml and lid, used.		125,00 EUR
	1004668	Anesthesia mask for dogs, material PVC, diameter 136 mm, length 140 mm, opening of elastic rubber cover 54 mm, NEW		39,00 EUR
	1004431	Anesthetic mask for mice and small rodents: Perspex dome with separate intake and exhaust of the anesthetic gas. a latex membrane with small hole is clamped on the cylinder with an O-ring, PVC hose for inspiration and exhalation, NEW. This item is exclusively intended for veterinary use.		260,00 EUR
	1001601	Anaesthesia mask for mice and small rodents, acrylic glass, latex, diameter: 48mm, diameter opening: 24mm,, incl. gas mains adapter, new - This item is intended solely for veterinary medicine -		25,00 EUR
	1001603	Anaesthesia box for rodents made from acrylic glass measuring 300 x 150 x 160 mm, removeable cover with handle, 2 standard connections at the front sides (diameter 15 mm internal/external 22 mm), neu		179,00 EUR
	1004624	3-way stopcock for fresh gas hose, NEW; incl. overpressure valve NEW (this valve must be installed as safety, if accidentally both outputs of the 3-way valve should remain closed when the flowmeter is open.) In the picture you can see this article integrated into another anesthesia machine.		398,00 EUR
	1001604	DRÄGER/AESULAP Intubation equipment, unused, partly in original packaging, in very good condition, comprising: 1 battery handle, 4 laryngoscope spatulas in various sizes, with 6 stainless steel tube adapters, 1 arterial clamp, 1 pair of surgical scissors, 2 replacement bulbs, 1 mouth gag amongst other items, incl. fabric bags and metal containers L 350 x B 160 x H 75 mm, used		79,00 EUR
	1001605	DRÄGER/AESULAP Intubation equipment, second-hand but in very good condition. Comprising: 1 battery handle, 2 laryngoscope spatulas, with 6 stainless steel tube adapters. One arterial clamp, one pair of surgical scissors, two replacement bulbs, one mouth gag amongst other items inc. 2 fabric bags and green metal storage container, second-hand		69,00 EUR
	1004796	VET test lung: Max. Volume 6 liters at 20mbar Suitable for testing a horse respirator or for training. With device for clamping to a standard rail, used		480,00 EUR
	1005052	SIEMENS silicone patient hose, length 30cm, diameter 22mm, used		25,00 EUR






















1005219	Silicone patient hose, length 110cm with socket 22mm on both sides, inner diameter 19mm, autoclavable, new		39,00 EUR
1004586	DRÄGER 8605088: Flexible joint arm for breathing bag incl. connection block for Cosy 2, second-hand		89,00 EUR
1005152	AMBU PEEP 20 valve, adjustable PEEP valve for Ambu hand-held resuscitator, as new.		69,00 EUR
1004786	DRÄGER MX08221-01: Evita S2 SET 2Y, new		28,00 EUR
1004787	DRÄGER MX08220-02: Evita 2/S2/2DURA/4/XL SET 1Y, new		28,00 EUR
1004790	DRÄGER 8405000-10: Evita Medicament Nebulizer (pneu) NEW		125,00 EUR
1004791	DRÄGER MX08995 Primus ServSet 6 J (plus 2 hrs. labor for installation)		2.290,00 EUR
1004792	DRÄGER MX08284 Savina ServSet 6 J (plus 2 hrs. labor for installation)		790,00 EUR
1004793	DRÄGER MX08756 Oxylog 3000 ServSet 6 J (plus 3 hrs. labor for installation)		1.270,00 EUR
1003967	DRÄGER patient component for Cicero and Cato, used		450,00 EUR
1003966	DRÄGER (Ventilator) piston pump for Cicero, used		280,00 EUR

























1004906	DRÄGER patient part for anesthesia unit Julian, second-hand		980,00 EUR
1004955	DRÄGER M33841 Bellows for Julian		185,00 EUR
1004956	DRÄGER bellows container for Julian, used		195,00 EUR
1003043	DRÄGER AGS imaging system, used		280,00 EUR
1004651	DRÄGER Cosy 1 circle system, used		980,00 EUR
1005255	DRÄGER 8605534, Power supply unit for heater COSY 2.6, second-hand		580,00 EUR
1001620	DRÄGER circuit component incl. soda lime container, inspiration and expiration valve, second-hand		750,00 EUR
1004674	DRÄGER pressure relief valve (with reversing lever) for circuit component ISO, used		98,00 EUR
1001607	Leck-Such spray to detect leaks in the anaesthesia device, new		5,00 EUR
1001621	DRÄGER soda lime container, second-hand		75,00 EUR
1004184	DRÄGER stand and lid for scrubber tank, second-hand		15,00 EUR











1002792	DRÄGER Drägersorb 800 plus soda lime for anaesthesia devices, canister with 5 Litres		42,50 EUR
1004710	DRÄGER CLIC Absorber 800+ (1.2 L), NEW		25,00 EUR
1001623	DRÄGER expiration valve, second-hand		45,00 EUR
1004725	DRÄGER expiration valve for DRÄGER Evita 2, 4, XL, used		395,00 EUR
1002123	DRÄGER sight glass for inspiratory and expiratory valves, NEW		21,00 EUR
1001625	DRÄGER valve separator for inspiration and expiration valve, 1 piece, new		27,00 EUR
1004628	SET with reusable tubing, NEW, consisting of: DRAEGER patient tube silicone 110 cm, 22 mm, Y-piece angled with LUER Lock connection, DRAEGER breathing bag 2.0 liter made of chloroprene, DRAEGER iso-spout, DRAEGER patient tube made of silicone, 50 cm, 22 mm		212,00 EUR
1003908	Schlauchkonnektor, zur Befestigung eines 11mm Patientenschlauchs an ein Dräger Inspirations-/Expirations - Ventil mit 22mm Aussendurchmesser, neu		6,50 EUR
1004900	Dräger MP00301 Anesthesia tubing set BASIC, latex free disposable tubing set consisting of 2 patient tubes breathing bag with 60cm tube, Y-piece. NEW		8,50 EUR
1001629	DRÄGER patient hose, silicone, 50 cm, 22 mm, new		51,00 EUR
1002299	DRÄGER patient hose, silicone, 110 mm, 22 mm, new		58,00 EUR

1005262	Patient tube made of silicone, length 150cm with socket 22mm on both sides, tube inner diameter 19mm, autoclavable, new		49,00 EUR
1005239	Silicone patient hose for children, length 110cm with socket 22mm on both sides, inner hose diameter 11mm, autoclavable, new		39,00 EUR
1004939	Ulmer children's set made of silicone. Consisting of 3 hoses approx. 110cm, D=9.5mm, children's Y-piece angled, breathing bag 0.5 litre, 1x ISO double nozzle, NEW		169,00 EUR
1005260	Breathing bag made of silicone 0.5 litre, sleeve 22mm, autoclavable, new		38,00 EUR
1005258	Breathing bag made of silicone 1.0 litre, sleeve 22mm, autoclavable, new		38,00 EUR
1004759	Breathing bag 2.3 liters, silicone, used		15,00 EUR
1004695	Breathing bag 5.0 liters, black rubber, NEW		19,00 EUR
1004766	Ambu bag for adults, with peep valve, used		45,00 EUR
1005151	LAERDAL Resuscitator Silicone Adult, hand resuscitator for adults, autoclavable, possibility to connect an oxygen reservoir, second-hand.		85,00 EUR
1001631	DRÄGER ISO nozzle, second hand		9,80 EUR
1001547	DRÄGER Y pieces (1x straight, 1x curved) and 2 x 7 pieces tube connectors (material: steel) (size 4, 5, 6, 7, 8, 9, 10), unused - in original container		98,00 EUR

1004893	Y-piece children straight, for use with patient tubes with a diameter of 11mm, new		15,00 EUR
1004874	Y-piece angled, for adults. Suitable for patient tubes with 22mm diameter, NEW		18,00 EUR
1005257	DRÄGER Y-piece angled with Luer-Lock connection for gas sample tubing. Suitable for patient tubes with 22 mm diameter, used		19,00 EUR
1004410	Tube adapter with LuerLock connection for attaching the CO2 sample tube to the Y-piece, second-hand		9,20 EUR
1005256	Hose connector straight, 1 outlet 22mm, 1 outlet 15mm, new		4,90 EUR
1004815	DRÄGER parking bracket for circular part (with wall support) With the aid of this system, a DRÄGER circular part can be attached to a standard rail		98,00 EUR
1004937	DRÄGER parking bracket for circular part: With the help of this system a DRÄGER circular part can be fixed to a standard rail		48,00 EUR
1002954	Kuhn System / Jackson Rees valveless, half-open anaesthesia system. Particularly suitable for small animals due to the low free volume and low respiratory resistance, can be autoclaved, new		142,00 EUR
1004186	ADAPTER for DRÄGER mixed gas hose / Jackson Rees system, NEW		23,00 EUR
1004518	Coaxial anesthesia tubing system with rotatable connection, variable expiration tubing, Luer lock connection for gas sampling (CO2 measurement) Length 1.6 m, NEW		12,50 EUR
1001642	DRÄGER anaesthetic gas exhaust air hose, rubber spiral hose incl. bypass for ambient air (prevention of a vacuum), 3 metre, second-hand		25,00 EUR






1003835	DRÄGER M33298 AGS 3 m aspiration hose, as new		49,00 EUR
1004402	DRÄGER AGFS/AGSS anesthesia gas exhaust plug connector straight , EN-ISO, NEW		103,00 EUR
1004807	AGFS/AGSS anesthesia gas exhaust plug connector, 45 degrees angled, EN-ISO, NEW		103,00 EUR
1005241	AGFS-AGSS anaesthetic gas connector, angled 45 degrees, EN-ISO, second-hand		55,00 EUR
1001643	DRÄGER standard plug NGA GD for anaesthetic gas exhaust air hose, straight or angulate (please note when ordering!), second-hand		45,00 EUR
1004510	Anaesthetic gas suction adjustable (0.5-20 litres/min), delivery includes 12V power supply unit, exhaust hose 180 cm, new		890,00 EUR
1004768	SIEMENS Evac 180: Anesthesia gas delivery system for SERVO 900		145,00 EUR
1001644	Exhaust air hose, PVC, product sold by the metre, new		2,00 EUR
1003190	Anaesthetic gas filter for anaesthesia machines, activated carbon filter for connection to the exhaust air hose. The filter must be replaced if the weight increases by 200g. This article is intended exclusively for veterinary medicine and export. NEW		45,00 EUR
1001646	Large animal circuit, circuit part for large animals, (material: stainless steel, chrome-plated copper), complete with respiratory lime container, overpressure valve, 2 patient tubes, Y-piece, 4 different adapters for tracheotubes, including airway pressure measurement, NEW, all parts are compatible with the MATRX VML Large animal anesthesia machine circuit		4.900,00 EUR
1001647	Patient hose length: 175 cm, internal diameter (at attachment) 50 mm, for horse circuit component, 1 pair, new -This item is intended solely for veterinary medicine -		295,00 EUR







1001654	Respiratory bag for circuit component, 30 litres, horse, new. - This item is intended solely for veterinary medicine -		425,00 EUR
1005228	Bellows for large animal breather Eickemeyer LAVC-2000D, used, was only used for testing purposes.		450,00 EUR
1001655	Excess pressure valve for circuit component, horse, new - This item is intended solely for veterinary medicine -		380,00 EUR
1003419	Metal Y-piece for horse patient tubes, with Luer lock (for connection of the anaesthesia gas measuring device), new, - This article is only intended for veterinary medicine -		298,00 EUR
1001657	metal Y-piece for horse patient tubes, with Luer lock (for connection of the anaesthesia gas measuring device). This item fits to the tube adapter made from metal, new, - This item is intended solely for veterinary medicine -		298,00 EUR
1004134	Holder for horse hoses with Y-piece, NEW		165,00 EUR
1001663	Cap nut made from metal for Y piece, horse, new, - This item is intended solely for veterinary medicine -		82,00 EUR
1001651	Plastic tube adapter, 1 set (4 piece) for plastic Y piece, horse, new. - This article is intended solely for veterinary medicine -		240,00 EUR
1003518	Tube adapter, diameter 22 mm, new		80,00 EUR
1003519	Tube adapter, diameter 27 mm, new		80,00 EUR
1003520	Tube adapter, diameter 32 mm, new		80,00 EUR








1003521	Tube adapter, diameter 40 mm, new		80,00 EUR
1005105	Tube adapter universal, diameter 40mm with sealing ring and union nut. With this adapter you can connect your Y-piece of the horse circle part with endotracheal tubes of the brand SurgiVet/Smith, new		189,00 EUR
1001662	Cover lid (dome for inspiration and expiration valve on horse circuit component, NEW, - This item is intended solely for veterinary medicine -		122,00 EUR
1003803	Valve disk for MATRX horse circle part, diameter: 65mm		68,00 EUR
1001664	Soda lime component for circuit component, horse, new - This item is intended solely for veterinary medicine -		455,00 EUR
1003820	Befestigungsring für Schauglas am Kreisteil Pferd, neu		88,00 EUR
1003821	Sealing ring for sight glass on circular part horse, outer diameter 103mm, inner diameter 98mm, seal thickness 3.4mm, new		18,00 EUR
1004519	DRÄGER M13386 Sealing ring for Dräger scrubber tank, 2 seals are needed for the scrubber tank (1x top, 1x bottom), NEW		34,00 EUR
1004885	Dichtring für Dräger Atemkalkbehälter, für den Atemkalkbehälter werden 2 Dichtungen benötigt (1x oben, 1x unten), NEU ACHTUNG: Dieses Teil ist KEIN DRÄGER Originalteil und darf nur in der Tiermedizin eingesetzt werden		18,00 EUR
1004520	Sealing ring for in- and expiration valve, NEW		4,50 EUR











**Anaesthesia equipment**

1004527	DRÄGER Primus inclusiv Infinity Gamma XL and vaporizers: mobile anaesthesia workstation with E-Vent Plus, electronic gas mixing, electronic ventilator, heated patient system, open platform concept and extensive interfaces for modern data communication. Fresh gas flow rate 0 - 18 l/min, ORC function, operating modes: man., spontaneous, IPPV, PCV, stroke volume 20 - 1400 ml, tidal volume. 1 - 1400 ml, respiration frequency 3 - 80/min, monitoring: O <sub>2</sub> , N <sub>2</sub> O, CO <sub>2</sub> , Hal., Enfl., Isofl., Sevofl., Desfl., MV and VT, respiration frequency, peak pressure, plateau pressure, average pressure, PEEP. PARAMAGNETIC GAS MEASUREMENT, Dimensions: H 137 x W 80 x D 80, weight 115 kg. incl. DRÄGER Infinity Gamma XL, 2 DRÄGER Vapor 2000 (Isoflurane and Sevoflurane), used		12.900,00 EUR
1003876	DRÄGER Fabius Tiro. Mobile anaesthesia device. Supplied with compact circuit system, colour display, ventilation modes: Volume controll, Pressure Support, MAN, SPONT, incl. DRÄGER Infinity Delta (EKG, SPO <sub>2</sub> , NibP, Press, Temp) und DRÄGER SCIO Oxy four (O <sub>2</sub> , CO <sub>2</sub> , N <sub>2</sub> O, Halothan, Sevofluran, Isofluran, Enfluran, Desfluran), used		11.900,00 EUR
1004146	DRÄGER Fabius Tiro. Mobile anaesthesia unit with flow meter for O <sub>2</sub> , N <sub>2</sub> O and Air, delivery incl. compact breathing system. Ventilation modes: IPPV, MAN, SPONT, PCV. used		9.800,00 EUR
1004557	DRÄGER Fabius CE, mobile anaesthesia workstation, with dual vaporisor support, electric ventilation, breathing circuit, integrated monitoring: insp. O <sub>2</sub> conc., breathing rate, tidal volume, AMV, mean and plateau pressure, peak respiratory pressure and PEEP, fresh gas flow, second-hand Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511, inclusiv: DRÄGER PM 8050 respiratory gas monitor, airway pressure, expir. breathing volume, insp. O <sub>2</sub> conc., anaesthesia conc. in the bypass procedure for halothane, enflurane, isoflurane, desflurane, sevoflurane, CO <sub>2</sub> , DRÄGER Vaporiser 2000, Sevoflurane, with DRÄGER fill system, GE Dash 4000 multi-parameter monitor, battery operated, ECG, SPO <sub>2</sub> , temp. NIBP, 2x invasive blood pressure, very good condition, second-hand		10.900,00 EUR
1004189	DRÄGER Fabius CE, mobile anaesthesia workstation, with dual vaporisor support, electric ventilation, breathing circuit, integrated monitoring: insp. O <sub>2</sub> conc., breathing rate, tidal volume, AMV, mean and plateau pressure, peak respiratory pressure and PEEP, fresh gas flow, second-hand Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511		6.900,00 EUR






1003366	<p>DRÄGER Fabius GS, anaesthesia workstation with electronically controlled and operated ventilator, 45-minute battery back-up in the event of general power supply failure. ventilation mode: manual, spontaneous, IPPV, PCV. Dimensions: Width: 89.5cm, x height: 130 cm, x depth: 82 cm, electronic measurement of fresh gas flow with semi-closed COSY breathing system, integrated monitoring: insp. O2 conc., breathing rate, tidal volume, AMV, mean and plateau pressure, peak respiratory pressure and PEEP, fresh gas flow, In addition this Fabius GS is provided with the following monitoring: Dräger Infinity Gamma XL with the following functions: ECG, NIBP, SpO2, CO2, NIBP, PRESS, Temp and 5 anaesthesia gases (halothane, enflurane, isoflurane, desflurane and sevoflurane), inclusiv DRÄGER Vaporiser 2000, Sevoflurane, with DRÄGER fill system, DRÄGER support, DRÄGER Vapor 2000, Isoflurane, with DRÄGER-FILL filling unit, used</p>		11.900,00 EUR
1004123	<p>Standard rail (50 cm) with wall mounting, NEW, for mounting DRÄGER Titus etc.</p>		198,00 EUR
1000735	<p>DRÄGER Sulla 909V anaesthesia unit, blue model, Ventilog 2, ORC, low-flow, chassis with drawer and writing tablet, second-hand</p>		4.200,00 EUR
1001450	<p>DRÄGER Sulla 808V anaesthesia machine, blue model, on original chassis with drawer unit, writing surface and bottle support, incl. Ventilog 2, second-hand The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		3.980,00 EUR
1005055	<p>DRÄGER Sulla 808 V mobile anesthesia machine, without drawers, with ventilator Ventilog 2 used, surcharge for bottle holder (10l) (Art. No. 1004512 and 1004511), used The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs are EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		2.680,00 EUR
1004686	<p>DRÄGER Sulla 800 mobile anaesthesia device, with Ventilog 2 ventilation device, vaporizer 2000 isoflurane, and circuit component 8 ISO, second-hand. Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		2.680,00 EUR










1001454	DRÄGER Sulla 800 mobile anaesthesia unit, second-hand. Surcharge for cylinder holder (10l) (art. no. 1004512 and 1004511) The standard inlet pressure for anaesthesia units is approx. 5 bar. If you wish to operate your unit with an oxygen concentrator, we have to convert the unit to concentrator operation. The costs are EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		2.100,00 EUR
1001449	DRÄGER Sulla 808 V mobile anaesthesia device with drawer, with Ventillog 2 ventilation device, second-hand. Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511		2.950,00 EUR
1001451	DRÄGER Sulla 808 mobile anaesthesia machine, second-hand. Surcharge for gas cylinder (10 l) mount, art. no. 1004512 and 1004511 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		2.300,00 EUR
1001453	DRÄGER Sulla 800 V anaesthesia device, Ventillog, trolley/chassis with writing surface, second-hand. Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		2.450,00 EUR
1001468	DRÄGER Tiberius, simple, robust anaesthesia device with drawer unit, O2 flush, second-hand Surcharge for gas cylinder (10 l) mount, art. 1004512 or 1004511) The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.750,00 EUR
1001459	DRÄGER Trajan 800 anaesthesia device, space-saving anaesthesia device for wall mounting (fixed or adjustable by means of quick support), second-hand. - This item is intended solely for veterinary medicine and export -		1.850,00 EUR
1001458	DRÄGER Trajan anaesthesia device, space-saving anaesthesia device for wall mounting (fixed or adjustable by means of quick support), second-hand. - This item is intended solely for veterinary medicine and export - The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.550,00 EUR

- |         |  |   |               |
|---------|--|---|---------------|
| 1005058 | <p>DRÄGER Fabian, wall anesthesia machine, incl. Circular part with inhalation and expiration valve, soda lime container, OHMEDA TEC 5 VET isoflurane, pull-out writing surface, used. This article is intended for veterinary medicine and export only The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs are EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>   |    | 2.100,00 EUR  |
| 1003852 | <p>Dr. MÜLLER anesthesia unit for horses, COMPLETE PRICE, our anesthesia unit with ventilation consists of the following components: DRÄGER TRAJAN 808 Anesthesia unit with flow meter for O<sub>2</sub>, N<sub>2</sub>O, suitable for central gas supply and/or bottle operation with O<sub>2</sub> Flush (240L/min.), with interchangeable Vapor mount, on braked trolley with bottle holder (please state if required!), used DRÄGER vapor 19.3 for Isoflurane, DRÄGER CIRCLE SYSTEM for foals and small animals, including soda lime container, inspiration and exhalation valve, used; LARGE ANIMAL CIRCLE SYSTEM (material: stainless steel, chrome-plated copper), complete with soda lime container, pressure releaf valve, hoses, Y-piece, 4 different adapters for tracheal tubes DR. MÜLLER horse ventilator, pressure-controlled ventilator for horses based on the bag-in-bottle system,with modified BIRD MARK 7 ventilator (sec. hand): output 450 l/min, required compressor output min. 500 l/min (ventilation unit that generates pressure and vacuum in an acrylic glass cylinder). Manual, assisted or controlled ventilation. The device is mobile, extremely robust (4-wheel aluminum trolley with acrylic glass bin), can be easily dismantled, -intended solely for veterinary medicine - information: MyClean® FD: Disinfectant for Dr. Müller horse ventilator (100g contain: 2,5 g Polyhexamethylenbiguanidhydrochlorid, 0,8 g Didecylmethylammoniumchlorid) We recommend an annual service of the complete machine. This concerns in particular the BIRD Mark 7 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p> |    | 11.590,00 EUR |
| 1004896 | <p>DR. MUELLER Laboratory anaesthesia device: Portable table anaesthesia device mounted on a stainless steel carrying frame (WxHxD: 35x31x17cm). Delivery with flowmeter for O<sub>2</sub> and Air (each 0,1 litre -1litre), 2-way valve, overpressure valve, 2 fresh gas hoses 2m with quick coupling and DRÄGER Vapor 2000 for isoflurane (for filling a filling support is necessary, price 89,00 EURO), used The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs are EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>  |  | 2.990,00 EUR  |

1004573	Dr. Müller laboratory animal anaesthesia device: Portable table-top anaesthesia device mounted on a stainless steel supporting frame (WxHxD: 35x31x17 cm). Supplied with Vet. med vapour for isoflurane, NEW (a filling piece is required for filling, price EURO 89.00) flow meter (0.1 litres - 1 litre), two-way valve, pressure relief valve, used The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		2.340,00 EUR
1004173	Dr. Müller laboratory animal anaesthesia device: Portable table-top anaesthesia device mounted on a stainless steel supporting frame (WxHxD: 35x31x17 cm). Supplied with Vet. med vapour for isoflurane, NEW (a filling piece is required for filling, price EURO 93.00) flow meter (0.1 litres - 1 litre), new The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.640,00 EUR
1003830	DR. MÜLLER portable anesthesia unit with flow meter for O2, for hanging in standard rails or placement on the tabletop, incl. stainless steel carrying frame, used. This item is exclusively intended for veterinary medicine and export The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		2.040,00 EUR
1003297	Dr. Müller wall anesthesia unit with flow meter for O2 0,8l-6l, isoflurane vapor (new) and circle part with scrubber tank, in- and ex-spriation valve, used The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.680,00 EUR
1003825	Univet Porta, portable laboratory animal anesthesia unit. Flow meter for O2 and nitrous oxide. Ohmeda TEC 5 Sevofluranvapor, used, This item is exclusively intended for veterinary use The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.650,00 EUR



1001473	<p>DRÄGER Romulus 19 mobile anaesthesia device with drawers and writing surface, second-hand. Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511,- This item is intended solely for veterinary medicine and export- The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		1.850,00 EUR
1004873	<p>DRÄGER Romulus 19 mobile anesthesia unit with drawers and writing tablet, second-hand, surcharge for bottle holder (10l) (Art. No. 1004512 and 1004511), No high-pressure connection (supply only via pressure-reducer or central gas supply) No O2 deficiency signal, no nitrous oxide barrier This item is exclusively designed for veterinary medicine and export- The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs are EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		1.600,00 EUR
1001474	<p>DRÄGER Tiberius 19 mobile anaesthesia device with drawers and writing surface, second-hand. - This item is intended solely for veterinary medicine and export- Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		1.650,00 EUR
1001472	<p>DRÄGER Sulla 19 mobile anaesthesia device, second-hand. - This item is intended solely for veterinary medicine and export- Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		1.790,00 EUR
1002554	<p>DRÄGER Fabius: Mobile anaesthesia machine with FlowMeter for O2, Air and N2O, double holders for Vapor, drawer unit, second-hand Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.</p>		2.100,00 EUR












1001466	DRÄGER Remus A mobile anesthesia machine with double vapor holder, used. Surcharge for bottle holder (10l) (Item No. 1004512 and 1004511) The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		2.100,00 EUR
1003988	Heyer Kemperhof with ABV-A respirator on Siemens trolley for operation with Selectatec Vapor, used. This item is exclusively intended for veterinary use and export The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.900,00 EUR
1003897	HOYER MCM 590: anesthesia unit with integrated ventilation and DATEC Cardiacap 2 monitor, used The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.980,00 EUR
1003770	DRÄGER swiveling drawer block, used		180,00 EUR
1004291	Support frame (stainless steel) for DRÄGER Trajan. This support frame is used to produce a portable anaesthesia device from a Trajan device, new		280,00 EUR
1003831	Fitting for hanging the wall-mounted anesthesia unit on a standard rail, used		75,00 EUR
1001464	Quick-mount bracket for Trajan anaesthesia device (all models) to be attached to a standard rail, new		198,00 EUR
1001463	Quick-mount wall bracket for Trajan anaesthesia device (all models), new - This item is exclusively designed for being used in the field of veterinary medicine or for export purposes -		76,00 EUR
1001465	Support frame (stainless steel) for wall anaesthesia devices. This support frame is used to produce a DRÄGER portable anaesthesia device from a wall-mounted anaesthesia device, new		230,00 EUR

























1004450	VÖLKER mobile anaesthesia unit with isoflurane vapour, flowmeter and circular section, incl. trolley and bottle holder, second-hand The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.		1.390,00 EUR
1004600	DRÄGER 8601765-06: Netzteil für DRÄGER Julian, Uin: 90-264 VAC, Pout: 600 Watt, gebraucht		1.350,00 EUR
1000757	Complete check of your anesthesia device incl. Vapor without electronic monitoring (plus repairs and spare parts according to cost estimate) - This item is exclusively intended for veterinary use		209,00 EUR
1003911	In-Home Tutoring - Anesthesia training in your veterinary practice/clinic by Dr. med. vet. Christine Lendl, CertVA: Veterinary specialist for anesthesia for zoo, penned and wild animals. Switching to inhalation anesthesia - a new anesthesia unit - extension and intensification of your surgical practice equipment? We offer you further training tailored to your needs in your own rooms by a well-known veterinary specialist for anesthesiology! In a half-day training session in your practice or clinic you will learn the fundamentals of anesthesia and analgesics according to your requirements, their combination with anesthesia regimes, the changeover from injection to inhalation anesthesia, the technical principles of inhalation anesthesia, selection and use of the various patient systems, performance and monitoring of anesthesia as well as the termination of inhalation anesthesia. The training can be supplemented with ventilation, care of equipment and peri-operative medicine. We will be glad to assist you and		400,00 EUR
1004344	NARKOVET CONSULTING® offers advice & training on the acquisition of new equipment, the development of a competent staff and the preparation of practice-tested work protocols in the field of anaesthesia, analgesia and perioperative intensive care. She trains veterinarians and veterinary assistants in the proper operation and use of modern anaesthesia, monitoring and ventilation equipment as well as the associated monitoring. For multi-module training courses, NARKOVET CONSULTING® works together with the European School for Advanced Veterinary Studies (ESAVS). Prices and information on request		0,01 EUR
<b>Oxygen bubblers</b>			
1001614	DRÄGER AIR-bubbler (to moisten the airways), without pressure reducer (for wall connection) with integrated flow meter, second-hand		120,00 EUR
1002984	DRÄGER O2 pressure reducer with three outflows: Flow meter with Perlator, 6 bar. Ejector, used		280,00 EUR












**Vapors  
(anaesthetic  
vaporisers)**

1003258	DRÄGER O2 pressure reducer with FlowMeter and humidifier, used		98,00 EUR
1000775	DRÄGER Vapor 19.1 rededicated for isoflurane -with open filling device (no safety filler neck is required for these vapors), used, -this item is intended exclusively for veterinary medicine-		580,00 EUR
1005088	A.E. Vapor Isoflurane VET. Anaesthetic vaporiser for isoflurane with cagemount connection. A safety filler neck is required for filling-This article is intended exclusively for veterinary medicine- new		980,00 EUR
1001482	DRÄGER Vaporiser 19.3 for Isoflurane (original DRÄGER Isofluran Vaporiser). second-hand This article is intended exclusively for veterinary medicine and export		950,00 EUR
1001488	DRÄGER Vaporiser 19.3, Sevoflurane, with quick fill system, second-hand This item is exclusively designed for being used in the field of veterinary medicine or for export purposes		980,00 EUR
1004089	DRÄGER Vapor 2000, Isoflurane, with DRÄGER-FILL filling unit used		1.250,00 EUR
1001483	DRÄGER Vaporiser 2000, Isoflurane, with safety filling device, second-hand		1.100,00 EUR
1004663	DRÄGER Vapor 2000, Isoflurane, with funnel fill system (open filling device), used		790,00 EUR
1005110	NorVap Jupiter Sevoflurane Vapor with Selectatec connection and Quick-Fill filling device, as new. This vapor is intended exclusively for veterinary medicine or export.		780,00 EUR
1004687	DRÄGER Vapor 2000, halothane, with funnel fill system -open filling device-, used		690,00 EUR
1001491	DRÄGER Devapor, vaporiser for desflurane, with integrated heating, used		890,00 EUR









1000774	DRÄGER Vapor for Desflurane, with integrated heating, used		980,00 EUR
1001492	Vaporiser refitting 19.1/19.for isoflurane, check your vaporiser subsequently refitted for isoflurane, second-hand. - This item is intended solely for veterinary medicine -		298,00 EUR
1001493	Technical inspection of vaporiser - all manufacturers (halothane, enflurane, isoflurane, desflurane), test for leaks and correct output		127,00 EUR
1004640	Technical inspection of vaporiser - all manufacturers (sevoflurane), test for leaks and correct output		155,00 EUR
1003260	DRÄGER sealing plate for vaporiser filling device, used - This item is suitable for human medicine -		129,00 EUR
1001494	Sealing plate for DRAEGER vaporiser filling device, new - This item is exclusively designed for being used in the field of veterinary medicine or for export purposes -		129,00 EUR
1002405	Vapor-Stand (fixed): This item acts as a stand for a DRÄGER Vapor and is to be operated seperately from the anaesthesia device, NEW -This item is intended solely for veterinary medicine-		360,00 EUR
1004556	DRÄGER Vapor adapter (plug-in system), AUTO EXCLUSION SYSTEM, NEW		455,00 EUR
1002270	DRÄGER (fixed) vaporiser support for two vaporisers (switchable), second-hand		180,00 EUR
1001497	(Fixed) vaporiser support, new - This item is intended solely for veterinary medicine -		298,00 EUR
1005212	Vapor holder tripod for Dräger Vapore 19.1 and 19.3. With this holder a Dräger Vapor can be attached to a round tube such as an infusion pole. This product is not an original Dräger product and may therefore only be used in veterinary medicine, used		375,00 EUR

1000795	DRÄGER Vapor-Parkhalterung, zur Befestigung an einem Narkosegerät oder der Montage an der Wand, gebraucht.		49,00 EUR
1001499	DRÄGER vaporiser park support for 2 vaporizer, second-hand		49,00 EUR
1004563	SELECTATEC holder for DATEX TEC 3, 4, 5 and 7, NEW		198,00 EUR
1001515	SELECTATEC adapter unit. This adapter can be used to adapt DRÄGER vaporisers to TEC systems, new. - This item is intended solely for veterinary medicine -		179,00 EUR
1005087	Cagemount adapter stainless steel, male/female for mounting on TEC 3 vaporisers, hose connection 6mm, new		89,00 EUR
1001506	OHMEDA Tec 5 vaporiser for isoflurane, second-hand		680,00 EUR
1004814	OHMEDA Vapor Tec 5, for isoflurane (converted enflurane vapor), second-hand This item is exclusively designed for veterinary medicine and export-		580,00 EUR
1001507	OHMEDA Tec 5 vaporiser for sevoflurane with quick fill system, second-hand		680,00 EUR
1000783	OHMEDA Vapor Tec 5, for Halothane, used		680,00 EUR
1003333	DATEX TEC 7: Vapour for halothane with quick fill device, used		780,00 EUR
1001509	OHMEDA Tec 7 vaporiser for enflurane with quick fill system, in mint condition		680,00 EUR












1002806	OHMEDA Vapor Tec 7, for Sevoflurane, with quick closure, used		980,00 EUR
1004460	OHMEDA Vapor Tec 7, for Isoflurane, with quick closure, good condition		980,00 EUR
1003884	PENLON Sigma Delta for Sevoflurane, with Selectatec mount, used		620,00 EUR
1001512	SIEMENS VAPORISER for isoflurane, second-hand		380,00 EUR
1003161	SIEMENS VAPORISER (for KION) for isoflurane, used		580,00 EUR
1003321	SIEMENS VAPOUR (for KION) for sevoflurane, used		780,00 EUR
1002759	BLEASE Sevoflurane Vaporiser with Quick Fill device, for Selectatec-Support, second-hand		490,00 EUR
1004009	DATEX Aladin Vapor for Sevoflurane, used		1.200,00 EUR
1004514	DATEX Aladin Vapor for Isoflurane, used		1.200,00 EUR
1004515	DATEX Aladin Vapor for Enflurane, used		1.200,00 EUR
1001516	Safety filling supports (plastic) for isoflurane, new		93,00 EUR












1001517	DRÄGER safety filling supports for isoflurane, material: metal, used		175,00 EUR
1004278	DATEX-Ohmeda Bottle Adapter Isoflurane, new		95,00 EUR
1004763	DATEX-Ohmeda Bottle Adapter Sevoflurane, new		95,00 EUR
1004764	Datex-Ohmeda Bottle Adapter Halothane, new		95,00 EUR
1003434	DRÄGER fill adapter isoflurane quick fill, NEW		75,00 EUR
1004654	DRÄGER fill adapter sevoflurane DRÄGER fill, NEW		75,00 EUR
1001519	DRÄGER safety filling supports for halothane, material: metal, second-hand		127,00 EUR
1004499	DRÄGER safety filling supports for halothane, material: metal, new		160,00 EUR
1001518	Safety filling supports (plastic) for halothane, second-hand		45,00 EUR
1004498	Safety filling supports (plastic) for halothane, new		89,00 EUR
1001520	Safety filling supports (plastic) for enflurane, used		45,00 EUR

**Pressure reducers  
und central gas  
supply**












1001521	DRÄGER safety filling supports for enflurane, material: metal, second-hand		127,00 EUR
1001522	Safety filling supports (plastic) for sevoflurane, with special adapter for Abbott bottle with quick fill system, new		125,00 EUR
1003755	ABBOT Sevoflurane discharge adapter for Vapor, NEW		69,00 EUR
1001570	DRÄGER OMS 800 oxygen deficiency signal with nitrous gas (retrofit kit), second-hand. - This item is intended solely for veterinary medicine and export-		190,00 EUR
1004362	GCE Mediselect O2 pressure reducer, adjustable flow from 0.05 L / min - 1L / min, new		59,00 EUR
1001571	GREGGERSEN pressure reducer for oxygen with quick coupler, new		283,00 EUR
1004878	HIGH-PRESSURE SEAL O2, sealing ring for oxygen pressure reducer (This seal is used at the connection to the oxygen cylinder), NEW		7,50 EUR
1004947	DRÄGER pressure reducer O2 pin index for DRÄGER Primus and Pallas, second-hand		450,00 EUR
1004445	DRÄGER O2 switch NIST with O2 quick coupling, second-hand		98,00 EUR
1003892	Pressure reducer O2 Pin Index, British system, NIST, used		145,00 EUR
1004946	DRÄGER pressure reducer N2O, pin index, with quick coupling for the British gas system, second hand		145,00 EUR











1005047	DRÄGER pressure reducer for oxygen with double oxygen withdrawal box, ejector for emergency suction, used		150,00 EUR
1001572	DRÄGER pressure reducer for oxygen, second-hand		127,00 EUR
1002134	DRÄGER Oxett: Support frame for a small O2 bottle (2 litres) with 3-fold distribution box for O2, second-hand		220,00 EUR
1001575	DRÄGER pressure reducer for nitrous gas (standard in the field of human medicine), second-hand		127,00 EUR
1003891	Pressure reducer N2O, Pin Index, with hose 1m, NIST, second-hand		98,00 EUR
1004439	DRÄGER pressure reducer N2O, pin index, NIST, second-hand		150,00 EUR
1004215	DRÄGER 8603462-01 Pressure reducer N2O Pinindex for DRÄGER Primus and Palla, second-hand		180,00 EUR
1004381	DRÄGER 8603463 N2O pressure reducer for DRÄGER Primus, as new, second-hand		280,00 EUR
1004359	Pressure reducer f. CO2 German standard, NEW The following accessories can be obtained from OLYMPUS: OLYMPUS MAJ-590: Insufflation hose OLYMPUS -591: Suction hose		520,00 EUR
1004553	DRÄGER quick release coupling (female) for O2, used		98,00 EUR
1003867	DRÄGER quick release coupling (female) for N2O, used		98,00 EUR

1002439	DRÄGER Pressure Supply Display for O <sub>2</sub> ,N <sub>2</sub> O and Air incl. supply hoses, second-hand		180,00 EUR
1001576	GREGGERSEN pressure reducer for compressed air with quick coupler (standard in the field of human medicine), new		283,00 EUR
1004440	DRÄGER pressure reducer for compressed air with NIST, used		198,00 EUR
1003887	Central O <sub>2</sub> gas supply hose, British system, with quick release coupling, used		56,00 EUR
1003889	Central gas supply N <sub>2</sub> O hose, British system, NIST, used		25,00 EUR
1003893	Quick release O <sub>2</sub> coupling, double, with Y-piece, British system, used		98,00 EUR
1004512	Surcharge for bottle holder on anaesthesia unit for 1 pc. O <sub>2</sub> bottles (10l), NEW		198,00 EUR
1004511	Surcharge for bottle holder on anaesthesia unit for 2 pcs. O <sub>2</sub> bottles (10l), NEW		248,00 EUR
1001589	Support for anaesthesia bottle (O <sub>2</sub> and N <sub>2</sub> O): Stainless steel cross for firm attachment of bottle base, NEW		85,00 EUR
1001590	Bracket for anaesthesia bottle (O <sub>2</sub> and N <sub>2</sub> O): Support fork with adjustable strap to attach the upper part of the bottle to the anaesthesia device or to the wall, etc. NEW		85,00 EUR
1003623	DRÄGER gas bottle wrench, used		12,00 EUR

## Flowmeters










1004465	DRÄGER A-cone switch, second-hand		480,00 EUR
1005149	DRÄGER O2 flowmeter, 2-16 litres/min., for rail mounting, outlet plug, connection with angled plug for ZGV socket, second-hand.		190,00 EUR
1004496	DRÄGER Flow meter O2 for standard rail with hose and, sec. hand		98,00 EUR
1004830	DRÄGER tube flowmeter O2, rail unit, with tube and O2 ZGV connection, incl. nebuliser bottle		184,00 EUR
1005051	DRÄGER O2 flow meter double for operation on an oxygen withdrawal can, each 1-15 liters / min, used		160,00 EUR
1003840	In-line air flow meter, rail mounted unit, NIST connection, NEW		198,00 EUR
1004995	GREGGERSEN Kolibri, O2 rail flowmeter with NIST connection, adjustable from 0.1L-6L/min, new		180,00 EUR
1004552	DRÄGER O2 distributor with pressure reducer, flowmeter (1 to 15l/min), direct outlet and ejector, second-hand		280,00 EUR
1004158	DRÄGER Flowmeter for air, with quick coupling for central gas supply, second-hand		75,00 EUR
1005044	DRÄGER pressure reducers for oxygen and tube flow meters. Adjustable flow: 1 liter-15 liters, used		130,00 EUR
1004881	Dwyer RMA Series Flowmeter Air, 4-fold 1-5 litres/min, second-hand		150,00 EUR

**Compressors,  
concentrators**


1004619	Ohmeda Quantiflex: flowmeter for O2 and N2O with O2 flush, used		380,00 EUR
1004509	Conversion of your anesthetic machine to operation with an O2 concentrator (0,59 bar bei 6l/ min). After this transformation the operation with gas cylinders or central gas supply will be possible but the O2 lack signal will be out of order.		89,00 EUR
1003618	DeVilbiss oxygen concentrator, output 5L/min, dimensions: HxWxD: 62.2x34.2x30.4 cm, weight: 16.3kg, noise level: 40 dbA (50Hz), flow rate (continuously variable): 0.5-5L/min, oxygen concentration: 93% +/- 3%(0.5-5L/min), operating voltage: 230V/50Hz, power consumption: with 2L approx. 290 Watt / with 5L approx. 312 Watt, outlet pressure 0.59bar, flow rate 5 liter with 93Vol. %, ACCESSORIES (not to be used in connection with an anesthesia machine!), including bubble humidifier and 15 m connection hose, NEW		790,00 EUR
1004313	DeVilbiss oxygen concentrator, output 10L/min, dimensions: HxWxD: 62.2x34.2x30.4 cm, weight: 16.3kg, noise level: 40 dbA (50Hz), flow rate (continuously variable): 0.5-5L/min, oxygen concentration: 93% +/- 3%(0.5-5L/min), operating voltage: 230V/50Hz, power consumption: with 2L approx. 290 Watt / with 5L approx. 312 Watt, outlet pressure 0.59bar, flow rate 5 liter with 93Vol. %, ACCESSORIES (not to be used in connection with an anaesthesia machine!) including bubble humidifier and 15 m connection hose, NEW		1.490,00 EUR
1004811	DeVilbiss connection fitting for oxygen for O2 concentrator DeVilbiss, plastic, new		9,80 EUR
1004403	SILENT compressor for compressed air, operating pressure: 8 bar, boiler capacity 100 l, weight 75 kg, suction capacity: 370 l/min, OIL-FREE, NEW Use exclusively in ANIMAL MEDICINE		745,00 EUR
1004866	SILENT compressor for compressed air, for operating veterinary anaesthesia equipment with ventilation (e.g. Dräger Sulla 808V, Sulla 800V). Oil-free compressor, whisper-quiet, pressure tank 40 litres, motor power 750W. 2 pressure gauges (for container and working pressure). Delivery incl. pressure reducer, quick coupling and water separator, NEW		450,00 EUR
1005273	PLANET-AIR L-S50-15, whisper-quiet compressor with 15-litre tank, suction capacity 50 l/min, pressure reducer and water separator, second-hand -This item is intended exclusively for veterinary medicine and export.		650,00 EUR

**OP and  
investigation  
lights**











**Bulbs**

1001401	Illuminant for HANAULUX Hamburg, 24V / 40W, new		18,00 EUR
1001402	Halogen illuminant for HANAULUX Oslo and Heidelberg, 22,8V40W, new		15,50 EUR
1001403	Halogen illuminant for HANAULUX London, 22,8V50W, new		16,90 EUR
1001404	Halogen illuminant for HANAULUX Universal and London, old model, 24V/50W, new		24,90 EUR
1001405	Halogen illuminant for HANAU Hamburg, 24V / 40W, new		33,00 EUR
1001406	Illuminant for HANAULUX Boston and Universal, 24V50W, new		14,90 EUR
1001407	Illuminant for HANAULUX Univeral, 230V50W, new		24,90 EUR
1001408	Halogen illuminant for HANAULUX Amsterdam, 22.8V75W, new		27,90 EUR
1001409	Halogen illuminant for HANAULUX 2000 series, 22.8V50W, new		16,00 EUR
1001410	Halogen illuminant for MARTIN Chromophare C570D, 24V150W, new		20,90 EUR












**OP lights**

1002419	Berchtold Chromophare D650 (22.8 V 150 W) with satellite D500 (22.8 V 150 W). Field adjustment 17-28 cm electrical brightness adjustment, suitable for room heights over 3 metres, used		4.900,00 EUR
---------	---	---	--------------











1005202	BERCHTOLD Chromophare D530plus. Mobile operating light. 22,8V, 150W halogen, field adjustment 17 - 28cm, electronic brightness control, approx. 80.000Lux, diameter of lamp body 52cm, second-hand.		1.850,00 EUR
1004318	DRÄGER Sola Premium Duo, consisting of DRÄGER Sola 500 (250 Watt) and DRÄGER Sola 700 (300 Watt). From room height 330 cm, The lamps are NEW and in their original packaging.		4.500,00 EUR
1000657	DR. MACH Triaflex operating light, 3 x 40 Watt, wall-mounted, focusable, 2 brightness levels, incl. sterilisable handle sleeve, second hand		980,00 EUR
1001380	HANAULUX Hamburg OR lamp, 4x40W halogen, wall mount or ceiling mount (please specify together with ceiling height when ordering!), second-hand		890,00 EUR
1000656	DR. MACH Trigenflex R96, 3 x 50 watts, approx. 75,000 lux, focusable. Diameter of lamp body 48cm, light field: 15 - 26 cm, ceiling or wall mounting, with sterilisable handle sleeve, incl. transformer, second-hand		1.390,00 EUR
1001386	Martin ML 701 surgical light. 24V/150Watt , approx. 80,000Lux. Focusable. Incl. sterilisable handle and transformer, good condition, used. For ceiling mounting from a ceiling height of 260cm (please specify ceiling height when ordering).		1.850,00 EUR
1001375	HANAULUX Oslo operating theatre lamp, 4 x 40 W halogen, wall or ceiling installation (please state with ceiling height when ordering), focusable, ceiling height from 260 cm, second-hand		1.190,00 EUR
1001376	HANAULUX Oslo Duo, 4 x 40 watt in each, ceiling installation, focusable, ceiling height from 300 cm, second-hand		2.100,00 EUR
1001377	HANAULUX London operating theatre lamp, 3 x 50 W halogen, ca. 75,000 Lux, wall or ceiling installation (please state with ceiling height when ordering), ceiling height from 260 cm, second-hand		1.080,00 EUR
1001378	HANAULUX London Duo operating theatre lamp, 2 London light fixtures (3 x 50 W halogen, ca. 75,000 Lux), with general ceiling installation, ceiling height from 300 cm, second-hand		1.980,00 EUR
1000635	HANAULUX 2003 operating theatre lamp, 3 x 50 W, ca. 75,000 Lux, wall or ceiling installation (please state with ceiling height when ordering), ceiling height from 260 cm, second-hand		1.650,00 EUR

1000636	HANAULUX 2003 duo operating theatre lamp for ceiling installation (minimum height of 3.00 m). 2 light fixtures (3 x 50 W halogen, ca. 75,000 Lux), incl. sterilisable grip sleeve and new transformer (state ceiling height when ordering), ceiling height from 300 cm required		3.100,00 EUR
1002365	HANAULUX 2004, four-beam operating light 22.8 V/V50 W (ca. 80,000 Lux), with wall or ceiling fixing, ceiling height from 300 cm, used		1.980,00 EUR
1001383	HANAULUX 2004 Duo, 2 x four-ray operating light 24 V 50 W (2 x ca. 80000 Lux), with joint ceiling fixing, ceiling height from 300 cm, used		3.500,00 EUR
1002735	HANAULUX 2005 duo (2x ca. 75.000 Lux) 22.8V/50W. Ceiling height from 300 cm, second-hand		4.400,00 EUR
1001385	HANAULUX 2007 (ca. 130,000 Lux) operating theatre lamp, with satellite 2003 (ca. 100,000 Lux) with general ceiling installation, ceiling height from 330cm, second-hand.		3.500,00 EUR
1004177	MAQUET blue 130, surgical light for ceiling mounting, 3 x 90 Watt, focusable, ceiling height from 280cm, second-hand		2.980,00 EUR
1005232	DERUNGS Triango 100-3 Duo. Duo operating light for ceiling mounting . Light output: 2x 100000lux/1m, colour temperature adjustable, light field focusable 18-28cm, dimmable, 3 different light colours selectable. Please state the room height when ordering: Distance between floor and concrete ceiling (if you have a false ceiling, please also specify the distance to the concrete ceiling). NEW		5.590,00 EUR
1005230	DERUNGS Triango 100-3. Operating light for ceiling mounting or optionally for wall mounting. Light output: 100000lux/1m, adjustable colour temperature, focusable light field 18-28cm, dimmable, 3 different light colours selectable. When ordering, please state whether you want ceiling or wall installation (for ceiling installation, we also need information on the room height and whether there is a false ceiling). NEW		4.190,00 EUR
1005231	DERUNGS Triango 100-3 F. Operating light on rolling stand. Light output: 100000lux/1m, adjustable colour temperature, focusable light field 18-28cm, dimmable, 3 different light colours selectable. NEW		4.530,00 EUR
1004625	Ceiling canopy, diameter 60 cm, height 15.5 cm, diameter of the opening 14.0 cm, NEW; With retaining ring for fixing the canopy. Please specify the diameter of your ceiling flange when ordering!		360,00 EUR



















1004626	Ceiling canopy, diameter 45 cm, H 19 (9) cm, diameter of the opening 5,3 cm, NEW; With retaining ring for fixing the canopy. Please specify the diameter of your ceiling flange when ordering!		240,00 EUR
1001390	Ceiling anchorage ring for all operating theatre lamps, new incl. nuts and drilling template		82,00 EUR
1001393	HANAULUX sterilisable handle for 2000 series operating theatre lamps, NEW		68,00 EUR
1001395	HANAULUX sterilisable handle for 2007 operating theatre lamps, made from metal, used		89,00 EUR
1001394	MARTIN Chromophare sterilisable handle for operating theatre lamps, NEW		54,00 EUR
1004517	Sterilisable handle for HANAULUX London and Oslo operating theatre lamps (aluminium), NEW This item is exclusively designed for being used in the field of veterinary medicine or for export purposes		89,00 EUR
1003528	MARTIN sterilizable handle for Martin ML series, new		63,00 EUR
1004771	MARTIN sterilisable handle for Martin ML series, second-hand.		30,00 EUR
1005233	DERUNGS sterilisable handle for Triango 100 and Triango 100 duo operating lights, new		38,00 EUR
<b>Transformers and UPS</b>	1004265 Transformer 24 V DC 751W, for OP luminaires, NEW Only for mounting under the ceiling canopy! A housing is required for surface mounting.		310,00 EUR
	1004264 Transformer 24 V DC 600W, for OP luminaires, NEW Only for mounting under the ceiling canopy! A housing is required for surface mounting.		210,00 EUR

## Investigation lights












1004263	Transformer 24 V DC 320W, for OP luminaires, NEW Only for mounting under the ceiling canopy! A housing is required for surface mounting.		120,00 EUR
1001364	HANAULUX Heidelberg examination light, wall or ceiling mounting (please state together with ceiling height when ordering), 40 watt halogen, ca. 25,000 Lux, integrated transformer, used		920,00 EUR
1001367	HANAULUX Blue 30, for ceiling mounting, 50 watt, ca. 30,000 Lux, used		980,00 EUR
1002851	DERUNGS reading light with flexible gooseneck arm to attach to DRÄGER Julian side rail and similar, second-hand		80,00 EUR
1005011	Dr. MACH Soloflex, examination lamp for ceiling or wall mounting (please specify room height when ordering !), 24V, 50 Watt, light field focusable, second-hand		490,00 EUR
1003477	Dr. MACH Soloflex, examination lamp for ceiling or wall mounting (please specify room height when ordering !), 24V, 50 Watt, second-hand		450,00 EUR
1001370	BERCHTOLD Chromophare C200 examination lamp 22.8V / 55W, approx. 25,000Lux for wall or optional ceiling mounting, room height from 250cm, good condition, second-hand		690,00 EUR
1002383	DR. MACH 130, examination light for wall mounting, 22,8V/50 Watt halogen, approx. 30000 Lux, focusable with integrated transformer , second-hand		750,00 EUR
1005234	DERUNGS Visiano 20, examination lamp LED for ceiling mounting. Light output 60.000Lux/0,5m. Light field dimmable in 4 steps, two adjustable colour temperatures (3500K/4500K), light field diameter 21 cm /1meter, NEW		1.720,00 EUR
1005012	HANAULUX universal examination lamp, wall-mounted or ceiling-mounted (please specify together with ceiling height when ordering!), 50 watt, 24 volt, approx. 20,000 lux, second-hand		490,00 EUR












## Operating tables and patient trolleys

<b>Operating tables</b>	1004485	MAQUET Alphamaxx 1133.12B3, mobile operating table SELF DRIVEN; HIGH MECHANICAL LOAD UP TO 450 kg, with one-piece leg plate, all functions electrically operated, including cabled control unit, supplied with new cushions, used		13.800,00 EUR
	1005263	MAQUET 1150 Alphamaquet Plus (1150 02 D0), transportable operating table system with undivided leg plate, cable remote control for adjustment of back plate, height, Trendelenburg / Antitrendelenburg bearing, inclination, including charger, table top suitable for X-rays, holder for X-ray cassette, table top can be moved longitudinally on its base, upholstery used without damage, second-hand		11.900,00 EUR
	1005270	MAQUET 1140.00 Betamaquet, mobile operating table column system (the table column can be transported in the room together with the table top or without it by means of a carriage, as it is not firmly anchored in the floor but can be moved and set down freely), incl. transporter 1140. Separate transport of the table top possible. Electric adjustment of height, lateral tilt left/right, tilt forwards/backwards, battery operation. Manual adjustment of the back plate, split and bendable leg plates, adjustable height: 72-108 cm, width of the table top: 50 cm, dimensions of the base: 50x100 cm, removable leg plates. Operation is possible with mains connection as well as with battery, second-hand		8.900,00 EUR
	1004576	MAQUET Alphastar 1132 01A0, mobile operating table with one-piece leg plate, all functions electrically adjustable: inclination, Trendelenburg, tilt left / right, tilting back plate and leg plate, zero position length 2095 mm, battery and mains operation, control via cable control unit, incl. used upholstery in good condition, second-hand		9.800,00 EUR
	1001737	MAQUET 1140 60A transporter without table plate, good condition, second-hand		980,00 EUR
<b>Accessories for operating tables</b>	1004656	SCHMITZ 101.601 uro-adaptor for extension of the seat section, fully radiolucent, with side rails, second-hand		780,00 EUR
	1005201	SMITH EQ-5000 Level 1 Equator, convective patient warmer. Adjustable temperature range: 36°C, 40°C, 44°C. Delivery incl. 2-metre supply hose, second-hand. We recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: vet@moeckundmoeck.de		590,00 EUR

1005240	COVIDIEN WarmTouch. Patient warming system. Adjustable temperature in 4 steps (32°C, 38°C, 43°C, 45°C Boost). Delivery includes transport frame and air hose, used. We recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: vet@moeckundmoeck.de		740,00 EUR
1005227	SMITHS EQ-5000 Level 1 Equator, convective patient warmer. Adjustable temperature range: 36°C, 40°C, 44°C. Delivery incl. supply hose 2meter, used. Delivery includes 5-foot chassis. We recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: vet@moeckundmoeck.de		670,00 EUR
1005247	SMITH Level 1, convective patient warmer with adjustable temperature levels: 37°C, 40°C, 43°C. Delivery incl. 2 metre supply hose, second-hand. We recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: vet@moeckundmoeck.de		650,00 EUR
1003744	MAQUET 1130.70AC Maquetmatic, without accessories, as illustrated, used		1.600,00 EUR
1005021	MAQUET 1005.01B0, gynaecological leg holder according to Göpel, suitable to pair of leg plates Maquetmatic, second-hand		650,00 EUR
1001671	Anaesthesia arm, length adjustable on both sides, NEW		380,00 EUR
1003765	AMATECH Allen L 5000 surgical unit, 1 piece, used		320,00 EUR
1005195	MAQUET 1130.55A0, table top extension with integral foam padding, suitable for tables of the series Maquet 1130, second-hand		380,00 EUR
1005196	ESCHMANN table top extension for operating tables with lateral standard rails and gynaecological cut-out, second-hand		395,00 EUR

1004754	TRUMPF 4544770 Extension accessories, second-hand		680,00 EUR
1005135	TRUMPF 7019121, undivided leg plate 88x50cm, suitable for Trumpf Titan. On request, a new pad can be delivered for an extra charge, used		590,00 EUR
1005019	TRUMPF fastening element for the standard rail. Rotatable mounting bracket with fixable tube holder, used.		110,00 EUR
1004718	MAQUET 1007.04AC positioning unit for spinal surgery, with transport frame, suitable for Maquet 1150, used		2.980,00 EUR
1001741	MAQUET gynaecology leg support for affixing to the normal rail by radial clamps, 1 pair, used		550,00 EUR
1003364	MAQUET 1003.33 knee supports 1 pair, used		780,00 EUR
1004661	MAQUET 1002.36AO steady arm with upholstery, one pair, used		360,00 EUR
1004963	MAQUET 1338, side support without upholstery, second-hand		98,00 EUR
1003432	MAQUET Arm strap		180,00 EUR
1003077	Radial clamp, 2 pcs. for fixation of Göpel leg supports to the unit rail, NEW		385,00 EUR
1004135	DRÄGER/MAQUET clamp for fixing a unit to the standard rail, second-hand		45,00 EUR

1004136	Screw clamp for fastening a unit to the standard rail, second-hand		39,00 EUR
1004356	TRUMPF 1218832: Pulling spindle unit, second-hand		560,00 EUR
1005020	MAQUET 1002.24 Bracket for handcuff. The actual handcuff is not included in the price, used		75,00 EUR
1005022	MAQUET 1003.49BO, spar extension, used.		165,00 EUR
1005024	MAQUET guide rail 75cm, for mounting attachments or extension accessories, second-hand.		225,00 EUR
1004335	MAQUET plate extension, second-hand		250,00 EUR
1004341	MAQUET 1150.59 Mounting clamp for standard rail, second-hand		190,00 EUR
1004962	MAQUET head calotte with new pad. The head calotte is attached to the mounting joint of the operating table. The pad is attached to the head calotte with a Velcro fastener. Used		395,00 EUR
1004353	MAQUET 1003.50AO Lower leg counter traction bar incl. pad roll, inclinable by means of toothed rim joint, second-hand		490,00 EUR
1004355	MAQUET connecting rod, used to extend the head plate or head dome		120,00 EUR
1004369	MAQUET 1004.86Bo / 07: Connector, second-hand		175,00 EUR

1004370	MAQUET 1150.36Ao / 05: head plate holder, second-hand		980,00 EUR
1004376	MAQUET 1001.97AO / 02: Foot plate incomplete		45,00 EUR
1004377	MAQUET 1004.92A: Connector, second-hand		78,00 EUR
1004378	MAQUET 1003.49AO / 02: Connector, second-hand		78,00 EUR
1004749	MAQUET 1002.18AO. Teleskop-Stütze für Hand- OP-Tisch 100163KO, gebraucht		78,00 EUR
1004750	MAQUET 100.35AO, Rotationszugbügelkloben, gebau cht		128,00 EUR
1004751	MAQUET 1003.5100 Führungsholm mit Universalstütze, sehr guter Zustand, gebraucht		280,00 EUR
1004752	MAQUET 1003.34AO Rotationskippkloben, gebraucht		168,00 EUR
1004753	MAQUET Ballonpumpe, gebraucht		35,00 EUR
1002138	SCHAERER Additional plate for operating table S-200, second-hand		380,00 EUR
1004110	Cushion for hand surgery table, NEW		280,00 EUR



1004419 MAQUET Telemaquet-IR 3110.30B5, used 320,00 EUR



1004929 MAQUET 1009.80.C foot switch, used 450,00 EUR



1002996 MAQUET waste collector 17 x 27 cm with unit rail clamp claw, used 35,00 EUR



1005023 Maquet 1002.5800, slide rail piece with braked ball joint adjustment for head plates or head calottes, second-hand 275,00 EUR



#### Patient trolleys

1004775 MAQUET 4748 Transmobil, patient transporter with lateral patient protection plates (foldable), 4 lockable swivel castors, back plate can be set up, upholstery used in good condition, X-ray capable table top, used 2.690,00 EUR



#### Osteosynthesis

#### X-ray fluoroscopy

**x-ray machines** 1005286 Veterinary X-ray unit with universal X-ray table, stainless steel frame, plastic table top, cassette drawer, longitudinally adjustable stand column. Dimensions: 140x60x64cm Height 207cm Siemens Polymobil II generator, output: 40-90kV (adjustable in 16 steps), 0.32mAs-200mAs (adjustable in 29 steps), second-hand 3.400,00 EUR



**Accessories for x-ray machines** 1001700 PLANILUX iris light 100, used for the accurate detection of minor contrasts, for the improved assessment of over-exposed X ray images with high light density, max. 70.000 CD/qm, continuous dimming, precision iris diaphragm of 5 - 100 mm in diameter, continuous collimation to the smallest image areas, in mint condition 180,00 EUR











1005274 X-ray fixation cot, V-shaped, adjustable positioning cot for small animals made of wood, second-hand -This article is intended exclusively for veterinary medicine and export- 75,00 EUR



1005275 LUMAX X-ray protective gloves, pair, Pb-equivalent 0.5, universal size, flexible material, second-hand 110,00 EUR





	1005277	LUMIX X-ray cassette 30x40cm with green-sensitive reinforcement film 400, second-h and		49,00 EUR
	1005278	LUMAX X-ray cassette 24x30cm with green-sensitive reinforcement foil 400, second-h and		40,00 EUR
	1005279	LUMAX X-ray cassette 18x24cm with green-sensitive reinforcement foil 400, second-hand		35,00 EUR
	1005283	Wall holder for X-ray aprons, foldable, 26x15cm, second-hand		49,00 EUR
<b>Mobile c-arms</b>	1004502	Siemens Siremobil ISO-C, mobile C- arm. Fluoroscopic output 40KV to 110KV 0.2bis 8.9mA. Image intensifier 23cm. , Resolution 1.4LP/mm CRT display, Memoscope C- Mod SUB, Subtraction, Roadmap, Landmark, 900 Image Speicher. second-hand		12.800,00 EUR
	1003589	OEC 7700: Mobile C-arm with double monitor, 23 cm BV, digital image storage, incl. black and white video printer, used		11.800,00 EUR
<b>Accessories for mobile c-arms</b>	1004305	Test specimen for constancy test, NEW		980,00 EUR
	1004456	IBA measuring case IBArad/flu-A: test case Primus A, dosimeter Dosimax plus with semiconductor detector DEDX IBA measuring case IBAcn incl. photometric measuring head LX-LS, NEW		5.500,00 EUR
<b>Other</b>	1002744	DeSoutter CCS Clean Cast System: Plaster Cast Saw with electric suction device (mobile), used.		980,00 EUR
	1004957	Instrument container made of stainless steel, with lid 60 x 23 x 7 cm, used		25,00 EUR
	1004696	HAEBERLE 16580 Kombi-Set, grey-white, to mount on dispenser, NEW		28,00 EUR

1004697	HAEBERLE picbox 5-fold, NEW		89,00 EUR
1004698	HAEBERLE picbox 3-fold, NEW		99,00 EUR
1004588	VBM Tourniquet 2 x 500: Operating pressure: 2-10 bar, two outlets, with tripod, basket and 4 m pressure hose and wall connector, without cuffs, used		570,00 EUR
1005264	MarMed operating and treatment table with pump-hydraulic height adjustment 71-111cm, stainless steel top 130x60cm with beaded edge and drain, 8 fastening hooks under the table top, lifting force approx. 100kg, top can be tilted almost vertically, second-hand		790,00 EUR
1005281	SOEHNLE Professional 8320, battery-operated small animal scale, graduation in 10g steps, weighing range up to max. 15kg. Dimensions: 74x5.3x31.5cm(WxHxD), foldable, second-hand		69,00 EUR

## Ultra sound devices



### Ultrasound systems












1004970	GE IC5-9 vaginal probe for GE Voluson 730 Pro, good condition but the rubber sheath is discolored, used		1.490,00 EUR
1005197	BARD BARDSCAN IIS, portable ultrasound scanner on trolley, with fold-out touch screen, ultrasound probe 3.5/5.0MHz, including printer. Mains and battery operation. 1 battery is included in the scope of delivery, second-hand		490,00 EUR




## Ultrasonic cleaners

### Ultrasonic nebulizers

### Investigation

1004691	KaWe Trömmel reflex hammer, heavy, NEW		19,00 EUR
1004721	SOLENI Classic Comfort 14: mobile massage and treatment table, electric height and tilt adjustment, second-hand		1.280,00 EUR

	1004774	Digital clinical thermometer with extended measuring range 28.0 to 42.9°C, NEW 10 pcs.		29,00 EUR
<b>Video printer and recorder</b>	1002128	Sony CVPM1E. Colour video printer, second-hand, good condition.		290,00 EUR
	1003275	Sony colour printing pack UPC-2010, colour ribbon and paper for 200 printouts, compatible with SONY printers 2000 series		198,00 EUR
	1004413	SONY UP-55MD, color video printer for ultrasound and endoscopy, incl. remote control, second-hand		980,00 EUR
	1004271	SONY UP-850: Video Graphic Printer b/w, second-hand		350,00 EUR
	1004357	SONY UP 1200 PM Mavigraph Color Video Printer, used		590,00 EUR
	1002536	SONY UP 1800EPM Mavigraph Colour Video Printer with digital display, second-hand		580,00 EUR
	1003919	SONY UP-2850P Color video printer, used		680,00 EUR
	1001770	SONY UP-2800P Color Video Printer, second-hand		680,00 EUR
	1002535	SONY UP 3030P Mavigraph Colour Video Printer, second-hand		680,00 EUR
<b>Dental treatment</b>	1001779	CAVITRON inserts for ultrasound plaque removal, pointed, new		93,00 EUR

1004279	KAVO 640 C Super-Torque: turbine, used		345,00 EUR
1004407	SATELEC Suprasson ZEG handpiece, used		145,00 EUR
1004446	Dental X-ray equipment: please see at X-ray equipment		0,01 EUR

## General Terms and Conditions of Business

of Dr. Wilfried Müller GmbH, Lechstraße 13 - 15, 86931 Augsburg, Germany

Dr. Wilfried Müller GmbH (hereinafter referred to as "Dr. Wilfried Müller GmbH") offers its customers merchandise for sale over the Internet through its own online shop.

1.

### Preliminary remarks

1. The deliveries, services and offers of Dr. Wilfried Müller GmbH are provided exclusively on the basis of these General Terms and Conditions.
2. These General Terms and Conditions only apply for companies, legal entities under civil law or special funds under civil law.
3. Terms and conditions of business or purchase of the purchaser opposing or differing from these General Terms and Conditions are herewith explicitly rejected, unless Dr. Wilfried Müller GmbH explicitly agrees with their application in writing. Such an agreement shall only apply to the respective individual case, and not to previous or future services.
4. These General Terms and Conditions shall also apply, when Dr. Wilfried Müller GmbH is aware of the terms and conditions of the purchaser contradicting or differing from these Terms and Conditions, but still performs the delivery to the customer without reservation.
5. Departures from these Terms and Conditions of Sale require the explicit written consent of Dr. Wilfried Müller GmbH.

2.

### Object of the agreement

1. The object is solely the supply of the merchandise or services which are defined in the description of the respective article. Drawings, images, measurements and other performance data are only binding after written confirmation by Dr. Wilfried Müller GmbH.
2. Contractual services which are not included in the article description, in particular those which are published for marketing purposes, only form part of the contractual subject matter if Dr. Wilfried Müller GmbH confirms these in writing.
3. Only the features included in the order confirmation shall be accepted as standard properties of the purchase object.

3.

### Entry into contract

a.

#### Onlineshop [www.drwmuellergmbh.de](http://www.drwmuellergmbh.de) or [www.med-geraete.de](http://www.med-geraete.de)

1. Dr. Wilfried Müller GmbH presents the articles in the online shops [www.drwmuellergmbh.de](http://www.drwmuellergmbh.de) / [www.med-geraete.de](http://www.med-geraete.de) in a non-obligatory and non-binding manner; this means that it is an invitation to the customers to submit an order with a binding offer to Dr. Wilfried Müller GmbH. The same applies to articles in newsletter promotions or in catalogues.
2. The customer is obliged to supply its contact details in full and to confirm its status as a company during the order process on [www.drwmuellergmbh.de](http://www.drwmuellergmbh.de) or [www.med-geraete.de](http://www.med-geraete.de).
3. The customer submits a binding order for the articles in the virtual shopping cart by confirming (clicking) the "Chargeable order" button at the end of the order process. Dr. Wilfried Müller GmbH shall promptly confirm receipt of the order by e-mail. This confirmation of receipt does not in itself represent an acceptance of the order; however, it may be combined with the declaration of acceptance.
4. In the event that the customer notices input errors, it has the option to correct them by using the browser functions (e.g. the Back button) before sending the order.
5. The text of the agreement is saved, and can be sent by e-mail upon request.
6. Dr. Wilfried Müller GmbH is entitled to accept the contract offer submitted with the customer's order within 14 days of receipt. The acceptance can be declared in form of a specific notification or the delivery of the merchandise.

b.

#### eBay

1. Dr. Wilfried Müller GmbH also offers merchandise in auctions on the Internet platform "eBay".
2. During an auction a purchase contract comes into effect with the party who has made the highest bid in the auction. The amount of the purchase price is determined by this highest bid. A bid expires when another bidder makes a higher bid during the auction.
3. If an offer also has the "Buy now" option, a purchase contract for the acquisition of the article comes into effect regardless of the end of the offer time and without conducting an auction based on the "Buy now price", when a purchaser clicks on the "Buy now" button and then confirms it. The purchaser can exercise this option, as long as no offer has been made for the article or the reserve price set by Dr. Wilfried Müller GmbH has not yet been reached. The purchaser can also obtain offers with the "Buy now" option as a guest without the need to register on eBay. In this case the purchaser accepts the offer by clicking on the "Buy now" button, selecting the "Buy as a guest" option on the following login page under "Continue without registration", and then completing the subsequent payment process.
4. If Dr. Wilfried Müller GmbH activates a "Suggested price" function for the article, the price for this article is negotiable. If the purchaser suggests a price, this becomes a binding offer. If Dr. Wilfried Müller GmbH declares its agreement with the suggested price, this serves as the corresponding binding acceptance by Dr. Wilfried Müller GmbH.
5. Dr. Wilfried Müller GmbH sends confirmation and details about the purchase process after entering into the contract.

c.

#### Other (e.g. sale on site, telephone enquiry)

1. Dr. Wilfried Müller GmbH is entitled to accept a customer order by sending an order confirmation or by delivering the merchandise within a period of two weeks. However, any confirmation of receipt Dr. Wilfried Müller GmbH that may send is not a binding acceptance of the order.
2. Offers, even those made in the name of Dr. Wilfried Müller GmbH, are subject to change and non-binding. Any contract is not considered legally binding and entered into, until Dr. Wilfried Müller GmbH has confirmed it in writing or executed it.
3. All agreements entered into between the vendor and the customer when concluding the contract require the written form. The vendor's members of staff are not entitled to agree changes or amendments to the contract without respecting the written form. If there be typing or printing errors in the vendor's order confirmation, or if the price be based on technical transmission errors, the vendor is entitled to contest these. Previously made payments shall be promptly refunded to the customer.

4.

### Prices, payment, offsetting

1. The prices quoted by Dr. Wilfried Müller GmbH are posted as net prices on the offer page in connection with the order process together with VAT at the applicable rate (currently 19%) and excluding packaging costs. VAT is not charged for deliveries to European countries, if the order contains a valid VAT ID. For deliveries outside of Europe no VAT is levied. The customer must observe the national tax regulations.
2. The purchaser must also bear the costs of delivery. Due to the variety of the offered merchandise these cannot be immediately determined. Dr. Wilfried Müller GmbH shall forward the information on the specific costs of delivery, as soon as it knows these.
3. Unless agreed otherwise in writing, the invoices of Dr. Wilfried Müller GmbH are immediately due for payment without a discount.
4. Deliveries to other countries are only made after advance payment. In this case, the ordered merchandise is only dispatched after receipt of the invoiced amount.
5. The purchaser falls into arrears, as soon as the contractually agreed payment term is exceeded.
6. When the purchaser falls into arrears, Dr. Wilfried Müller GmbH is entitled to charge interest amounting to 9% above the base interest rate from this moment on. In the event that Dr. Wilfried Müller GmbH proves higher damages caused by the delay, it shall be entitled to claim for it.
7. In the event that the purchaser defaults on the payment of deliveries or services, which are based on the same legal relationship, Dr. Wilfried Müller GmbH shall be entitled to demand advance payment and to withhold any merchandise which has not yet been dispatched or services which have not yet been provided - without

**Dr. Wilfried Müller GmbH - Lechstr. 13-15 - D-86931 Prittriching**

Phone: 0049-(0)8206-9624-0 | Fax: 0049-(0)8206-9624-29

E-Mail: [info@med-geraete.de](mailto:info@med-geraete.de) | Internet: <http://www.med-geraete.de>

13.03.2024

Page 119 of 121

- prejudice to any other rights.
8. The purchaser is only entitled to offset, if the counterclaims are legally binding, recognised by Dr. Wilfried Müller GmbH or if they are not disputed. Counterclaims arising from the same contractual relationship are excluded from this exclusion of offsetting.
9. The purchaser is only entitled to exercise the retention right, if its counterclaim is based on the same contractual relationship.
- 5.

### Delivery time, transfer of risk

- Compliance with the delivery obligation is conditional on the timely and proper fulfilment of the purchaser's obligations. The defence of breach of contract is reserved. The customer undertakes to provide correct and complete address details. In the event that incorrect details incur additional shipment costs, e.g. for a second shipment, these shall be borne by the customer.
  - In the event that the purchaser is in default of payment or of acceptance, the risk of an accidental transfer or an accidental impairment of the merchandise is transferred to the purchaser at the time, when the purchaser enters into default of payment or acceptance.
  - The agreed delivery times are always considered non-binding, unless explicitly agreed otherwise in writing. The agreed delivery time is observed, if Dr. Wilfried Müller GmbH communicates the production and readiness for collection of the merchandise to the purchaser by the end of the delivery time, unless discharge of debt at creditor's domicile or at the debtor's domicile have been exceptionally agreed.
  - Dr. Wilfried Müller GmbH is entitled to make part deliveries and perform partial services, unless acceptance of such a part delivery or part service cannot be reasonably expected from the purchaser.
  - In the event that the merchandise is sent to this address or to another address provided by the purchaser on the purchaser's request, the risk of the accidental loss or impairment of the merchandise is transferred to the customer, as soon as the merchandise has been handed over by the transport agent or have left the warehouse of Dr. Wilfried Müller GmbH for shipment. This applies irrespective of whether the merchandise is shipped from the place of fulfilment and which party bears the transport costs.
  - Operational disruptions, disruption to power or raw material suppliers, transport disruptions, insofar as these events were unforeseeable, as well as strikes, lock-outs, government decrees and cases of force majeure free the party affected for the duration of the disruption and to the extent of its effect from obligations of delivery or acceptance. In the event that the delivery acceptance is delayed by more than three months, both parties are entitled to withdraw from the contract with regard to the amount affected by the disruption to delivery or acceptance.
  - For deliveries to foreign countries, the purchaser shall be obliged to comply with the applicable customs, import and other regulations. Any costs for customs, import, storage and any other fees incurred in connection with the delivery to a foreign country must be borne by the purchaser; any taxes incurred must also be paid by the purchaser.
- 6.

### Retention of ownership

- Dr. Wilfried Müller GmbH reserves the ownership of the supplied merchandise to secure all claims against the purchaser which may arise from the business relationship. In the event of breach of contract by the purchaser, in particular arrears of payment, Dr. Wilfried Müller GmbH is entitled to reclaim the merchandise. If Dr. Wilfried Müller GmbH reclaims the merchandise, this automatically signifies its withdrawal from the contract. Dr. Wilfried Müller GmbH is entitled to sell the merchandise after it has been reclaimed; the income from this sale is to be offset against the purchaser's liabilities - less reasonable sales costs.
  - The ownership of Dr. Wilfried Müller GmbH extends to new products created by processing of the merchandise subject to retention of title. In this case, the processing is performed for Dr. Wilfried Müller GmbH as the manufacturer. In the event that objects not belonging to Dr. Wilfried Müller GmbH are processed, combined or mixed with merchandise belonging to the latter, the Dr. Wilfried Müller GmbH shall obtain co-ownership in the proportion the invoiced value of the merchandise subject to retention of title to the invoiced values of third part materials.
  - The purchaser assigns all claims from the sale of merchandise subject to the retention of title including bills and cheques to Dr. Wilfried Müller GmbH to secure the corresponding claims pursuant to item 6.1. Dr. Wilfried Müller GmbH accepts this assignment. When selling merchandise of which Dr. Wilfried Müller GmbH is the co-owner, the assignment is limited to the share of the claims corresponding to its co-ownership.
  - As long as the purchaser is prepared and capable to fulfil its obligations vis-à-vis Dr. Wilfried Müller GmbH, it may use the merchandise of which Dr. Wilfried Müller GmbH is the owner or co-owner in the normal course of business and itself collect claims assigned to Dr. Wilfried Müller GmbH. Dr. Wilfried Müller GmbH is obliged not to collect the claims, as long as the purchaser fulfils its payment obligations with the collected proceeds, is not in arrears, and in particular no application for the initiation of insolvency proceedings has been filed or payments have been suspended. If the latter is the case, Dr. Wilfried Müller GmbH can demand that the purchaser discloses the assigned claims and their debtors to Dr. Wilfried Müller GmbH, provides all details necessary for the collection, hands over the corresponding documentation and communicates this assignment to the debtors.
  - The purchaser can only agree assignments as collateral, pledges and assignments of claims, including by means of a sale of claims, only with the prior written consent of Dr. Wilfried Müller GmbH.
  - If the value of the securities exceeds the secured claims by more than 10%, Dr. Wilfried Müller GmbH shall release the securities of its choice upon the purchaser's request.
- 7.

### Rights of the purchaser due to faults

- The mere presentation of the articles in the online shop is classed as the simple functional specification, but not as a guarantee for the consistency of the articles. Declarations of guarantee of a third party, e.g. manufacturer's warranties, remain untouched.
  - The purchaser must promptly after receipt of the merchandise, whether it complies with the contractually agreed properties and is suitable for the intended purpose.
  - In case of obvious faults, Dr. Wilfried Müller GmbH must be promptly notified in writing thereof. Notification of hidden faults must be promptly given in writing after they have been discovered. Otherwise, the delivery is regarded as accepted.
  - In the event that the merchandise is faulty the purchaser has the option to demand the removal of the defect in the form of subsequent fulfilment. There is no right to subsequent delivery when no replacement object is available, and Dr. Wilfried Müller GmbH cannot procure it from a third party. If the subsequent fulfilment in the meaning of § 440 S. 2 BGB (German Civil Code) fails, the purchaser is entitled to choose between withdrawing from the contract and reducing the purchase price. This also applies if Dr. Wilfried Müller GmbH seriously and finally refuses to provide the subsequent fulfilment. In the event that the fault is minor, the purchaser shall have the right to a reduced price.
  - In the event that a trader's discount is granted, the warranty is excluded with the exception of intent, gross negligence and personal injury.
  - A right of return within 14 days applies for items offered within the 'repository' section, starting from the date of receipt of the goods. Any further warranty rights are excluded. You will be deemed to have met the deadline if you dispatch goods within the 14-day period.
  - With the exception of cases of intent, gross negligence and personal injury, warranty claims regarding the delivered merchandise lapse within one year,
    - unless otherwise agreed in writing.
    - Endoscopes shall have a warranty period of six months, which differs from the above period of one year.
    - For resellers a warranty period of 14 days shall apply. (A reseller is a party which intends to resell the merchandise at the time of entry into the contract with the vendor).
- The term starts at the time determined in § 199 BGB. It starts at the latest at the end of the maximum periods determined in § 199 para. 3 and para. 4 BGB.
- In the event that the purchaser or a third party conducts improper repair work or makes changes to the merchandise, they cannot assert any warranty claims for defect for the merchandise or the resulting consequences.
  - Any transport, travel, work and material costs for the purpose of subsequent fulfilment are not borne by Dr. Wilfried Müller GmbH, insofar as these increase, because the merchandise has been transported to another place than the premises of the purchaser after their delivery, unless this transport is required for its intended use.
  - There is no entitlement for warranty claims for insignificant deviations from the agreed characteristics, insignificant impairment of usability, natural wear & tear or damage, which occur after the transfer of risk due to incorrect or negligent use, excessive use, faulty assembly or due to external influences, which could not be expected according to the contract. There is also no entitlement for material defect claims, if the purchaser does not observe the instructions on handling, maintenance and inspections as well as care for the merchandise in the instruction manual.
  - Warranty is excluded for batteries, lights and wearing parts.
- 8.



## Liability

1. All claims for compensation of the purchaser vis-à-vis Dr. Wilfried Müller GmbH are excluded irrespective of their legal basis, unless Dr. Wilfried Müller GmbH or its vicarious agents have acted purposefully or with gross negligence or at least have breached important contractual duties with minor negligence. Important contractual duties are those which protect significant legal entitlements of the contractual partner relevant to this contract, which must be granted to them under the contract according to its content and purpose. Contractual duties are also those duties, which make orderly execution of the contract possible in the first place, whereby the contractual partner relies and should be able to rely on compliance therewith.
2. In the case of gross negligence or a breach of an important contractual duty due to minor negligence, compensation is limited to the typical and foreseeable damage.
3. The liability for any damage, which is not incurred on the object of the contract, is excluded except for cases of intent and gross negligence.
4. The above is without prejudice to liability due to fraudulent behaviour of Dr. Wilfried Müller GmbH and for accepted warranties, for claims according to the Product Liability Act (ProdHaftG) as well as damages due to injury to life, limb or health.
5. In the event that the liability of Dr. Wilfried Müller GmbH is excluded or limited, this also applies to the personal liability of the members of staff, employers, representatives and vicarious agents of Dr. Wilfried Müller GmbH.
6. For all claims for compensation or reimbursement of fruitless expenditure in case of contractual or non-contractual liability, which are asserted against Dr. Wilfried Müller GmbH - except for cases of intent, gross negligence or personal injury - a period of limitation of one year applies. The term starts at the time determined in § 199 BGB. The term starts at the latest at the end of the maximum periods determined in § 199 para. 3 and para. 4 BGB.

9.

## Saving of the text of the contract

Dr. Wilfried Müller GmbH saves the text of the contract. Upon request Dr. Wilfried Müller GmbH shall send the order confirmation together with all order details and the General Terms and Conditions to the email or address provided by the purchaser.

10.

## Closing provisions, applicable law, place of jurisdiction

1. This contract is governed by German law under exclusion of the UN Sales Law, unless written agreements to the contrary have been made.
2. The exclusive place of jurisdiction for all disputes in connection with the business relationship is the domicile of Dr. Wilfried Müller GmbH. The latter is, however, entitled to initiate legal proceedings against the purchaser at its general place of jurisdiction. In the event that the purchaser does not have a general place of jurisdiction in Germany or in another EU member state, the place of jurisdiction for all disputes arising from and in connection with this contract is the domicile of Dr. Wilfried Müller GmbH.
3. The place of fulfilment is the domicile of the vendor.
4. If one or several of the provisions in these terms and conditions are or become ineffective, incomplete or in need of amendment, the effectiveness of all other contractual provisions shall not thereby be affected.
5. The contractual parties are in this case obliged to agree on a provision which comes closest to the financial purpose of the invalid provision or to fill any identified omission.

Issued: June 2016