

DRÄGER Tiberius, simple, robust anaesthesia device with drawer unit, O2 flush, ...

??? ?????: 1001468

1.750,00 €
????? ??????



??? ??????

DRÄGER Tiberius, simple, robust anaesthesia device with drawer unit, O2 flush, second-hand Surcharge for gas cylinder (10 l) mount, art. 1004512 or 1004511) The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.

-> ????? ??? ?????? ??? ?????? <-

