

**DRÄGER Fabius: Mobile anaesthesia machine with FlowMeter for O2, Air and N2O, double ...**

??? ?????: 1002554

**2.100,00 €**  
????? ??????



??? ??????

DRÄGER Fabius: Mobile anaesthesia machine with FlowMeter for O2, Air and N2O, double holders for Vapor, drawer unit, second-hand Surcharge for gas cylinder (10 l) mount art. No. 1004512 and 1004511 The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.

-> ????? ??? ??????? ??? ??????? <-

