

**DR. MÜLLER portable anesthesia unit with flow meter for O2, for hanging in standard ...**

??? ?????: 1003830

**2.040,00 €**  
????? ??????



??? ??????

DR. MÜLLER portable anesthesia unit with flow meter for O2, for hanging in standard rails or placement on the tabletop, incl. stainless steel carrying frame, used. This item is exclusively intended for veterinary medicine and export. The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.

-> ????? ??? ?????? ??? ?????? <-

