

VÖLKER mobile anaesthesia unit with isoflurane vapour, flowmeter and circular section, ...

??? ?????: 1004450

1.390,00 €
 ?????? ??????



??? ??????

VÖLKER mobile anaesthesia unit with isoflurane vapour, flowmeter and circular section, incl. trolley and bottle holder, second-hand The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.

-> ????? ??? ?????? ??? ?????? <-

