

Dr. Müller laboratory animal anaesthesia device: Portable table-top anaesthesia device ...

??? ?????: 1004573

2.340,00 €
????? ??????



??? ??????

Dr. Müller laboratory animal anaesthesia device: Portable table-top anaesthesia device mounted on a stainless steel supporting frame (WxHxD: 35x31x17 cm). Supplied with Vet. med vapour for isoflurane, NEW (a filling piece is required for filling, price EURO 89.00) flow meter (0.1 litres - 1 litre), two-way valve, pressure relief valve, used The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs amount to EUR 89.00. This conversion is only possible for devices used in veterinary medicine.

-> ????? ??? ?????? ??? ?????? <-

