

??? ?????: 1005058

????? ???? ??



DRÄGER Fabian, wall anesthesia machine, incl. Circular part with inhalation and expiration valve, soda lime container, OHMEDA TEC 5 VET isoflurane, pull-out writing surface, used. This article is intended for veterinary medicine and export only. The standard inlet pressure for anaesthesia units is approx. 5 bar. If you would like to operate your device with an oxygen concentrator, we must convert the device to concentrator operation. The costs are EUR 89.00. This conversion is only possible for devices used in veterinary medicine.

-> ???? ??? ?????? ??? ?????? <-

