

**DeVilbiss oxygen concentrator, output  
5L/min, dimensions: HxWxD:  
62.2x34.2x30.4 cm, ...**

??? ?????: 1003618

**790,00 €**  
????? ??????



??? ??????

DeVilbiss oxygen concentrator, output 5L/min, dimensions: HxWxD: 62.2x34.2x30.4 cm, weight: 16.3kg, noise level: 40 dbA (50Hz), flow rate (continuously variable): 0.5-5L/min, oxygen concentration: 93% +/- 3%(0.5-5L/min), operating voltage: 230V/50Hz, power consumption: with 2L approx. 290 Watt / with 5L approx. 312 Watt, outlet pressure 0.59bar, flow rate 5 liter with 93Vol. %, ACCESSORIES (not to be used in connection with an anesthesia machine!), including bubble humidifier and 15 m connection hose, NEW

-> ????? ??? ?????? ??? ?????? <-

