

DRÄGER Babylog 2000: Ventilator for premature and newborn infants with a bodyweight of ...

??? ?????: 1000805

1.980,00 €
????? ??????



??? ??????

DRÄGER Babylog 2000: Ventilator for premature and newborn infants with a bodyweight of up to 6 kg, with integrated gas supply mixer function - oxygen and compressed air. For the following ventilation modes: controlled ventilation IPPV, intermittent mandatory ventilation IMV, spontaneous ventilation under continuous positive airway pressure CPAP With monitoring of positive airway pressure PAW, electricity supply to the monitoring module and the gas supply, second-hand

-> ????? ??? ?????? ??? ?????? <-

