

DRAEGER Babylog 8000: respirator for premature, newborn and Children up to 20kg body ...

??? ?????: 1000806

4.900,00 €
????? ??????



??? ??????

DRAEGER Babylog 8000: respirator for premature, newborn and Children up to 20kg body weight, controlled ventilation (IPPV), synchronized ventilation (SIPPV), intermittent mandatory ventilation (IMV), synchronized intermittent mandatory ventilation (SIMV), spontaneous breathing under positive airway pressure (CPAP), monitoring for insp. O₂ concentration FiO₂, airway pressure Paw, Flow V, minute volume MV, respiratory volume VT, dimensions basic device: (WxHxD) 212x280x390mm, 14.5kg, used

-> ????? ??? ?????? ??? ?????? <-

