

**DRÄGER Babylog 8000 plus:  
Ventilator for premature babies,  
newborn infants and children ...**

??? ?????: 1001536

**5.200,00 €**  
????? ??????



??? ??????

DRÄGER Babylog 8000 plus: Ventilator for premature babies, newborn infants and children up to 20 kg in bodyweight, controlled ventilation (IPPV), synchronised ventilation (SIPPV), intermittent mandatory ventilation (IMV), synchronised intermittent mandatory ventilation (SIMV), spontaneous respiration with continuous positive airway pressure (CPAP), monitoring for insp. O<sub>2</sub> concentration FiO<sub>2</sub>, positive airway pressure PAW, Flow V, minute volume MV, respiratory volume TV, dimensions of basic device: (WxHxD) 212x280x390 mm, 14.5kg, second-hand

-> ????? ??? ??????? ??? ??????? <-

