

**SIEMENS Servo 900 C ventilator,
ventilation pattern: volume control rt
,volume control ...**

??? ?????: 1001554

1.980,00 €
????? ??????



??? ??????

SIEMENS Servo 900 C ventilator, ventilation pattern: volume control rt ,volume control with sighs, pressure-controlled, pressure-supported, SIMV, SIMV + pressure support, CPAP manual, min. vol. 0,5-40 l/min, respiratory rate: 5-120 breaths/min, similar display: ventilatory pressure, exp. minute vol., dig. display: respiratory rate, O2 concentration, insp. tidal volume, exp. tidal volume, exp. minute volume, peak, intermittent, moderate, max. working pressure: 120 mbar, 50x23x32cm, weight : 19kg , from bodyweight of 10 kg , second-hand

-> ????? ??? ?????? ??? ?????? <-

