

**DRÄGER Oxylog 2000 transport ventilator with accessories, time-controlled, constant ...**

??? ?????: 1001559

**2.450,00 €**  
????? ??????



??? ??????

DRÄGER Oxylog 2000 transport ventilator with accessories, time-controlled, constant volume for patients with a tidal volume from 1000 ml ventilation mode: IPPV \ SIMV \ CPAP. Variables 1:E, also IRV, with integrated PEEP valve 0-15 mbar. Air-Mix \ No-Air-Mix switch, integrated expiratory volume measurement with display, second-hand

-> ????? ??? ?????? ??? ?????? <-

