

DRÄGER Ventilog C: Anaesthesia ventilation device, pneumatic driver, time-controlled ...

??? ?????: 1003341

1.480,00 €
????? ??????



??? ??????

DRÄGER Ventilog C: Anaesthesia ventilation device, pneumatic driver, time-controlled ventilation rate: 6 - 60 / min I:E ratio: 1:1, 1:2, 1:3 Inspiration flow: 20 to 80 L/min at 20 mbar stroke volume 30 to 150 ml (bellows for children) 100 to 1400 ml (bellows for adults) minute volume max 25 L/min, compressed gas supply up to 2 bar (O₂ or oil-free compressed air), 16 L/min, used. The small bellows are not included in the price, only for operation with COSY SYSTEM, used

-> ????? ??? ?????? ??? ?????? <-

