

**MELAG Euroklav 23VS+, chamber size
45 cm x diameter 25 cm, 22 liter, load
capacity 4 kg ...**

??? ?????: 1003751

3.650,00 €
????? ??????



??? ??????

MELAG Euroklav 23VS+, chamber size 45 cm x diameter 25 cm, 22 liter, load capacity 4 kg of instruments, 1 kg textiles, dimensions (WxHxD): 42.5 x 48.5 x 69 cm, the unit fits on a 60 cm tabletop, weight 45 kg, power consumption 2600 Watt, used

-> ????? ??? ?????? ??? ?????? <-

