

**WEBECO A65, steam steriliser with
integrated steam generator.
Cylindrical chamber to ...**

??? ?????: 1005297

10.900,00 €
????? ??????



??? ??????

WEBECO A65, steam steriliser with integrated steam generator. Cylindrical chamber to accommodate 1 sterilisation unit (STE). Equipped as standard with 5 programmes (textile programme, programme for thermolabile products, quick programme, special programme, prion programme) and 2 test programmes. Usable space 340x340x 600mm, 75 litres usable space. Weight 330kg, second-hand

-> ????? ??? ?????? ??? ?????? <-

