

**OLYMPUS CLV 160, xenon cold light source 300W, good condition, new xenon lamp, used ...**

??? ?????: 1001962

**2.900,00 €**  
 ????? ?????



??? ?????

OLYMPUS CLV 160, xenon cold light source 300W, good condition, new xenon lamp, used The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.

-> ????? ??? ?????? ??? ?????? <-

