

# OLYMPUS CLV-U20, Xenon cold light source 300 W with insufflation, only (OES mode, no OES ...

??? ?????: 1001963

**980,00 €**  
????? ??????



??? ??????

OLYMPUS CLV-U20, Xenon cold light source 300 W with insufflation, only (OES mode, no OES 100 or OES 200), second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.

-> ????? ??? ?????? ??? ?????? <-

