

**OLYMPUS CLV-S 20 cold light source
, 220 V/300 W,xenon, second-hand
The cold light ...**

??? ?????: 1002333

680,00 €
????? ??????



??? ??????

OLYMPUS CLV-S 20 cold light source , 220 V/300 W,xenon, second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.

-> ????? ??? ?????? ??? ?????? <-

