

OLYMPUS CLV-S40 Visera: Xenon cold-light source for rigid endoscope, 300 watt, used ...

??? ?????: 1002474

980,00 €
????? ??????



??? ??????

OLYMPUS CLV-S40 Visera: Xenon cold-light source for rigid endoscope, 300 watt, used The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.

-> ????? ??? ?????? ??? ?????? <-

