

OLYMPUS CLV-S30, xenon cold light source for rigid endoscopy, new xenon lamp, new lamp ...

??? ?????: 1002646

780,00 €
????? ??????



??? ??????

OLYMPUS CLV-S30, xenon cold light source for rigid endoscopy, new xenon lamp, new lamp life meter, used The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.

-> ????? ??? ?????? ??? ?????? <-

