

OLYMPUS CLV-180: Cold light source, xenon 300 Watt, including new bulb, second-hand The ...

??? ?????: 1004163

3.900,00 €
 ????? ?????



??? ?????

OLYMPUS CLV-180: Cold light source, xenon 300 Watt, including new bulb, second-hand The cold light source is equipped with a new xenon illuminant. The average service life is approx. 500 hours. The service life indicator of the illuminant is currently not available. This is only used for approximate determination of the service life.

-> ???? ??? ?????? ??? ?????? <-

