

Wisap CO2 Aqua-Purator 1611. The device provides continuous flushing of the abdomen and ...

??? ?????: 1002131

480,00 €
 ?????? ??????



??? ??????

Wisap CO2 Aqua-Purator 1611. The device provides continuous flushing of the abdomen and aqua dissection during a diagnostic laparoscopy. The accessory (bottle pin for NACL) is not included in the delivery and must be ordered directly from the manufacturer. Technical data: continuously adjustable pressure up to 750 mmHg. NACL instillation quantity: selectable 1 litre or Hi-Flo 3 l/min, second-hand

-> ????? ??? ?????? ??? ?????? <-

