

**OLYMPUS OTV-SI: video processor
with integrated cold light source 150
Watt halogen ...**

??? ?????: 1003939

2.900,00 €
????? ??????



??? ??????

OLYMPUS OTV-SI: video processor with integrated cold light source 150 Watt halogen (MD-151), for connection of video endoscopes: ENF-V2, ENF-VT2, ENF-VQ, or of camera heads: OTV-S7H-N, OTV-S7H-1N. OTV-S7H-1D, used

-> ????? ??? ?????? ??? ?????? <-

