

**ERBE Erbotom T175VET, monopolar and bipolar, HF power: monopolar cutting 175 watts, ...**

??? ?????: 1000499

**1.100,00 €**  
????? ??????



??? ??????

ERBE Erbotom T175VET, monopolar and bipolar, HF power: monopolar cutting 175 watts, monopolar coagulation 100 watts, bipolar 50 watts, both monopolar and bipolar functions can be activated without foot switch, branded ATTENTION: This device is intended exclusively for use in veterinary medicine. ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODES, 3 DIFFERENT INSERTS)

-> ????? ??? ?????? ??? ?????? <-

