

ERBE ICC 200 VET without EndoCut function (button present on operating panel, but ...

??? ?????: 1000505

1.980,00 €
????? ??????



??? ??????

ERBE ICC 200 VET without EndoCut function (button present on operating panel, but function not connected). Monopolar 200 watt, bipolar 120 watt. Very good condition, used. ALL HF SURGICAL UNITS ARE SUPPLIED WITH ACCESSORIES (NEUTRAL PLATE CABLE FOR SINGLE NEUTRAL ELECTRODES, MONOPOLAR AND BIPOLAR ELECTRODE, 3 DIFFERENT INSERTS). - This item is exclusively designed for being used in the field of veterinary medicine or for export purposes -

-> ????? ??? ??????? ??? ??????? <-

