

BERCHTOLD Elektrotom 621: HF surgical device with the following functions: Monopolar type ...

??? ?????: 1000516

2.390,00 €
????? ??????



??? ??????

BERCHTOLD Elektrotom 621: HF surgical device with the following functions: Monopolar type of cutting current with low sloughing (max 200 W) and moderate sloughing (max 150 W), monopolar haemostatic type of cutting current (polypectomy, papillotomy) (max 100 W), monopolar contact-coagulation (max 120 W), forced coagulation (max 150 W) and spray coagulation (max 90 W). Bipolar coagulation (max 95 W). Digital performance display, service display, memory (last setting saved at shut-down), used ALL HF SURGICAL DEVICES ARE DELIVERED WITH ACCESSORIES (NEUTRAL PLATE, MONO- AND BI-POLAR ELECTRODE, 3 DIFFERENT INSERTS)

-> ????? ??? ?????? ??? ?????? <-

