

ERBE APC 300, Argon Beamer for argon plasma coagulation. Only for use in combination with ...

??? ?????: 1002689

1.590,00 €
????? ??????



??? ??????

ERBE APC 300, Argon Beamer for argon plasma coagulation. Only for use in combination with ERBOTOM ACC 450, 451, ICC 200 Argon Coag, 300, 350. For argon plasma coagulation (coagulation method in which argon gas which has been ionized by HF current is passed onto the tissue to be coagulated), without pressure reducer, used.

-> ????? ??? ?????? ??? ?????? <-

