

**DRÄGER Photo therapy 800, with chassis, for the light irradiation of premature or newly ...**

??? ?????: 1001357

**590,00 €**  
????? ??????



??? ??????

DRÄGER Photo therapy 800, with chassis, for the light irradiation of premature or newly born infants in the incubator or in the infrared-heated thermal bed to reduce the serum-bilirubin concentration (phototherapy), frequency: 50 Hz, output: 230W, current: 2.2 A, lamp housing 54 x 25.5 x 20.5cm, second-hand

-> ????? ??? ?????? ??? ?????? <-

