

**DRÄGER Isolette 8000plus, mobile incubator for newborns. Continuous monitoring of ...**

??? ?????: 1005288

**4.500,00 €**  
????? ??????



??? ??????

DRÄGER Isolette 8000plus, mobile incubator for newborns. Continuous monitoring of central and peripheral skin temperature. Manual or automatic humidity control. Height adjustment from 133cm to 153cm. Width 104cm, depth 76cm. Weight approx. 98kg. Access flaps on front and back, 6 access openings. Mattress size: 38x74x3cm, used.

-> ????? ??? ?????? ??? ?????? <-

