

BRAUN Perfusor Compact, very good condition, delivery rate 0.1-99.9 ml/h, bolus rate 800 ...

??? ?????: 1001335

690,00 €
 ?????? ???????



??? ???????

BRAUN Perfusor Compact, very good condition, delivery rate 0.1-99.9 ml/h, bolus rate 800 ml/h, new battery, second-hand
 This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please!
 For retroactive changement (after sale) we have to charge EUR 55,-

-> ????? ??? ??????? ??? ??????? <-

