

BRAUN Perfusor fm, very good condition, delivery rate 0.1 to 200 ml/hour, with new ...

??? ?????: 1001337

490,00 €
 ?????? ???????



??? ???????

BRAUN Perfusor fm, very good condition, delivery rate 0.1 to 200 ml/hour, with new battery and integrated mains component, second-hand This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please! For retroactive changement (after sale) we have to charge EUR 55,-

-> ????? ??? ??????? ??? ??????? <-

