

**BRAUN Perfusor Compact S, volume:
0.01 - 200 ml/hour, time selection
0.01-99h, max. ...**

??? ?????: 1002673

750,00 €
????? ??????



??? ??????

BRAUN Perfusor Compact S, volume: 0.01 - 200 ml/hour, time selection 0.01-99h, max. infusion rate: 200 ml/h, max. bolus rate: 1200 ml/h, new battery, used This perfusor is ready for use with B. BRAUN perfusor syringes. If you want to use different syringe types let us know please! For retroactive changement (after sale) we have to charge EUR 55,-

-> ????? ??? ?????? ??? ?????? <-

