

CAREFUSION Alaris plus GH, syringe pump, compatible with many standard disposable ...

??? ?????: 1005117

690,00 €
????? ??????



??? ??????

CAREFUSION Alaris plus GH, syringe pump, compatible with many standard disposable syringes 50ml (Braun Omnifix, BD, Codan, Fresenius, Terumo). Max. Infusion rate with 50ml. Syringes 0.1ml/h-1200ml/h. Max. Bolus rate 10ml/h-1200ml/h. System accuracy +/-2%. Rechargeable NIMH battery to operate the device for approx. 4 hours, used

-> ????? ??? ?????? ??? ?????? <-

