

BAXTER EVO IQ LVP VET, volumetric infusion pump for mains or battery operation. Rate ...

??? ?????: 1005170

590,00 €
????? ???????



??? ???????

BAXTER EVO IQ LVP VET, volumetric infusion pump for mains or battery operation. Rate 0.1-2000ml/h, bolus rate 0.1-2000ml/h. We recommend the use of B.Braun Intrafix Safeset infusion lines, used. This article is exclusively intended for use in veterinary medicine or export.

-> ????? ??? ??????? ??? ??????? <-

