

**FRESENIUS Injectomat Agilia.
Syringe pump for battery and mains
operation, incl. ...**

??? ?????: 1005251

640,00 €
????? ??????



??? ??????

FRESENIUS Injectomat Agilia. Syringe pump for battery and mains operation, incl. retaining claw for infusion stand. Approx. 50 usable syringe types. Delivery rate 0.1 to 99.9ml/h in 0.1ml steps. Bolus rate 50-1200ml, second hand

-> ????? ??? ?????? ??? ?????? <-

