

**FRESENIUS Injectomat TIVA Agilia.
Syringe pump for battery and mains
operation, incl. ...**

??? ?????: 1005291

980,00 €
????? ??????



??? ??????

FRESENIUS Injectomat TIVA Agilia. Syringe pump for battery and mains operation, incl. holding claw for infusion stand. Approx. 50 usable syringe types. Delivery rate 0.1 to 99.9ml/h in 0.1ml increments. Bolus rate 50-1200ml. Please note that the medication database must be programmed by Fresenius field service. second hand

-> ????? ??? ?????? ??? ?????? <-

