

DRÄGER IRIS: anaesthetic monitor to measure the concentration continuously and monitor ...

??? ?????: 1001817

580,00 €
????? ??????



??? ??????

DRÄGER IRIS: anaesthetic monitor to measure the concentration continuously and monitor the inhalation anaesthetics halothane, enflurane and isoflurane in the fresh gas supply from inhalation anaesthesia devices. Measurement range 0 to 9.9 Vol. %, measurement limits can be adjusted individually, second-hand

-> ????? ??? ?????? ??? ?????? <-

