

DRÄGER Oxydig: O2 measurement and warning device used in oxygen tents, incubators, ...

??? ?????: 1001821

620,00 €
????? ??????



??? ??????

DRÄGER Oxydig: O2 measurement and warning device used in oxygen tents, incubators, ventilators and anaesthetic devices, etc. Measuring range 0-100 Vol.% O2, spiral sensor cable can be extended to a maximum length of 1.5 m , dimensions of the display instrument 170 x 93 x 65 mm, weight 0.5 kg, sensor dimensions: diameter 40 mm, H 45 mm, alarm is activated visually and acoustically when the individual adjustable warning limits are exceeded, second-h and

-> ????? ??? ?????? ??? ?????? <-

