

**DR. MACH Trigenflex R96, 3 x 50
watts, approx. 75,000 lux, focusable.
Diameter of lamp ...**

??? ?????: 1000656

1.390,00 €
????? ??????



??? ??????

DR. MACH Trigenflex R96, 3 x 50 watts, approx. 75,000 lux, focusable. Diameter of lamp body 48cm, light field: 15 - 26 cm, ceiling or wall mounting, with sterilisable handle sleeve, incl. transformer, second-hand

-> ????? ??? ??????? ??? ??????? <-

