

SMITH EQ-5000 Level 1 Equator, convective patient warmer. Adjustable temperature range: ...

??? ?????: 1005201

590,00 €
????? ??????



??? ??????

SMITH EQ-5000 Level 1 Equator, convective patient warmer. Adjustable temperature range: 36°C, 40°C, 44°C. Delivery incl. 2-metre supply hose, second-hand. We recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: vet@moeckundmoeck.de

-> ????? ??? ?????? ??? ?????? <-

