

SMITH Level 1, convective patient warmer with adjustable temperature levels: 37°C, ...

??? ?????: 1005247

650,00 €
????? ??????



??? ??????

SMITH Level 1, convective patient warmer with adjustable temperature levels: 37°C, 40°C, 43°C. Delivery incl. 2 metre supply hose, second-hand. We recommend the use of convective warming pads from Moeck&Moeck. These are compatible with almost all patient warming systems. Moeck&Moeck also offers special warming pads for veterinary medicine. Please request a quotation directly from Moeck&Moeck: vet@moeckundmoeck.de

-> ????? ??? ?????? ??? ?????? <-

