

## **DRÄGER Fabius GS, anaesthesia workstation with electronically controlled and operated ...**

Article number: 1003366

**11.900,00 €**  
net price



### **Product description**

DRÄGER Fabius GS, anaesthesia workstation with electronically controlled and operated ventilator, 45-minute battery back-up in the event of general power supply failure. ventilation mode: manual, spontaneous, IPPV, PCV. Dimensions: Width: 89.5cm, x height: 130 cm, x depth: 82 cm, electronic measurement of fresh gas flow with semi-closed COSY breathing system, integrated monitoring: insp. O<sub>2</sub> conc., breathing rate, tidal volume, AMV, mean and plateau pressure, peak respiratory pressure and PEEP, fresh gas flow, In addition this Fabius GS is provided with the following monitoring: Dräger Infinity Gamma XL with the following functions: ECG, NIBP, SpO<sub>2</sub>, CO<sub>2</sub>, NIBP, PRESS, Temp and 5 anaesthesia gases (halothane, enflurane, isoflurane, desflurane and sevoflurane), inclusiv DRÄGER Vaporiser 2000, Sevoflurane, with DRÄGER fill system, DRÄGER support, DRÄGER Vapor 2000, Isoflurane, with DRÄGER-FILL filling unit, used

-> [Click here to come to the site.](#) <-

