

**DeVilbiss oxygen concentrator, output
10L/min, dimensions: HxWxD:
62.2x34.2x30.4 cm, ...**

Article number: 1004313

1.490,00 €
net price



Product description

DeVilbiss oxygen concentrator, output 10L/min, dimensions: HxWxD: 62.2x34.2x30.4 cm, weight: 16.3kg, noise level: 40 dbA (50Hz), flow rate (continuously variable): 0.5-5L/min, oxygen concentration: 93% +/- 3%(0.5-5L/min), operating voltage: 230V/50Hz, power consumption: with 2L approx. 290 Watt / with 5L approx. 312 Watt, outlet pressure 0.59bar, flow rate 5 liter with 93Vol. %, ACCESSORIES (not to be used in connection with an anaesthesia machine!) including bubble humidifier and 15 m connection hose, NEW

-> [Click here to come to the site.](#) <-

