

## **DRAEGER Babylog 8000: respirator for premature, newborn and Children up to 20kg body ...**

Article number: 1000806

**4.900,00 €**  
net price



### **Product description**

DRAEGER Babylog 8000: respirator for premature, newborn and Children up to 20kg body weight, controlled ventilation (IPPV), synchronized ventilation (SIPPV), intermittent mandatory ventilation (IMV), synchronized intermittent mandatory ventilation (SIMV), spontaneous breathing under positive airway pressure (CPAP), monitoring for insp. O<sub>2</sub> concentration FiO<sub>2</sub>, airway pressure P<sub>aw</sub>, Flow V, minute volume MV, respiratory volume VT, dimensions basic device: (WxHxD) 212x280x390mm, 14.5kg, used

-> [Click here to come to the site.](#) <-

