

**DRÄGER Ventilog 2 ventilation device,
incl. accessories, powered with
oxygen or ...**

Article number: 1001528

980,00 €
net price



Product description

DRÄGER Ventilog 2 ventilation device, incl. accessories, powered with oxygen or compressed air, for controlled ventilation, constant volume, time-controlled, tidal volume 50 - 150 ml (small bellows) and 150 - 1600 ml (large bellows), frequency 6-60 min, used, Inspirat: flow 20-80 l/min, compressed gas supply: 2-6 bar, compressed gas consumption: 16 l/min, used, selectable I : E ratio = 1:1, 1:2, 1:3, used. The small bellows are not included in the price

-> [Click here to come to the site.](#) <-

